

**Vitrebond™ Light Cure Glass Ionomer Liner/Base  
Adper™ Adhesive Systems  
Filtek™ P60 Posterior Restorative  
Sof-Lex™ Discs, Strips and Finishing Brush**

Prepare tooth and isolate.

**Apply liner/base if desired:**

- Mix a level scoop of Vitrebond powder and a drop of Vitrebond liquid on a mixing pad.
- Apply a thin layer of the liner/base to dentin surfaces using a ball applicator.
- Light-cure for 30 seconds.

**Adhesive Placement:**

- Apply adhesive following directions for 3M™ ESPE™ Adper™ adhesive of choice.

**Place Restorative:**

- Place Filtek P60 restorative in increments less than 2.5mm.
- Light-cure each increment for 20 seconds.

**Finish and Polish:**

- Finish occlusal surface using a Sof-Lex Finishing Brush.
- Finish proximal surfaces with Sof-Lex Extra Thin Discs and Strips.



**Vitrebond™ Light Cure Glass Ionomer Liner/Base**  
**Adper™ Adhesive Systems**  
**Filtek™ P60 Posterior Restorative**  
**Sof-Lex™ Discs, Strips and Finishing Brush**

---



**Check occlusion:**

- Check lateral and centric occlusion.
- Adjust if necessary.