

# Silmet Luting Guide

## Glass Ionomer Cement

## RMGI

## Aesthetic Resin Cement

## Adhesive Resin Cement

## Self-Adhesive Resin Cement

Advantages	<ul style="list-style-type: none"> <li>High fluoride release, rechargeable</li> <li>Low chemical bond to tooth</li> <li>Adhere in wet environment</li> <li>Low film thickness</li> <li>Low technique sensitivity</li> </ul>	<ul style="list-style-type: none"> <li>Fluoride release, rechargeable</li> <li>Low- medium chemical and micromechanical bond to tooth</li> <li>Less technique sensitive than resin cement.</li> </ul>	<ul style="list-style-type: none"> <li>Highest bond strength</li> <li>Minimal shade shift over time' if light-cured</li> <li>Highly aesthetic</li> </ul>	<ul style="list-style-type: none"> <li>No etching to tooth structure is needed.</li> </ul>	<ul style="list-style-type: none"> <li>No etch or primer</li> <li>Easy to use</li> <li>Low postoperative sensitivity</li> <li>Less technique sensitivity</li> <li>Easy cleanup</li> </ul>
Disadvantages	<ul style="list-style-type: none"> <li>Sensitivity may occur if tooth is over-dried</li> <li>Time(24 hours) is required to maximum strength</li> <li>Low strength</li> <li>Water sensitive during setting phase</li> </ul>	<ul style="list-style-type: none"> <li>Moisture-sensitive technique</li> <li>Sensitivity may occur if tooth is over-dried.</li> </ul>	<ul style="list-style-type: none"> <li>Most technique sensitive</li> <li>Moisture sensitive.</li> <li>.Highest chance of postoperative sensitivity if used with total-etch bonding</li> </ul>	<ul style="list-style-type: none"> <li>Limited availability of shades</li> <li>May require oxygen inhibition gel.</li> </ul>	<ul style="list-style-type: none"> <li>Can have shade shift over time</li> <li>Not as strong as adhesive resin cements.</li> </ul>
Indications	<ul style="list-style-type: none"> <li>Metal and ceramic-metal restorations</li> <li>High -strength ceramic(zirconia) crowns and bridges</li> <li>Posts (metal)</li> <li>Laboratory composites</li> </ul>	<ul style="list-style-type: none"> <li>Metal and ceramic-metal restorations</li> <li>High -strength ceramic(zirconia) crowns and bridges</li> <li>Laboratory composites</li> <li>Implant restorations</li> </ul>	<ul style="list-style-type: none"> <li>All- ceramic crowns and veneers in aesthetic zone</li> </ul>	<ul style="list-style-type: none"> <li>All- ceramic crowns ,onlays ,inlays, bridges</li> <li>Metal or ceramic-metal crowns/bridges</li> <li>High-strength ceramic(zirconia) crowns ,onlays ,inlays, bridges</li> <li>Maryland bridges(metal wings)</li> <li>Posts ( cast metal, ceramic, fiber- reinforced resin).</li> </ul>	<ul style="list-style-type: none"> <li>All- ceramic crowns ,onlays ,inlays, bridges</li> <li>All -Metal or ceramic-metal crowns/bridges</li> <li>High -strength ceramic(zirconia) crowns ,onlays ,inlays, bridges</li> <li>Posts (meta and fiber).</li> </ul>
Contraindications	<ul style="list-style-type: none"> <li>All- ceramic restorations</li> </ul>	<ul style="list-style-type: none"> <li>Thin all ceramic restoration</li> <li>Posts</li> </ul>	<ul style="list-style-type: none"> <li>Opaque all-ceramic restoration</li> </ul>	<ul style="list-style-type: none"> <li>All- ceramic veneers</li> </ul>	<ul style="list-style-type: none"> <li>Ceramic veneers</li> <li>Crown or bridge with poor retention</li> <li>Resin -bonded bridges.</li> </ul>
Silmet's Brand	<b>PROGLASS™ One</b>	<b>PROGLASS™ Plus</b> <b>PROGLASS™ Plus CEM</b>	None	<b>PROLink™ CEM</b>	<b>PROLink™ CEM Plus</b>
Similar to:	Fuji 1™,	Fuji Cem 2™, Relyx Luting™ , Nexus RMG™I, Riva Cem	None	Maxcem Elite™, Relyx Unicem 2™,	RelyX ultimate™ , Panavia™ V5, G-Cem one™