**ACTIVA™ PRESTO™****PRODUCT DESCRIPTION**

**ACTIVA PRESTO™** is an esthetic, light cure material that contains calcium, phosphate, and fluoride in a durable, wear-resistant resin matrix indicated for all restorative procedures. The material is stackable, does not slump, and holds its shape prior to light curing.

**ACTIVA PRESTO** is tougher than traditional composites. It contains a patented, shock-absorbing, rubberized-resin component that resists chipping and fracture. This increases flexural strength while maintaining the esthetics, hardness, compressive strength, and durability required of dental composites.

**ACTIVA PRESTO** is a urethane-based resin that contains no Bisphenol A, no Bis-GMA, and no BPA derivatives.

**INDICATIONS**

Recommended as a filling material for pits, root surface cavities, and Class I, II, III, IV and V restorations where there is no pulpal involvement.

**CONTRAINDICATIONS**

Not indicated for direct placement on the exposed pulp.

**INSTRUCTIONS FOR USE**

1. Isolate and prepare tooth to receive a restoration. Bevel and round off enamel margins. In Class V lesions, bevel or undercut enamel.

2. Place appropriate pulp protection, if indicated.
3. Etch and apply a bonding agent following the manufacturer's instructions for use.
4. **ACTIVA PRESTO** is a dentin and enamel replacement material. Place a 19-gauge applicator tip on the syringe and apply in 2mm increments, light curing for 20 seconds between each layer.
5. Finish and polish as usual.

**STORAGE AND HANDLING**

- Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.
- Shelf life of unopened product: 2 years from date of manufacture.
- Re-cap immediately after use.

**CAUTION**

Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves. Do not exceed manufacturer's recommended curing time for the light you are using.

Note: Multi-dose syringes should be encased in a fresh protective barrier for each patient. For additional information, refer to: <http://www.fda.gov/Medical Devices/Products and Medical Procedures/Dental Products/ucm4044 72.htm>.

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