
Cartridge Delivery

- ① Bite registration Tip
- ② Mixing Tip
- ③ Dispensing Gun

Distributed in USA/Canada exclusively by:
Kulzer
300 Herkules Way
South Bend, IN 46614-2517
1-800-431-1785

Caution: Federal law restricts this device to sale by or on the order of a dentist.

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Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
Made in Romania

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Flexitime® Bite

Instructions for use GB

Scannable polyvinyl siloxane for bite registration, superior hardness, very short time in mouth.

Description

Flexitime® Bite is an addition-cross-linking polyvinyl siloxane with a very high final Shore D hardness of 40. Its surface properties and optimized color enable it to be optically recorded in CAD/CAM systems to portray the opposing dentition. Flexitime® Bite can be scanned in the Kulzer scanner and many other CAD/CAM systems. It complies with the specifications in DIN 13903.

- Uses**
• Standard bite registrations in the end bite position
• Portrayal of opposing dentition
• Coating of bite forks
• Key material for needle point registration
• Production of small model segments

Instructions for use

Please visually inspect for damage before use. Do not use damaged products. Mixing tips and bite registration tips are for single use only.

Dispense a small amount of material from a 50 ml cartridge before attaching a mixing tip. After use, leave the mixing tip on the cartridge to keep it closed.

Immediately after use, remove the mixing tip and attach a new one.

When bite registration materials based on polyvinyl siloxane are touched with latex gloves, retraction solutions or freshly inserted composite fillings, it can change the curing behavior.

Flexitime® Bite should be at room temperature when used (23°C (73°F)). Lower temperatures, for example when the materials are stored in a refrigerator can extend the working time, and higher temperatures can shorten it.

Hardened polyvinyl siloxanes are chemically resistant. Clothes can be permanently stained if the contact material is not completely cured.

Flexitime® Bite is used with Mixing Tip green (1:1). It can also be combined with the Bite Registration Tip.

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