
Cartridge Delivery

- 1 Bite registration Tip
- 2 Mixing Tip
- 3 Dispensing Gun 2

Distributed in USA/Canada exclusively by:
Kulzer, LLC
300 Hereaus Way
South Bend, IN 46614-2517
1-800-431-1785

Caution: Federal law restricts this device to sale by or on the order of a dentist.

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Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
Made in Romania

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Flexitime[®] Bite

Instructions for use GB

Scannable polyvinyl siloxane for bite registration, superior hardness, very short time in mouth.

Description

Flexitime Bite is an addition-cross-linking polyvinyl siloxane with a very high final Shore D hardness of 40. Its surface properties and optimized color enable it to be optically recorded in CAD/CAM systems to portray the opposing dentition. Flexitime Bite can be scanned in the Kulzer scanner and many other CAD/CAM systems. It complies with the specifications in DIN 13903.

Uses

- Standard bite registrations in the end bite position
- Portrayal of opposing dentition
- Coating of bite forks
- Key material for needle point registration
- Production of small model segments

Instructions for use

Please visually inspect for damage before use. Do not use damaged products. Mixing tips and bite registration tips are for single use only. Dispense a small amount of material from a 50 ml cartridge before attaching a mixing tip. After use, leave the mixing tip on the cartridge to keep it closed. Immediately before the next use, remove the mixing tip and attach a new one.
When bite registration materials based on polyvinyl siloxane are touched with latex gloves, retraction solutions or freshly inserted composite fillings, it can change the curing behavior.
Flexitime Bite should be at room temperature when used (23°C (73°F)). Lower temperatures, for example when the materials are stored in a refrigerator can extend the working time, and higher temperatures can shorten it.
Hardened polyvinyl siloxanes are chemically resistant. Clothes can be permanently stained if the contact material is not completely cured.
Flexitime Bite is used with Mixing Tip green (1:1). It can also be combined with the Bite Registration Tip.

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¹⁾ Total working time (including mixing) at 23°C (73°F); shortened to 20s at oral temperature (95°C (95°F))
²⁾ Minimum time that the bite registration material must be in a closed mouth.

The data given are furnished to the best of our knowledge and are based on internal test results.

Disinfection

Flexitime Bite can be treated with an aqueous sodium hypochlorite solution (5.25% for up to 10 minutes*). The treatment does not affect the surface or dimensions. *DGZMK 3/93 V1.0
When using other disinfectant solutions suitable for silicones, follow the disinfection and processing instructions of the manufacturer.

Storage conditions

Do not store above 25°C (77°F)! Do not use after the expiration date (see the information on the package).

Warning and safety instructions

If the material contacts your eyes, rinse with copious amounts of water for several minutes. If symptoms continue, consult a doctor.
Do not inhale, swallow or ingest. If health problems arise after swallowing impression material, immediately contact a doctor. In rare cases, intestinal blockage may arise. Immediately consult a doctor if the product is inhaled.
Do not leave any residual material in the sulcus or oral cavity.
Flexitime Bite should not be used for conventional impressions since its high final hardness makes it very difficult to remove from the mouth. Undesired tooth extractions may result.
Hypersensitive persons may experience an allergic reaction to the product.
Only to be used by qualified dental personnel.

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