

# MARK3<sup>®</sup>

Premium Products. Affordable Prices.

# MARK3 EverCem Plus

Universal Resin Cement (Automix) with MDP

EverCem Plus from MARK3 now features MDP for enhanced bonding, especially to zirconia restorations. This dual-cure, self-etch, self-adhesive resin cement eliminates separate etching, priming, and bonding steps to save chair time. It is ideal for cementing ceramic, composite, metal restorations, and posts. Its low film thickness helps ensure a closer restoration fit with excellent bond strength.

## Features & Benefits

- **Contains MDP for strong, durable bond**
- Dual-cure for versatile performance
- Self-adhesive – no separate bonding agent in most cases
- Low film thickness for precise seating
- Easy automix syringe delivery
- Radiopaque for clear x-ray identification

### Performance Snapshot

Compressive strength: 275–300 MPa  
Flexural strength: 75–125 MPa  
Shear bond strength: 28–32 MPa  
Film thickness: 12 µm



## Indications for Use

- All-ceramic crowns, inlays & onlays
- Metal & ceramic-metal crowns/bridges
- High-strength zirconia restorations
- Metal or fiber posts

**COMPARES TO** RelyX Ultimate, Kerr NX3, Panavia SA Universal, G-Cem One, Premier ZR-Cem

**PACKAGING** Contains (1) 8 gm Automix Syringe & 6 mixing tips

**Efficient universal cementation with fewer steps and reliable bond strength.**

MARK3 EverCem Plus Universal Resin Cement