


NEW CLEANING SOLUTIONS


Authorized for sale in Canada



NON-AMMONIATED General Purpose Ultrasonic Solution

- Advanced formulation enhances your ultrasonic machines' cleaning action by combining nonionic surfactants, powerful chelators and corrosion inhibitors.
- Helps protect and extends the life of delicate instruments.
- Use directly in the ultrasonic tank or as a sonic conductor when using MARK3 Tartar, Stain and Permanent Cement Remover.
- The 10:1 concentrate makes storage easier.


#0092 1 Gallon Liquid



ULTRASONIC SOLUTION Tartar, Stain & Permanent Cement Remover

- Quickly removes tarter build-up, calculus deposits, food stains and tobacco stains from appliances such as partials, dentures and crown and bridges.
- No mixing necessary

#0094 1 Gallon Liquid



MULTI-ENZYME BLEND Evacuation System Cleaner

- Concentrated, multi-enzyme blend deodorizes as it cleans
- Dissolves organic fluids and debris
- Non foaming
- Daily use improves efficiency
- Neutral pH for safe use with amalgam separators
- Environmentally friendly
- Pleasant peppermint scent

#0090 1 Gallon Liquid