


Ultra-speed



## The intraoral film that delivers high-contrast and fine detail.

Ultra-speed dental film is a D-speed film that generates the high-quality, detailed images you need to make accurate diagnoses. Ideal for both manual and automated processing methods, Ultra-speed film is a versatile and easy-to-use solution for your practice.

### Enhanced image quality.

Because it accommodates such a wide range of processing and handling conditions, Ultra-speed film is a simple solution for your intraoral imaging needs. Exceptionally tolerant of process variations, this innovative film is designed to render excellent image quality.

### The standard in safety.

Ultra-speed film is also available with a ClinAsept barrier feature to help reduce the risk of infection. The built-in barrier protects the film packet from contaminants, while

its innovative design frees up your time and ensures peace of mind. Simply open the barrier envelope, place the film into the cup and you're ready to start processing.

Available in multiple sizes and formats to meet all of your intraoral imaging needs:

- Available in sizes 0-4.
- Select from flexible, easy-to-open Super Poly-Soft vinyl packets or our sturdier paper products for sizes 0 and 2.
- Bite-wing packets – paper packets with bite-wing tabs attached.
- ClinAsept barrier packets for added safety.
- For optimum results, process your radiographs using Carestream Dental chemicals.

## Ordering Information for Ultra-speed

| Size   | Description  | Code   | Packets/Package | REF No   |
|--------|--|--------|-----------------|----------|
| Size 0 | Super Poly-Soft packets                                      | DF-54  | 100             | 122 8840 |
| Size 0 | Super Poly-Soft packets (double film)                        | DF-53  | 100             | 122 8931 |
| Size 0 | Super Poly-Soft packets with ClinAsept barrier               | DF-54C | 75              | 144 5360 |
| Size 0 | Bite-wing paper packets                                      | DF-48  | 50              | 834 8658 |
| Size 1 | Paper packets  | DF-56  | 100             | 127 3747 |
| Size 1 | Paper packets (double film)                                  | DF-55  | 100             | 127 3721 |
| Size 2 | Super Poly-Soft packets                                      | DF-58  | 150             | 165 8194 |
| Size 2 | Super Poly-Soft packets (double film)                        | DF-57  | 130             | 175 3664 |
| Size 2 | Paper packets  | DF-58  | 150             | 149 1737 |
| Size 2 | Paper packets (double film)                                  | DF-57  | 150             | 149 1752 |
| Size 2 | Bite-wing paper packets                                      | DF-40  | 50              | 839 3043 |
| Size 2 | Super Poly-Soft packets with ClinAsept barrier               | DF-58C | 100             | 171 7131 |
| Size 2 | Super Poly-Soft packets with ClinAsept barrier (double film) | DF-57C | 100             | 898 7877 |
| Size 3 | Bite-wing paper packets                                      | DF-42  | 100             | 129 6771 |
| Size 4 | Occlusal paper packets                                       | DF-50  | 25              | 166 6163 |

For more information, call **800.933.8031** or visit [www.carestreamdental.com](http://www.carestreamdental.com).