

MARK3™
THE SMART ALTERNATIVE



Non-Splatter Formula Prophylaxis Paste

FEATURES & BENEFITS

- Contains 1.23% APF
- Maximizes stain removal
- Minimizes enamel loss and helps maintain an effective fluoride level
- Strengthens tooth enamel and prevents dental caries
- Disposable cups eliminate cross contamination
- Non-Splatter formula
- Great flavor with no bitter aftertaste

The MARK3 Line
was created to provide
High Quality Products
at an Affordable Price

order#	grit/flavor
1620Coarse - Assorted
1621Coarse - Bubble Gum
1622Coarse - Cherry
1623Coarse - Mint
1630Medium - Assorted
1633Medium - Mint

Each box Includes a prophy ring for ease of application.

200 per box

Compare to Nupro