

Silmet's Selection of Composites

(Micro-Hybrid, Nano Technology)

Features

- Reduced sensitivity to ambient light
- Light cured, radiopaque & color stable
- Low polymerization shrinkage
- High strength & wear resistance
- A variety of Enamel and Dentin shades to achieve the "Camellion Affect" and reach high esthetic results:

ProFil A1, A2, A3, A3.5, B1, B2, B3, C2, C3, incisal (P), OA2, OA3, OA3.5, OB2

ProFil Flow A1, A2, A3, A3.5, B1, B2, C2, C3

Presentation

ProFil

- Refill syringes: 4gr & 2gr syringes
- Unidose Capsules: Box of 20x0.25g + 33% Free Material (for waste control)
- Working Kit:
 - 4x4g syringes including total or self etch bonding
- Starter Kit:
 - 8x4g syringes including total or self etch bonding

ProFil Flow

- Refill syringes: 1g, 3g for the price of 2g
- Economy Kit: 4x1g syringes
- **New!** Self Leveling Viscosity (SLV) Kit:
 - 3x1g High viscosity +1g Low viscosity

SILMET
www.silmetdental.com

For more information,
call our Support Center:
1-800-972-7917/866-300-3657

	Universal Composites			Flowable Composites		
	ProFil™ Silmet	Herculite® XRV™ Kerr	Filtek™ Z250 3M ESPE	ProFil™ Flow Silmet	Revolution™ Formula 2 Kerr	Filtek™ Flow 3M ESPE
Weight of total inorganic filler (%)	78	59	60	60	55	65
Volumetric Shrinkage (%)	2.4	2.9	2.2	3.4	4.7	3.1
Depth of cure (mm)	1.5-2	2-3	2-2.5	1.5-2	2	2
Compressive strength (MPa)	315	380	400	295	360	320
Flexural Strength (MPa)	120	150	150	125	120	120

ProFil™ 

ProFil™ Flow 


Clinical evidence available:



Dental Products Report
(03/2012, Vol. 46)



The Dental Advisor
(Oct. 2009, Vol. 26 No. 8)



Implementing Sandwich
Technique with RMGI,
By Dr. Lev Nudelman (Oct. 2010)

Presentation



Special Offers*

Free LED



2+1: Buy 2 get 1 free



*Check availability with local distributor