



Polaroid

Wired Intraoral Camera

1 Features

-  **Auto focus**
Real Autofocus function shows clear images.
-  **Gyro mouse built-in**
Useable as a mouse.
-  **280° twist head**
Easy to get approach anywhere in oral cavity.
-  **Adjustable LED brightness**
Can adjust LED brightness to avoid light reflecting.
-  **OLED screen**
Bright and clear OLED screen indicates various device status.
-  **Software compatibility**
Support keyboard short-cut data changing program and twain driver.

2 Specifications

Camera (included head)	Operation voltage	5V
	Measurement	Approx. 202 x 30 x 28 mm
	Weight	Approx. 56g
	Resolution	1.3 Mega pixel sensor
	Sensor/Lens type	CMOS / Liquid lens
USB cable	Length	Approx. 3 M
Disposable cover	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

3 Minimum System Requirements

PC system:	2.0 Ghz CPU, 1GB RAM
OS:	Windows XP service pack 2
USB port:	USB 2.0

4 Integration

Integrates with almost ALL software's; or we provide a free Imaging software.

