

MARK3[®]

Premium Products. Affordable Prices.

MARK3 UniCore Universal

Micro-Hybrid Composite Unidose | .315gm | 20/pk |

MARK3's UniCore Universal Micro-Hybrid Composite delivers a versatile esthetics for anterior and posterior restorations. Designed with high-density inorganic filler loading, UniCore helps minimize polymerization shrinkage while supporting excellent marginal sealing, polishability, strength, and wear resistance.

Features & Benefits

- For anterior & posterior restorations
- High polishability with great esthetics
- Does not stick to instruments
- Optimal balance of physical properties
- High strength and wear resistance
- Ultra-low shrinkage



Performance Snapshot

Flexural strength: > 152 MPa

Water sorption : 1.8 - 14 ug/mm³

Water solubility: 1.2 - 3 ug/mm³

Curing time: Regular shades: 20 secs

Opaque shades: 30

Intended Use

- Class I - V restorations on posteriors & anteriors
- Veneering of discolored anteriors
- Splinting of loose teeth
- Reconstruction of traumatically affected anteriors
- Anchoring of orthodontic appliances

COMPARES TO

Filtek Supreme, Dentsply Esthet-X, Kerr Herculite XRV Ultra

PACKAGING

Contains (20) .315gm Compules

Esthetic composite performance with strength, polish, and easy handling.

MARK3 UniCore Universal