

ZEN-PX2

Portable Dental X-ray System

- Excellent Image Quality
- Elegant Design
- Ideal Output for Film & Digital sensor : 60kV, 2mA
- Easy to use

Compact Outside, Smart Inside

New era of Intra-Oral X-ray system is now in your hand.
ZEN-PX2 will offer you ease and comfort in your practice.




New era of Intra-Oral X-ray system is now in your hand.
ZEN-PX2 will offer you ease and comfort in your practice.

- Light and compact unit lets you bring your practice to anywhere.
- High performance Tube and HFG makes sharp and clear images.
- Advanced Radiation Shielding ensures less than 0.13% of allowed doses leak.
- Ideal to use with Digital sensor and chemical Film.

Technical Specification & Dimension

X-ray Generator	Type	High Frequency Inverter
	Power Output	180W
X-ray Installation	Tube Voltage	60 kV
	Tube Current	2mA
	Heat Capacity	8.5 KHU
X-ray Tube	Anode Type	Stationary
	Anode Angle	20°
	Focal Spot	0.8mm
X-ray Control	Exposure Time Set	56 Steps (0.01 ~ 2.0 sec)
Power Requirement		22.2 VDC
Inherent Filtration	Al equivalent	1.8 mmAl
Battery Capacity	Number of exam	100 per full charge
Weight	Main Body	5.18lb (2.35kg)
Dimension	Main Body	7.75" x 5.78" x 5.70" (19.7cm x 14.7cm x 14.5cm)

* Please do not hesitate to contact us if you need detailed specification of our products.



Components of ZEN-PX2

Main Unit

- X-ray tube assembly,
Beam limiting device,
Battery, Grip tie

Accessories

- Power adapter, Neck strap,
Remote control, Manual

● Compact Design

Camera type ZEN-PX2 offers convenient diagnosis with its portable ability and compact design. ZEN-PX2 also offers large battery capacity capable of over 100 exams on single charge.

● Graphic LCD

Large LCD display lets you monitor the process and setting.

● Grip Tie

Adjust length of this comfortable tie to fit your hand & be completely mobile.

● Operation Panel

X-ray operation have never been easier with easy to understand push button & optimized pre-settings.