

SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

Trade name: MARK3 UniCore Universal Composites Syringe and Compules, MARK3 UniCore Flow

Manufacturer: Cargus International, Inc.

Address: 135 North Route 9W, Congers, NY 10920

Telephone: 845-267-2600 **email:** info@cargus.com

Product use:

Intended Use: Dental product

Limitations on use: For use by dental professionals

Specific Use: Composite restorative material

SECTION 2: HAZARDS IDENTIFICATION**2.1. Classification of the substance or mixture**

Classification according to Regulation (EC) No 1272/2008

<i>Skin Irrit. 2</i>	H315	Causes skin irritation
<i>Eye Irrit. 2</i>	H319	Cause serious eye irritation
<i>Skin Sens. 1</i>	H317	May cause an allergic skin reaction

2.2 Label elements

Labelling according to Regulation (EC) No 1272/2008

The product is classified and labelled according to the CLP regulation.

Hazard pictograms

GHS07

Signal word Warning

Hazard-determining components of labelling

Bisphenol A glycidylmethacrylate

Triethyleneglycol dimetacrylate

Hazard Statements

H315 Causes skin irritation

H319 Causes serious eye irritation

H317 May cause an allergic skin reaction

H335 May cause respiratory irritation

Precautionary statements

- P280 Wear protective gloves/clothing/eye protection/face protection
P305 + P351+P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lens if present and easy removing clothes.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention
P501 Dispose of contents/container in accordance with national/international regulations.

Additional Information

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body, are exempted from the provisions of regulation (EC) No 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

Classification according to Regulation (EC) No 1272/2008: exempted

2.3 Other Hazards

Results of PBT and vPvB assessment

N/A

SECTION 3: INGREDIENTS

Chemical characterization: Mixture

Description: Mixture of substances listed below with nonhazardous additions

- CAS: 1565-94-2 Bisphenol A glycidyl dimethacrylate
CAS: 109-16-0 Triethyleneglycol dimetacrylate

SECTION 4: FIRST AID MEASURES

The following first aid recommendations are based on assumption that appropriate personal and industrial hygiene practices are followed.

Eye Contact: Flush eyes with large amounts of water. If signs/symptoms persist, get medical attention.

Skin contact: Note: Uncured material contact to skin directly may cause irritation or skin rash especially in sensitive individuals

Remove contaminated clothing and shoes. Immediately flush skin with large amounts of water. Get medical attention. Wash contaminated clothing and clean shoes before reuse.

Inhalation: Remove person to fresh air. If signs/symptoms persist, get medical attention.

If swallowed: Do not induce vomiting unless instructed to do so by medical personnel.

Give victim water. Never give anything by mouth to unconscious person. Get medical attention.



MARK 3	
Technical File	Composites
Safety Data Sheet	Rev.0
Title: SDS-Dental Composites	Page: 4 of 6
	Date: 20.10.2025

SECTION 5: FIRE FIGHTING MEASURES

5.1 Flammable properties

Auto ignition temperature Not data available

Flash point Not applicable

Flammable Limits – LEL Not applicable

Flammable Limits – UEL Not applicable

5.2 Extinguishing media

Water fog, CO₂, foam, dry chemical

Unusual Fire and Explosion Hazards: Upon combustion irritating gases may form.

SECTION 6: ACCIDENTAL RELEASE MEASURES

Accidental Release Measures:

Evacuate unprotected and untrained personnel from hazard area. The spill should be cleaned up by qualified personnel. Ventilate the area with fresh air.

Wipe up small spills with paper towels or dry cloth. Clean remainder with alcohol.

In the event of a release of this material, the user should determine if the release qualifies as reportable according to local, state and federal regulation.

SECTION 7: HANDLING AND STORAGE

7.1 Handling

Only by qualified staff in dental offices or laboratories

7.2 Storage

Avoid exposure to intense light and elevated heat, store in room temperature..

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Engineering controls

Use in a well-ventilated area.

8.2 Personal Protective Equipment (PPE)

8.2.1 Eye/Face Protection

Avoid eye contact.

The following eye protection(s) are recommended: Safety Glasses with side shields.

8.2.2 Skin Protection

Avoid skin contact. Wear lab coat.



MARK 3	
Technical File	Composites
Safety Data Sheet	Rev.0
Title: SDS- Dental Composites	Page: 5 of 6
	Date: 20.10.2025

8.2.3 Respiratory Protection

N/A

8.2.4 Prevention of Swallowing

Do not ingest. Wash hands after handling and before eating.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Specific Physical Form:	Paste
Odor, Color, Grade:	Slight acrylate odor, tooth-colored various shade
General Physical Form:	Solid
Autoignition temperature	No Data Available
Flash Point	Not Applicable
Flammable Limits - LEL	Not Applicable
Flammable Limits - UEL	Not Applicable
Boiling point	Not Applicable
Specific Gravity	>1
pH	N/A
Solubility in Water	Insoluble

SECTION 10: STABILITY AND REACTIVITY

Stability: Unstable > 40°C

Materials and Conditions to Avoid:

10.1 Conditions to avoid: Heat and light

10.2 Incompatible Materials: Peroxide and Amines

Hazardous Polymerization: Heat and light.

Hazardous Decomposition: None

SECTION 11: TOXICOLOGICAL INFORMATION

No data available

SECTION 12: ECOLOGICAL INFORMATION

ECOTOXICOLOGICAL INFORMATION: Not determined.

CHEMICAL FATE INFORMATION: Not determined.

SECTION 13: DISPOSAL CONSIDERATIONS

Waste Disposal Method:

EPA Hazardous Waste Number (RCRA): Not regulated

Dispose of contents/ container in accordance with the local/regional/national/international regulations.

Since regulations vary, consult applicable regulations or authorities before disposal.

SECTION 14: TRANSPORTATION INFORMATION

Land transport ADR/RID

Not classified as dangerous goods in the meaning of transport regulations

Sea transport IMDG

Not classified as dangerous goods in the meaning of transport regulations

Air transport ICAO-TI and IATA-DGR

Not classified as dangerous goods in the meaning of transport regulations

MARK 3[®]	
Technical File	Comosites
Safety Data Sheet	Rev.0
Title: SDS- Dental Composites	Page: 6 of 6
	Date: 20.10.2025

SECTION 15: Regulatory information

Classification according to EC-guidelines:

The product is a medical device according to EC-directive 93/42 EEC.

This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada

US FEDERAL REGULATIONS

This MSDS has been prepared to meet the U.S. OSHA Hazard Communication Standard, 29 CFR 1910.1200.

SECTION 16: OTHER INFORMATION

CAUTION: PRODUCT FOR PROFESSIONAL USE

The information on this safety sheet is based on presently available data and to our best knowledge for the correct handling of the product under normal conditions. Any use of this product in any way not indicated on this sheet or use of this product together with any other process/procedure will be exclusively under the user's responsibility.