

# MATERIAL AND SAFETY DATA SHEET

## Powder Free Latex Examination Gloves

### Supermax Healthcare Incorporated,

1899 Sequoia Dr.

Aurora, Illinois

60506 USA

Tel : 00-1-630-898-8886

Fax : 00-1-630-898-8855

Toll Free : 00-1-877-287-3542

### SECTION 1 COMPOSITION

These gloves are made of 100% Natural Rubber Latex providing strength, comfort, superior fit and grip due to its elastic properties.

### SECTION 2 PHYSICAL DATA

Meet with the requirements of ASTM D 3578 : 2000 Standard Specification for Rubber Examination Gloves and BS EN 455-1 , 2 & 3 : 2000 Medical Gloves For Single Use.

Width	95 ± 10 mm ( size medium)
Length	240 mm minimum
Thickness	0.08 mm minimum (providing tactile sensitivity)
Tensile Strength (unaged)	14 MPa minimum (providing superior strength)
Ultimate elongation (unaged)	700% minimum
Tensile Strength (aged)	14 MPa minimum
Ultimate elongation (unaged)	500% minimum
Watertightness	Substantially impermeable to water vapour and liquid water providing an excellent biological barrier. Double gloving is recommended for reduced risk.

### SECTION 3 HEALTH HAZARD INFORMATION

Biocompatibility data	Guinea Pig Sensitization (Buehler) - Did not indicate a potential for dermal irritation or allergic contact sensitization.  Repeated Insult Patch Test - Did not indicate a potential for dermal irritation or allergic contact sensitization.
-----------------------	--

