

maxill

TUBERCULOCIDAL • VIRUCIDAL[†] • BACTERICIDAL • FUNGICIDAL

mycolio™

DISINFECTANT WIPES

mycolio DISINFECTANT WIPES have a unique formulation of four active ingredients that work synergistically to create an unparalleled broad spectrum disinfectant ideal for daily use in dental offices. mycolio DISINFECTANT WIPES can be used on a multitude of hard surfaces to inactivate a broad spectrum of pathogenic microorganisms.

- Broad Spectrum Bactericide, Virucide, Tuberculocide and Fungicide
- Meets CDC guidelines for intermediate level surface disinfection
- Cleaner and disinfectant in one product; save money, time and storage space
- Compatible with a wide range of hard, non-porous surfaces
- Low lint, heavy duty, non-woven, pre-saturated wipes, 6" x 7" and 7.5" x 9.75"
- Packaged in recyclable HDPE tubs with self-closing dispenser lids
- Convenient refill rolls are available to reduce waste and save you money

1-minute contact time

Human Immunodeficiency virus type 1 (HIV-1)*

2-minute contact time

Vancomycin Resistant *Enterococcus faecalis* (VRE)
Methicillin Resistant *Staphylococcus aureus* (MRSA)
Duck Hepatitis B Virus (Surrogate for Human Hepatitis B Virus)*
Bovine Viral Diarrhea Virus (Surrogate for Human Hepatitis C Virus)*
Avian Influenza A (H3N2) Reassortant virus
Salmonella enterica

3-minute contact time

Mycobacterium bovis BCG
Pseudomonas aeruginosa
Staphylococcus aureus
Poliovirus Type 1
Trichophyton mentagrophytes



Cleaning Care for Dental Chairs

IN ADDITION TO DISINFECTING, A CLEANING PROTOCOL IS REQUIRED
Call 1-855-462-9455 for more information
www.maxill.com/us/chairecare



1, 2, 3 Minute
Contact Time

EPA Reg. No. 88897-1-90841

Designed with **You** in mind!

Product Availability

160 wipes tub	#61161
160 wipes refill roll	#61191
110 jumbo wipes tub	#61181
110 jumbo wipes refill roll	#61196

Visit us online for more information,
including SDSs and pricing:

www.maxill.com/mycolio

Made in Canada

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-61161BU
Rev.3: 2019-11-22

500 West Main Street, Cortland, OH 44410
www.maxill.com 1-855-462-9455