

Chroma Chameleon

Universal One-Shade Composite Restorative

Chroma Chameleon One Shade Light Cure Composite from MARK3 is a shade matching universal composite that easily blends and adapts to surrounding tooth structure. Used in anterior and posterior restorations, Chroma Chameleon matches all vita shades from A1-D3. Besides its permanent gloss and toothlike opalescence, Chroma Chameleon shows excellent physical properties and superior handling. Contains innovative fillers with optimal light diffraction technology. In rare cases, a blocker may be needed to mask tooth discoloration or stains. For class I-V restorations.

Features & Benefits

- **Advanced Shade Matching from A1-D3**
- Universal Composite for Anterior & Posterior Restorations
- Contains Innovative Spherical Fillers with Optimal Light Diffraction Technology
- Superior Polishability
- Permanent Gloss and Toothlike Opalescence
- Reduces chairtime and inventory for dentists.
- Radiopacity 200%



Performance Snapshot

Opacity: 65%
Flexural Strength: 121 MPa
Compressive Strength: 294 MPa
Modules of Elasticity: 9.5 Gpa

CHROMA BLOCK



Indications for Use (Block)

- Masking of the oral cavity in case of extended class III and IV cavities
- Masking of tooth discoloration.
- SKU #100-4850

COMPARES TO: Omnicroma from Tokuyama, Uni-Tone from Hexa, Transcend™ from Ultradent

PACKAGING: Contains (1) 4gm Syringe | SKU #100-4852