



PRODUCTS

See more innovative Silmet solutions on DentalProductShopper.com.



- High radiopacity
- "Flow-On-Demand," handling preferred by dentists
- Expanded shade range and stable color match
- Supreme polish retention

SILMET

ProFil Flow

ProFil Flow is a visible-light activated, radiopaque, flowable composite from Silmet. It features excellent flowable consistency characteristics because of unique nanotechnology additions, making it a versatile composite that can handle a variety of indications.

ProFil Flow offers high compressive and tensile strength for optimal wear resistance based on the advanced resin technology of ProFil.

ProFil Flow is indicated for use in Class 3, Class 5, smaller Class 4's, as well as a base/liner in Class 1 and Class 2 restorations. It also can be used to repair resin, porcelain, and acrylic

temporary materials. In addition, ProFil Flow is indicated for use as a pit and fissure sealant, an undercut blockout, and for restoration of minimally invasive cavity preparations.

ProFil Flow is available in 1-g and 3-g syringes in shades A1, A2, A3, A3.5, B1, B2, B3, and C2. The economy kit contains four 1-g syringes in a variety of shades.



SILMET LTD
 FREE Product Information
 CALL 800.520.1408 ext. 94209
 CLICK www.dps.li/a/46J-209
 CIRCLE Reader Service 209