

INSIGHT



## Minimize exposure. Maximize results.

The standard for intraoral film. With INSIGHT dental film, you'll get a true F-speed film that requires up to 60% less radiation per exposure.

- Delivers excellent and consistent image quality through exclusive T-Grain emulsion technology.
- Lowers patient and staff radiation exposure without sacrificing image resolution.
- Needs no continual equipment adjustments to accommodate processing conditions.
- Holds speed and contrast in variable chemistry conditions.
- Is available in all intraoral film sizes and formats, including ClinAsept barrier packets for added safety.

## Exposure Guidelines

| General D-speed film users |                           | General E-speed film users |                           |
|----------------------------|---------------------------|----------------------------|---------------------------|
| Automatic processing       | Reduce exposure time 60%* | Automatic processing       | Reduce exposure time 20%* |
| Manual processing          | Reduce exposure time 50%* | Manual processing          | Reduce exposure time 20%* |

\*Reduce exposure time only. Do not change mA or kVp settings. To fine tune image – check exposure recommendations inside box.

For best results, use with READYMATIC, READYPRO, or GBX processing chemicals. INSIGHT dental film processes as an F-speed film when processed in roller transport automatic processors; on the high-end of E-speed range when processed manually.

## Ordering Information for INSIGHT

| Size   | Description  | Code   | Films/Pkgs | REF No   |
|--------|--|--------|------------|----------|
| Size 0 | Super Poly-Soft packets                                      | IP-01  | 100        | 867 5332 |
| Size 0 | Super Poly-Soft packets (double film)                        | IP-02  | 100        | 128 0619 |
| Size 0 | Super Poly-Soft packets with ClinAsept barrier               | IP-01C | 75         | 120 0328 |
| Size 0 | Bite-wing paper packets                                      | IB-01  | 50         | 811 0785 |
| Size 1 | Paper packets  | IP-11  | 100        | 112 4981 |
| Size 1 | Paper packets (double film)                                  | IP-12  | 100        | 805 5402 |
| Size 2 | Super Poly-Soft packets                                      | IP-21  | 150        | 116 3401 |
| Size 2 | Super Poly-Soft packets (double film)                        | IP-22  | 130        | 179 8628 |
| Size 2 | Super Poly-Soft packets with ClinAsept barrier               | IP-21C | 100        | 153 9931 |
| Size 2 | Super Poly-Soft packets with ClinAsept barrier (double film) | IP-22C | 100        | 156 0390 |
| Size 2 | Paper packets  | IP-21  | 150        | 176 3960 |
| Size 2 | Paper packets (double film)                                  | IP-22  | 150        | 107 9086 |
| Size 2 | Bite-wing paper packets                                      | IB-21  | 50         | 180 7650 |
| Size 3 | Bite-wing paper packets                                      | IB-31  | 100        | 829 8929 |
| Size 4 | Occlusal paper packets                                       | IO-41  | 25         | 116 9143 |

For more information, call **800.933.8031** or visit **[www.carestreamdental.com](http://www.carestreamdental.com)**.