

The details
you need.



The comfort your
patients want.

Soft edges and high-contrast make
AGFA the patient's choice for comfort
and the dentist's choice for diagnosis.

- AGFA Dentus® intraoral film is soft, flexible, and features a radius edge to ensure patient comfort
- High contrast and maximum density with excellent grey graduation
- Intraoral films are fully protected from moisture contamination and can be disinfected



AGFA
Patient Preferred.

AGFA Dentus®

Intraoral Film

PATIENT PREFERRED

- Soft, flexible film features a radiused edge to ensure patient comfort
- Exclusively made with Softopac®, smooth edges and rounded corners eliminate irritation

OPTIMAL IMAGE SHARPNESS

- High contrast and high maximum density support the diagnosis due to excellent grey graduation
- Flexible but firm films help to reduce distortion

CONVENIENT

- Double sealed packets are fully protected from moisture contamination
- Can be immersed in disinfection liquid without worry
- Full offering to meet all your needs—adult PA (size 2), pedo (size 0), intraoral E/F and D Speed, extraoral, duplicating, and accessories

Extraoral Film

- Very high definition due to disc grain technique
- Very low Base-Plus-Fog Density
- Particular contrasting in various density ranges

AGFA Dentus M2 Comfort



Speed F/E 0.10 sec.



Speed D 0.20 sec.

Competitor Film



Speed F/E 0.10 sec.



Speed D 0.20 sec.

Prepared under identical condition to yield optimal density.
X-ray equipment: Shortcon, 70kV–7mA, head phantom – mandibular, filter 10 mm.
Automatic processor: 28°C, 6 minutes.



M2 COMFORT — INTRAORAL E SPEED	M2 COMFORT — INTRAORAL D SPEED	PANORAMIC
<p><i>Less Radiation! Shorter exposure time with E/F-Speed film</i></p> <p>SIZE 0</p> <p>65031124 M2-54 1 Film Per Packet (100 Count)</p> <p>65031256 M2-53 2 Films Per Packet (100 Count)</p> <p>SIZE 2</p> <p>65414615 M2-58 1 Film Per Packet (150 Count) Case Lot 24</p> <p>65414623 M2-57 2 Films Per Packet (150 Count)</p> <p>SIZE 4 (OCCLUSAL)</p> <p>65870955 M2-50 1 Film Per Packet (25 Count)</p>	<p><i>Traditional exposure time. No change to current technique or settings</i></p> <p>66009814 M2-D01 AGFA Dentus M2 Comfort D (2 x 3) 100 x 1 Sz-0</p> <p>66009815 M2-D02 AGFA Dentus M2 Comfort D (2 x 3) 100 x 2 Sz-0</p> <p>66009817 M2-D1 AGFA Dentus M2 Comfort D (3 x 4) 150 x 1 Sz-2 (Case Lot 24)</p> <p>66009818 M2-D2 AGFA Dentus M2 Comfort D (3 x 4) 150 x 2 Sz-2 (Case Lot 24)</p>	<p>66000311 Ortholux (5 x 12 in) 100 Films, Green Sensitive</p> <p>66000312 Ortholux (6 x 12 in) (15 x 30 cm) 100 Films, Green Sensitive</p>
AGFA DUPLICATING FILM	CEPHALOMETRIC	AGFA INTENSIFYING SCREENS (DENTUS GREEN)
<p>65038641 AGFA Dentus Duplicating Film (5 x 12 in, 100 Films)</p> <p>65038668 AGFA Dentus Duplicating Film (6 x 12 in, 100 Films)</p>	<p>66000314 Ortholux (8 x 10 in) 100 Films, Green Sensitive</p>	<p>65769847 Ortho Regular (15 x 30 cm) (6 x 12 in)</p>
		ACCESSORIES
		<p>66002011 Agfa Dentus Ortholux Cassette (15 x 30 cm)</p> <p>66002014 Agfa Dentus Ortholux Cassette (8 x 10 in)</p>

To learn more about AFGA call 1-877-258-5937.