

Section 1 – Identification

Product Name: Ionite-H Home Care 0.40% Stannous Fluoride Brush-on Gel
Active Ingredient: Stannous Fluoride – Tin (II) Fluoride
Manufacturer: Dharma Research, Inc. 5220 N.W. 72nd Avenue Miami, FL 33166
Information Contacts: (305) 482-9669
Toll Free: 1-877-833-3725
Emergency Phone Numbers: US & Canada 1 (877) 833-3725

Family: Fluoride Treatment
Product Use: Professional Dental Preventive Care
Product #: 56-00130, 56-00131

Section 2 – Hazards Identification

EMERGENCY OVERVIEW

This information is based on findings from related or similar materials.

- Will cause eye irritation.
- May cause skin irritation.
- Prolong exposure may cause respiratory system irritation.
- Ingestion of large amounts may cause nausea, vomiting, headaches, or dizziness.
- Please read entire MSDS for additional information.



Potential Health Effects, Signs and Symptoms of Exposure:

Primary Route of Entry Inhalation, skin, and ingestion.

Eye Vapors are irritating to the eyes. Splashes may cause redness of the eye.

Skin Repeated/prolonged contact may cause mild irritation of the skin. Symptoms include redness, burning, or itching.

Ingestion Swallowing small amounts during normal handling is not likely to cause harmful effects; swallowing large amounts may be harmful. Symptoms include irritation of the throat, nausea, dizziness, and upset stomach.

Inhalation If fine small particles are inhaled, irritation of the throat with a feeling of tightness in the chest may occur.

NOTE: Refer to Section 11, Toxicological Information for Details

Section 3 – Composition/Information on Ingredients

Chemical Identity	CAS Numbers	EINECS#	INCI Name	Exposure OSHA	Limits ACGIH	Carcinogen IARC/NTP/OSHA	%
				TWA/STEL	TWA/STEL		
Glycerol	56-81-5	200-289-5	Glycerin	(as respirable, mist) TWA = 5 mg/ m ³	(as mist) TWA = 5 mg/ m ³	Not listed	66 – 77%
Stannous Fluoride	7783-47-3	231-999-3	Stannous Fluoride	TWA = 2.5 mg/ m ³	TWA = 2.5 mg/ m ³	Not Listed	0-2%

N/E – None Established
 N/R – Not Reviewed
 N/DA – No Data Available
 N/A – Not Applicable

(items in parenthesis relate to 1999/45/EC)

Stannous Fluoride: Danger Symbol – GHS05 (C), GHS07 (Xn) Hazard Statement – H302 (R22), H313 (N/A), H315 (R38), H318 (R41)
 Precautionary Statements – P280 (S36/37/39), P305+351+338 (S26)

Glycerol: Danger Symbol – GHS07 (Xi) Hazard Statement – H315 (R38), H319 (R36) Precautionary Statement – P404 (S7), P403+ 235 (S3/9)

See Section 16 for Hazard and Precautionary Statement Key.

Section 4 – First Aid Measures

First Aid for Eye Flush with water for 15 minutes, including under eyelids. Get medical help if discomfort persists.

First Aid for Skin Wash thoroughly with soap and water. Remove contaminated clothing. Get medical help if discomfort persists.

First Aid for Ingestion If greater than normal dose is swallowed, do not induce vomiting. Drink large quantity of water or milk. Seek medical attention.

First Aid for Inhalation Move to fresh air. If having breathing difficulty, give oxygen. Seek medical attention if discomfort persists.

Section 5 – Fire Fighting Measures

Flash Point (°F/°C)	Flammable Limit (vol%)	Auto-ignition Temperature (vol%)
N/A	N/A	N/A

Method:

Extinguishing Media: Suitable extinguishing media for the surrounding fire should be used.

Fire Fighting Instructions: During a fire, irritating or highly toxic gases may be generated by thermal decomposition. If potential for exposure, wear complete personal protective equipment including self contained breathing apparatus, with full face operated in pressure demand. Fight fire from a safe distance/protected location.

Unusual Hazards: Evacuate all personnel.

Section 6 – Accidental Release Measures

Spill or Release Procedures Minor spills – Clean up immediately, avoid contact with skin and eyes. Mop area and clean with soapy water

Major spills – Clear area of personnel. Restrict access to area. Avoid contact with skin and eyes, and mop / squeegee area with soapy water. Collect contaminated water and dispose of properly.

See section 8 & section 12.

Section 7 – Handling and Storage

Handling Limit all unnecessary personal contact. Avoid breathing dust / vapors. Wear appropriate PPE.

Storage Store in a cool, well ventilated area away from heat, sparks and flame. Keep containers closed when not in use. Keep in original container provided by manufacturer.

Explosion Hazard None.

Section 8 – Exposure Controls / Personal Protection

Engineering Controls General exhaust is adequate under normal operating conditions.

Personal Protective Equipment

General To identify additional Personal Protective Equipment (PPE) requirements, it is recommended that a hazard assessment in accordance with the OSHA PPE Standard (29CFR1910.132), or European Standard EN166 be conducted before using this product.

Eye/ Face Protection Safety glasses / goggles or splash shields are required when handling. Ensure eye bath is on hand.

Skin Protection No special equipment for minor exposure (i.e. when handling small quantities), but is recommended to use impermeable clothing to prevent ANY contact with this product, such as gloves, apron, boots, or whole body suit.

Respiratory Protection In case of insufficient ventilation, wear suitable respiratory equipment.

Section 9 – Physical and Chemical Properties

Appearance	Odor & Odor Threshold	pH	Specific Gravity	Viscosity	% Volatile
------------	-----------------------	----	------------------	-----------	------------

Safety Data Sheet

0.4% STANNOUS FLUORIDE GEL

Colored Gel	Odor of Stated Flavor	N/A	(H2O = 1): N/A	N/A	N/A		
Boiling Point/ Freezing Point	Decomposition Temperature	Octanol/Water Partitioning Coefficient Log Po/w	Vapor Pressure:	Vapor Density	Evaporation Rate	Ignition	Solubility In Water (20°C)
N/A	N/DA	N/DA	N/A	N/A	N/A	N/A	N/A
Flash Point (°F/°C)		Flammable Limit (vol%)		Auto-ignition Temperature (vol%)			
N/A		N/A		N/A			

Section 10 – Stability and Reactivity

Stability: Stable	Incompatibility (Materials to Avoid): Strong Acids
Hazardous Decomposition Products: Toxic Fumes (HF)	Hazardous Polymerization: Will not occur
Conditions to Avoid: Heat and incompatible materials	

Section 11 – Toxicological Information

Acute Oral Toxicity N/DA	Acute Dermal Toxicity N/DA	Acute Inhalation Toxicity N/DA	Irritation – skin N/DA	Irritation – Eye N/DA
Sensitization N/DA		Mutagenicity N/DA		Sub-chronic Toxicity N/DA

Section 12 – Ecological Information

Ecotoxicological Information

Acute Toxicity to Fish N/DA	Acute Toxicity to Invertebrates N/DA	Acute Toxicity to Algae N/DA	Bioconcentration N/DA	Toxicity to Sewage Bacteria N/DA
---	--	--	---------------------------------	--

Chemical Fate Information

Biodegradability	N/DA. This material is not expected to significantly bioaccumulate.
Chemical Oxygen Demand	N/DA

Section 13 – Disposal Considerations

Dispose of in compliance with governmental regulation (state and federal).

Dispose of container and unused contents in accordance with federal, state and local requirements. For EU Member States, please refer to any relevant Community provisions relating to waste. In their absence, it is useful to remind the user that national or regional provisions may be in force.

Section 14 – Transport Information

DOT (49 CFR 172)	
Proper Shipping Name:	Non-regulated material
Identification Number:	N/A
Marine Pollutant:	No
Special Provisions:	N/A
Emergency Response Guidebook (ERG) #:	N/A
IATA (DGR):	Non-regulated material
Proper Shipping Name:	N/A
Class or Division:	N/A
UN or ID Number:	N/A

Safety Data Sheet

0.4% STANNOUS FLUORIDE GEL

Packaging Instructions:	N/A
Emergency Response Guidance (ICAO)#:	
IMO (IMDG):	
Proper Shipping Name:	Non-regulated material
Class or Division:	N/A
UN or ID Number:	N/A
Special Provisions & Stowage/Segregation:	N/A
Emergency Schedule (EmS)#:	N/A
Other Information:	N/A

Section 15 – Regulatory Information

US Federal Regulations

Clean Air Act: HAP/ODS	This product contains the following HAP's or ODS: <ul style="list-style-type: none"> NONE
Clean Water Act: Priority Pollutant	This product contains the following chemicals listed under the U. S. Clean Water Act Priority Pollutant and Hazardous Substance List: <ul style="list-style-type: none"> None
FDA: Food Packaging Status	This product has not been cleared by the FDA for use in food packaging and / or other applications as an indirect food-packaging additive.
Occupational Safety and Health Act	This product is considered to be hazardous under the OSHA Hazard Communication Standard. It's hazards are: <ul style="list-style-type: none"> Stannous Fluoride – CAS #7783-47-3 (Toxic, Irritant)
RCRA	This product contains the following chemicals considered to be hazardous waste under RCRA (40 CFR 261). <ul style="list-style-type: none"> None
SARA Title III: Section 302 (RQ)	This product contains no chemicals regulated under Section 302 as extremely hazardous substances.
SARA Title III: Section 302 (TPQ)	This product contains chemicals regulated under Section 304 as extremely hazardous chemicals for emergency release notification ("CERCLA" List): <ul style="list-style-type: none"> None
SARA Title III: Section 311-312:	This product is considered to be hazardous under the OSHA Hazard Communication Standard and is regulated under Section 311-312 (40 CFR 370). It's hazards are: <ul style="list-style-type: none"> Stannous Fluoride – CAS #7783-47-3 (Toxic, Irritant)
SARA Title III: Section 313:	This product contains the following chemicals which are subject to the reporting requirements of Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 and 40 CFR Part 372: <ul style="list-style-type: none"> None
TSCA Section 8(b): Inventory:	This product does contain chemicals listed on the TSCA inventory or otherwise complies with TSCA pre-manufacture notification requirements. <ul style="list-style-type: none"> Stannous Fluoride – CAS# 7783-47-3 Glycerol – CAS # 56-81-5
TSCA Significant New Use Rule:	None of the chemicals in this material have a SNUR under TSCA.

State Regulations

CA Right-to-Know Law: California No Significant Risk Rule:	None
MA Right-to-Know Law:	Glycerol – CAS # 56-81-5
NJ Right-to-Know Law:	Stannous Fluoride – CAS #7783-47-3 Glycerol – CAS # 56-81-5

Safety Data Sheet


0.4% STANNOUS FLUORIDE GEL

PA Right-to-Know Law:	Stannous Fluoride – CAS #7783-47-3 Glycerol – CAS # 56-81-5
FL Right-to-Know Law:	N/A
MN Right-to-Know Law:	Glycerol – CAS # 56-81-5

International Regulations

CDSL: Canadian Inventory (on Canadian Transitional List)	Stannous Fluoride – CAS #7783-47-3 Glycerol – CAS # 56-81-5
---	--

Labeling according to EC directives – 1272/2008 {CLP} AND 1999/45/EC (items in parenthesis relate to 1999/45/EC)

<p>European Community:</p> 	<p>For 0.4% Stannous Fluoride Gel (finished product):</p> <ul style="list-style-type: none"> DANGER SYMBOLS: GHS07 (Xi) – Warning / Irritant HAZARD STATEMENT: H316 (N/A), Causes mild skin irritation. H319 (R36), Causes serious eye irritation H335 (R37), May cause respiratory irritation. H303 (N/A), May be harmful if swallowed. PRECAUTIONARY STATEMENT: P102 (S2), Keep out of reach of children. P261 (S24), Avoid breathing dust/fume/gas mist/vapor sprays. P280 (S36/37/39), Wear protective gloves/clothing/eye protection/face protection. P309+314 (S45), If exposed or you feel unwell, get medical advice/attention.
--	--

Section 16 – Other Information

EU Classes and Risk / Safety Phrases for Referenced Ingredients (See Section 2):



(items in parenthesis relate to 1999/45/EC)

Danger Symbols:
GHS06 (T) – Toxic. GHS07 (Xi) – Warning / Irritants. GHS05 (C) – Corrosive. GHS07 (Xn) - Harmful

Hazard Statement: **H301 (R25)**, Toxic if swallowed. **H302 (R22)**, Harmful if swallowed. **H313 (N/A)**, Maybe harmful in contact with skin. **H315 (R38)**, Causes skin irritation. **H316 (N/A)** Causes mild skin irritation. **H319 (R36)**, Causes serious eye irritation. **H335 (R37)**, May cause respiratory irritation. **H336 (R67)**, May cause drowsiness or dizziness. **EUH032 (R32)**, Contact with acid liberates very toxic gas. **H303 (N/A)**, May be harmful if swallowed. **H318 (R41)**, Causes serious eye damage.

Precautionary Statement:
P102 (S2), Keep out of reach of children. **P260 (S23)**, Do not breathe dust/fume/gas/mist/vapors/spray. **P261 (S24)**, Avoid breathing dust/fume/gas mist/vapor sprays. **P280 (S36/37/39)**, Wear protective gloves/clothing/eye protection/face protection. **P309+314 (S45)**, If exposed or you feel unwell, get medical advice/attention. **P404 (S7)**, Store in a closed container. **P403+235 (S3/9)**, Store in a well ventilated place. Keep cool. **P305+351+338 (S26)**, If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do so. Continue rinsing.

Hazard Rating System (Pictograms)

<p>NFPA:</p> 	<p>HMIS:</p> 
---	---

MSDS Prepared by:	Ricardo Carles
Revision History:	09/06/15 Revision A

The information presented herein was obtained from sources considered to be reliable. However, this information is provided without any warranty, expressed or implied, regarding its correctness or suitability for consumers intended use and/or application. For this and other reasons, we assume no responsibility and expressly disclaim liability for loss, damage or expense arising out of any way connected with the handling, storage, use or disposal of the product. This SDS was prepared expressly for this product. Use the materials only as directed. If the product is used as a component of another product, the information contained within the SDS may not be applicable.

Dharma components are provided on an as is basis without warranties of any kind either expressed or implied. Dharma does not warrant the use or the results of use of the materials sold on an as is basis since they are intended for remanufacturing or repackaging. It is the sole responsibility of the user to examine and determine appropriate application and regulatory requirements associated with said Dharma components.