

Report No: BTK20200212002FDA		DA Date: Mar. 13, 2020	Page 1 of 4	
Applicant	:	Homar Bio-Technology (Guangzhou) Holding Co., Ltd.		
Address	:	No.28 Qiyuan Road, Huashan Town, Huadu District, Guangzhou 510880, Guangdong, P.R. China		

The following sample(s) was /were submitted and identified on behalf of the clients as:

Sample Name	:	Particulate Respirator KN95 Face Mask	
Trade Name	:	N/A	
Sample Model	:	KN95-PTMYKZ-01	
Sample Received Date	:	2020.02.12	
Testing Period	:	2020.02.12 to 2020.03.13	
Test Requested	:	As specified by client, for compliance with Food and Drug Administration Regulations for determining the amount of total extractives from polypropylene articles intended for repeated use.	
General test result	:	Approved	
Test Method	:	Please refer to next page(s).	
Test Result	:	Please refer to next page(s).	



Signed for and on behalf of



Dexter Chan/ Approved Signatory

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of BTK, this test report shall not be copied except in full and published as advertisement. BTK Physical & Chemical Lab.

Website: http://www.btktest.com.cn



Report No: BTK20200212002FDA

Date: Mar. 13, 2020

Page 2 of 4

Sample Description:

No.	Description	Name
1	Textile	Particulate Respirator KN95 Face Mask

Test Result: No.1

Test Method: As specified in FDA 21 CFR 177.1040

Test item	Simulant Used	Result	Limit
Acrylonitrile	Distilled Water	15ppm	80ppm

Test item	Simulant Used	Result	Limit
Nitrogen content	Distilled Water	18.1%	17.4%-19%

Simulant Used	Test Condition	Total nonvolatile extractives	Limit
deionized water	2 hours at reflux temperature	<0.01 mg/in ²	0.01 mg/in ²
8 % alcohol	2 hours at reflux temperature	<0.01 mg/in ²	0.01 mg/in ²
50 % alcohol	2 hours at reflux temperature	<0.01 mg/in ²	0.01 mg/in ²
n-heptane	2 hours at reflux temperature	<0.01 mg/in ²	0.01 mg/in ²

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of BTK, this test report shall not be copied except in full and published as advertisement. BTK Physical & Chemical Lab.

Guangzhou Bestek Testing Services Co., Ltd

Website: http://www.btktest.com.cn



Report No: BTK20200212002FDA

Date: Mar. 13, 2020

Page 3 of 4

Sample Photos:



General view for KN95-PTMYKZ-01



General view for KN95-PTMYKZ-01

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of BTK, this test report shall not be copied except in full and published as advertisement. BTK Physical & Chemical Lab.

Guangzhou Bestek Testing Services Co., Ltd

Website: <u>http://www.btktest.com.cn</u>



Report No: BTK20200212002FDA

Date: Mar. 13, 2020

Page 4 of 4



General view for KN95-PTMYKZ-01



General view for KN95-PTMYKZ-01

*** End of Report ***

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of BTK, this test report shall not be copied except in full and published as advertisement. BTK Physical & Chemical Lab.

Guangzhou Bestek Testing Services Co., Ltd

Website: http://www.btktest.com.cn