

MARK3[®]

Premium Products. Affordable Prices.

MARK3 UniCore Flow

Flowable Composite Syringe | 3gm | 1/pk

MARK3 UniCore Flow is a light-cured, radiopaque flowable composite designed for excellent handling, fast placement, and strong esthetic results. Its non-sticky formulation offers smooth flow without slumping, while high polish retention, stable color match, and advanced fluoride release make it ideal for efficient restorative procedures.

Features & Benefits

- High polish retention & stable color match
- No oozing or slumping
- Non-sticky formulation
- Low shrinkage with excellent esthetics
- Light-cured and radiopaque
- Fast and easy manipulation
- Advanced fluoride release

Performance Snapshot

Sensitivity to ambient light: Above 100 sec.
Weight of total inorganic filler: 78%
Volume of total inorganic filler: 59%
Depth of cure: 1.5 mm
Compressive strength: 230 MPa
Flexural Strength: 145 MPa
Water Absorption: < 40 µg / ㎖
Solubility: < 7.5 µg / ㎖



Intended Use

- Class III, Class V & small Class IV restorations
- Base/liner in Class I & II restorations
- Veneering of discolored anterior teeth
- Splinting, reconstruction & orthodontic anchorage

COMPARES TO Revolution by Kerr, Clearfil Majesty, Beautifil Flow, Venus Diamond Flow, Tetric EvoFlow

PACKAGING Contains (1) 3gm Syringe & (5) Applicator Tips

Smooth flowable handling with strong esthetic performance.

MARK3 UniCore Flow