

Proven solutions for every cementation need.

Resin Cements



RelyX™ Ultimate Adhesive Resin Cement

- Industry-leading bond strength
- Fast and easy—only two components (automix syringe and bottle)
- Versatile—for total-etch, selective-etch or self-etch procedures
- Fluorescent shades deliver world-class esthetics
- Ideal for CAD/CAM and glass ceramic restorations



RelyX™ Unicem Self-Adhesive Resin Cement

RelyX™ Unicem 2 Self-Adhesive Resin Cement

- World's most clinically proven self-adhesive cement
- Saves time by eliminating etching, priming and bonding steps
- Virtually no post-operative sensitivity
- Strong and moisture-tolerant
- Endo tip makes cementing posts easier and more reliable

Veneer Cements



RelyX™ Veneer Cement

- Customer-preferred shading system
- Light-cure-only system provides excellent color stability
- Try-in pastes offer excellent color match to final cured cement
- Easy to dispense, apply and seat; easy cleanup

Resin-Modified Glass Ionomer Cements



RelyX™ Luting Plus Cement RelyX™ Luting Plus Automix Cement

- Automix syringe offers precision delivery and eliminates measuring, hand mixing and loading
- Automix cement features tack light cure option for faster and easy excess removal
- Paste/paste formula in the Clicker™ Dispenser provides consistent mix ratios
- Sustained fluoride release for caries-prone patients
- Virtually no post-operative sensitivity

Conventional Glass Ionomer Cements



Ketac™ Cem Radiopaque Glass Ionomer Luting Cement Ketac™ Cem Aplicap™/Maxicap™ Glass Ionomer Luting Cement

- Chemical bond to enamel and dentin
- Sustained fluoride release
- Extremely low film thickness optimizes marginal integrity
- Radiopaque
- Capsule delivery offers simple handling, consistent working characteristics, and reliable quality

Temporary Cements



RelyX™ Temp NE Temporary Cement

RelyX™ Temp E Temporary Cement

- Strong, temporary adhesion
- Quick, easy excess removal
- Optimal fit
- Easy provisional removal
- Available in both eugenol and non-eugenol formulas

3M ESPE Dental Cements. Easy to choose. Easy to use.

Recommended indications*

CAD/CAM
Ceramic and Composite

Resin
Nano Ceramics

Metal/PFMs

Porcelain/Ceramic
Inlays, Onlays and Crowns

Veneers

Zirconia

Implants



Posts

Pediatric Crowns
Metal

Maryland Bridges
Metal and Composite

Temporary
cementation of

- Provisional restorations
- Crowns, bridges, inlays and onlays

RelyX™ Ultimate Adhesive Resin Cements	RelyX™ Unicem AND RelyX™ Unicem 2 Self-Adhesive Resin Cements	RelyX™ Veneer Cement	RelyX™ Luting Plus AND RelyX™ Luting Plus Automix Cements	Ketac™ Cem Radiopaque AND Ketac™ Cem Aplicap™/Maxicap™ Glass Ionomer Luting Cements	RelyX™ Temp NE Temporary Cement RelyX™ Temp E Temporary Cement
			  	 	
•	•				
•	•				
•	•		•	•	
•	•				
•	•	•			
•	•		•		
•	•		•		
•	•				
•	•		•	•	
•	•				
					•
Delivery Method      					

*For complete list of indications, please refer to product's Instructions for Use.

Customer Care Center: 1-800-634-2249

www.3MESPE.com/RelyX

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
U.S.A.
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

Please recycle. Printed in U.S.A.
70-2013-0495-6

3M, ESPE, Aplicap, Clicker, Ketac, Maxicap and RelyX are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
© 3M 2012. All rights reserved.