

# MARK3<sup>®</sup>

Premium Products. Affordable Prices.

# MARK3 UniCore Universal

Micro-Hybrid Composite | 4 gm Syringe

MARK3 UniCore Universal Micro-Hybrid Composite delivers a versatile restorative material with superior esthetics for anterior and posterior restorations. Its balanced physical properties, enhanced polishability, high strength, and wear resistance make it a reliable choice for everyday restorative procedures.

## Features & Benefits

- For anterior & posterior restorations
- High polishability with great esthetics
- Does not stick to instruments
- Optimal balance of physical properties
- High strength and wear resistance
- Ultra-low shrinkage

### Performance Snapshot

Flexural strength: > 152 MPa

Water sorption : 1.8 - 14 ug/mm<sup>3</sup>

Water solubility: 1.2 - 3 ug/mm<sup>3</sup>

Curing time: Regular shades: 20 secs  
Opaque shades: 30



### Intended Use

- Class 1-V posterior & anterior restorations
- Veneering of discolored anteriors
- Splinting of loose teeth
- Anchoring of orthodontic appliances

**COMPARES TO**

Filtek Supreme, Dentsply Esthet-X, Kerr Herculite XRV Ultra

**PACKAGING**

Contains (1) 4gm Syringe

**High-strength universal composite with smooth handling and esthetics.**

MARK3 UniCore Universal Syringe