

- Passive parallel-sided post for direct one-office visit post/core build-ups
- Ideal for single and multi-rooted post/core restorations
- Efficient and smooth cutting burs

FIBER POST

- Radio-opaque and light transmitting
- Optimal mechanical retention
- Elasticity close to dentin
- Enhanced physical properties
- Bonds with Bis-GMA cements and core build-up materials
- Aesthetic thanks to its translucent matrix
- 100% biocompatible (none epoxy based matrix)



METAL POSTS

- Available in Stainless Steel and Titanium alloy (can be bent if needed)
- Number coded for easy identification
- Serrated for optimum retention
- Vented to prevent cement to cause apical hydrostatic pressure



ORDER INFORMATION

Item number	Product	Quantity	#	Ø mm	Ø inches
-------------	---------	----------	---	------	----------

GLASS FIBER POST

09-300	Introduction Kit: 7 drills (1 of each size), 21 posts (3 of each size) and 1 vial of endo stops.				
09-301	Parallel Sided Vented Post	5	#3	0,90	0,36
09-302	Parallel Sided Vented Post	5	#4	1,00	0,40
09-303	Parallel Sided Vented Post	5	#4.5	1,14	0,45
09-304	Parallel Sided Vented Post	5	#5	1,25	0,50
09-304	Parallel Sided Vented Post	5	#5.5	1,40	0,55
09-305	Parallel Sided Vented Post	5	#6	1,50	0,60
09-306	Parallel Sided Vented Post	5	#7	1,75	0,70

TITANIUM POST

09-100	Introduction Kit: 7 drills (1 of each size), 21 posts (3 of each size) and 1 vial of endo stops.				
09-101	Parallel Sided Vented Post	10	#3	0,90	0,36
09-102	Parallel Sided Vented Post	10	#4	1,00	0,40
09-103	Parallel Sided Vented Post	10	#4.5	1,14	0,45
09-104	Parallel Sided Vented Post	10	#5	1,25	0,50
09-105	Parallel Sided Vented Post	10	#5.5	1,40	0,55
09-106	Parallel Sided Vented Post	10	#6	1,50	0,60
09-107	Parallel Sided Vented Post	10	#7	1,75	0,70

STAINLESS STEEL POST

09-001	Introduction Kit: 7 drills (1 of each size), 21 posts (3 of each size) and 1 vial of endo stops.				
09-002	Parallel Sided Vented Post	10	#3	0,90	0,36
09-003	Parallel Sided Vented Post	10	#4	1,00	0,40
09-004	Parallel Sided Vented Post	10	#4.5	1,14	0,45
09-005	Parallel Sided Vented Post	10	#5	1,25	0,50
09-006	Parallel Sided Vented Post	10	#5.5	1,40	0,55
09-007	Parallel Sided Vented Post	10	#6	1,50	0,60
09-008	Parallel Sided Vented Post	10	#7	1,75	0,70

DRILLS

09-200	Parallel Sided Drill	2	#3	0,90	0,36
09-201	Parallel Sided Drill	2	#4	1,00	0,40
09-202	Parallel Sided Drill	2	#4.5	1,14	0,45
09-203	Parallel Sided Drill	2	#5	1,25	0,50
09-204	Parallel Sided Drill	2	#5.5	1,40	0,55
09-205	Parallel Sided Drill	2	#6	1,50	0,60
09-206	Parallel Sided Drill	2	#7	1,75	0,70

PLACE YOUR ORDER WITH YOUR PREFERRED DEALER

For product information contact Michael Nordahl
 Phone: 330 221 7213 Email: michael@nordin-dental.com

NORDIN®
 SWISS DENTAL PRODUCTS
 OF DISTINCTION

Harald Nordin SA
 Z.I. La Foge – CH-1816 Chailly – Switzerland
 Tél.: +41 21 964 71 71 – E-mail: info@nordin-dental.com