

Powerful cleaning with minimal effort.

VibraKleen E²® Ultrasonic Cleaning Tablets
Dri-Clave® VK Ultrasonic Cleaning Solutions
VacuKleen E²® Evacuation System Cleaner

- Powerful ultrasonic energy quickly removes and breaks down hard to clean deposits like blood, bone, and food particles that chemical detergents can leave behind
- Powerful, yet gentle on your instruments
- Biodegradable and pH neutral
- Non-toxic, non-corrosive formula contains no chlorine or harsh alkaline to damage your instruments



**VK Solutions
& Tablets**

Your Complete Ultrasonic Cleaning System

Your Complete Ultrasonic Cleaning System

Product	Features	Benefits	Uses	Packaging	Product Numbers
VIBRAKLEEN E²® ULTRASONIC CLEANING SOLUTIONS					
VibraKleen E² Ultrasonic Cleaning Solutions	<ul style="list-style-type: none"> • Effervescent, enzymatic cleaning tablets • Non-corrosive and low-foaming • Biodegradable, non-toxic & pH neutral • No chlorine or harsh alkaline • Drop directly into ultrasonic cleaner • Compact, convenient package of 20 tablets dilutes to 10 gallons • Clean mint scent • Effective as evacuation cleaner • Power of protease enzymes 	<ul style="list-style-type: none"> • Breaks down hard-to-clean deposits in ultrasonic and evacuation systems • Easy on instruments • Environmentally friendly and safe to work with • Eliminates damage to stainless steel equipment • No messy mixing of liquid concentrates. Easy to use • Saves on counter and storage space. No lifting heavy liquids • Minimizes any odor from tank, pleasant to work with • One product for multiple uses. Savings of time and money • Effective, fast removal of blood, bone and food particles that chemical detergents leave behind 	<ul style="list-style-type: none"> • For general purpose ultrasonic cleaning, evacuation system cleaning and pre-soak 	Box of 20 Box of 80 2 tablets = 1 gallon	50037100 —Box of 20 (Case Lot 24) 50037101 —Clinic Pack—Box of 80 (Case Lot 12)
DRI-CLAVE® VK ULTRASONIC CLEANING SOLUTIONS					
VK-1 General Purpose Cleaner	<ul style="list-style-type: none"> • Non-ionic • Non-sudsing • 10:1 concentrated formulation • Color coded solution • Anticorrosive 	<ul style="list-style-type: none"> • Protects your equipment • Avoid messy cleaning • Economical gallon container • Mistake-free identification 	<ul style="list-style-type: none"> • Removes debris from instruments, burs, general equipment, ortho appliances, gold and porcelain • Use directly in ultrasonic tank 	Gallon—blue Liquid	50036202*
VK-3 Plaster and Stone Remover	<ul style="list-style-type: none"> • Color coded solution • Non-sudsing • Efficient cleansing time (5–10) minutes • Gallon container • No diluting necessary—use at full strength 	<ul style="list-style-type: none"> • Mistake-free identification • Protects your instruments • Avoid messy cleaning • Saves time • Economical 	<ul style="list-style-type: none"> • For removal of plaster and stone • Do not use with aluminum • Place in glass beaker inside ultrasonic cleaner 	Gallon—pink Liquid	50036206*
VK-4 Tartar and Stain Remover	<ul style="list-style-type: none"> • Efficient cleansing time (5–10) minutes • Gallon container • No diluting necessary—use at full strength 	<ul style="list-style-type: none"> • Saves time • Economical 	<ul style="list-style-type: none"> • Removal of tartar and stains • 5–10 minutes only • Place in glass beaker inside ultrasonic cleaner 	Gallon—clear Liquid	50036208*
VK-6 Temporary Cement Remover	<ul style="list-style-type: none"> • Anticorrosive • Efficient cleansing time (5–10) minutes • Gallon container • No diluting necessary—use at full strength • Color coded solution 	<ul style="list-style-type: none"> • Protects your instruments • Saves time • Economical 	<ul style="list-style-type: none"> • Removal of temporary cement including cavit and eugonal based temporary cement • 5–10 minutes only • Place in glass beaker inside ultrasonic cleaner 	Gallon—yellow Liquid	50036212*
VACUKLEEN E²® EVACUATION SYSTEM CLEANER					
VacuKleen E² Evacuation System Cleaner	<ul style="list-style-type: none"> • Highly effective powdered formula • Non-corrosive and low-foaming • Biodegradable, non-toxic & pH neutral • Convenient, easy-mix packets • Power of protease enzymes • Fresh, clean, scent of lemon • Simple to use 	<ul style="list-style-type: none"> • Removes hard-to-clean protein-based deposits that soapy solutions can leave behind • Easy on instruments • Environmentally friendly and safe to work with • Pre-measured packets—saves time and minimizes waste. Just add warm water • Effective and quick dissolving of blood, tissue, saliva and encrusted residue in evacuation lines • Real lemon extracts gives a clean smell that eliminates evacuation system odor and makes the operatories smell fresh • Add powder packet to warm water to make a low-foaming solution. Only mix what you need so you eliminate waste and storage problems 	<ul style="list-style-type: none"> • Cleaning of your evacuation equipment 	Box of 25 pouches Box of 60 pouches Box of 140 pouches One envelope makes one quart of solution	50036100 —Unit Dose Pouch for Gallons Clinic Pack—Box of 25 Pouches (Case Lot 24) 50036101 —Unit Dose Pouch for Quarts Intro Pack—Box of 60 Pouches (Case Lot 24) 50036102 —Unit Dose Pouch for Quarts Clinic Pack—Box of 140 Pouches (Case Lot 12)

* Each VK Solution Sold in Case Lots of Four (4) Gallons Only.

The following information does not replace or supercede the “directions for use” contained in product packaging. Please refer to the “directions for use” for complete instructions
 VibraKleen E²®, VacuKleen E²® and Dri-Clave® are registered trademarks of Heraeus Kulzer, LLC.
 Heraeus Kulzer, LLC 300 Heraeus Way, South Bend, IN 46614 Phone: (800) 431-1785 Fax: (877) 271-5211 www.heraeusdentalusa.com www.SmileByVenus.com