

Visit [www.gcamerica.com](http://www.gcamerica.com) for special offers and introductory pricing on products featured in this flyer and more!



## GI LUTING CEMENTS

Self Cure



### The First Paste-Paste Automix Resin Reinforced GI Cement

**GC FujiCEM™ 2** – Resin reinforced glass ionomer cement in paste-paste automix delivery for a consistent mix everytime.

**Indications:**

- Cementing metal, PFM, resin, high strength ceramic crowns, bridges, inlays and onlays
- Cementing posts

FujiCEM (original) is also available to those who prefer to hand mix.

Self Cure



### The Ultimate Physical Properties for Resin Reinforced GI Cementation

**GC Fuji PLUS®** – Resin reinforced glass ionomer luting cement available in capsule delivery and powder-liquid.

**Indications:**

- Cementing metal, porcelain fused to metal PFM, resin, high strength ceramic crowns, bridges, inlays and onlays
- Cementing posts

Self Cure



**GC Fuji I®** – Radiopaque luting cement with snap set. Available in capsule delivery and powder-liquid.

**Indications:**

- Cementing metal, PFM, bridges, inlays and crowns
- Cementing orthodontic bands
- Cementing stainless steel crowns
- Base or liner

Low shrinkage

Rechargeable fluoride

Improved translucency

**GC - the world leader in glass ionomer technology introduces several significant advances.**

Biomimetic

Moisture friendly

Bulk filling

**A class of material who's time has come!**

Biocompatible

Self-adhesive

Fluoride releasing

Hydrophilic

Virtually no sensitivity

	EQUIA	Fuji IX GP EXTRA	Fuji II LC	Fuji LINING LC	Fuji TRIAGE
<b>Restorations:</b>					
Class III restorations					
Class V, root surface caries, cervical erosion					
Class I and II pediatric restorations					
Long-term temporary restoration					
Class I and II in permanent teeth					
Core build-ups					
Bases and liners					
Sandwich technique					
Liner under composite or amalgam					
Base under amalgam, composite, porcelain					
Inlays and onlays					
Fissure sealant, root surface protection					
Protection for partially erupted molars					



New R&D Center

See inside for more information.



800.323.7063  
[www.gcamerica.com](http://www.gcamerica.com)  
[www.gcamerica.com/training](http://www.gcamerica.com/training)  
 © 2014 GC America Inc.



## Glass Ionomer Material Science Continues to Advance and Grow in Popularity Worldwide...

### Here's Why:

FEATURES	BENEFITS
Self-adhesive, excellent marginal seal	<b>No need for bonding agents</b>
Continual rechargeable Fluoride release	<b>Helps prevent secondary caries</b>
Biocompatible material with coefficient of thermal expansion like tooth structure	<b>Long lasting marginal integrity</b>
Fast and easy moisture tolerant technique	<b>Time saving/no need for isolation</b>
The new generation of glass ionomer	<b>Offering vastly improved esthetics and strength indicated for classes I, II, III, IV, &amp; V</b>

## GI RESTORATIVES

**Stronger, Smoother & Improves Aesthetics**

**NEW!**

**EQUIA™** – Self-adhesive aesthetic posterior restorative system. Bulk-filled two-step approach by combining EQUIA™ Fil and EQUIA™ Coat to improve hardness by 90%. Available in capsule delivery.

#### Indications:

- Class I, II and V restorations

Fast placement: Class I and II restorations can be placed in less than 3 minutes and 30 seconds.



Self Cure



## More Translucency, More Fluoride, Fast Setting

**GC Fuji IX GP®, Fuji IX GP® FAST, Fuji IX GP® EXTRA** – Packable posterior restorative. Self-curing. Available in eight Vita® shades in powder-liquid (GC Fuji IX GP only) and capsule delivery.

#### Indications:

- Class I and II restorations (non stress bearing)
- Core build-ups
- Pediatric / geriatric restorations
- Base or liner (sandwich technique)
- Intermediate restorative (IRM)
- Long term temporization

GC Fuji IX GP FAST sets twice as fast as GC Fuji IX GP, GC Fuji IX GP EXTRA has an even faster setting time, releases more fluoride and is extra translucent.

## The Standard in Resin Reinforced GI Materials

**GC Fuji II LC®** – Resin reinforced restorative with excellent aesthetics. Rapid light curing. Available in ten Vita® shades in capsule delivery and powder-liquid.

#### Indications:

- Class III and V restorations
- Cervical erosions
- Root surface lesions
- Pediatric restorations
- Base or liner (sandwich technique)

Light Cure



## GI BASE & LINERS

Light Cure



**GC Fuji LINING™ LC** - Radiopaque light-cured lining material. Available in paste-paste delivery and powder-liquid.

#### Indications:

- Base and liner under composite resins, amalgams and porcelains

## GI SEALANT / SURFACE PROTECTION

Self Cure



**GC Fuji TRIAGE®** – Innovative sealant and surface protectant. Works in moist environment. Self-curing. Available in capsule delivery in pink or white.

#### Indications:

- Partially erupting molars
- Provides lasting relief of hypersensitivity
- Temporary endodontic sealing - provides effective seal between endodontic visits
- Intermediate restorative for rampant caries - first step in overall treatment, stops progression to caries

Photos courtesy of Dr. Stanislaw Lassocinski



Glass Ionomer restoration done with GC Fuji IX GP® EXTRA.



Glass Ionomer restoration done with EQUIA.