

SILMET

Redefining Dental Adhesives for Superior Results!



ProLink™ Single Step Dental Adhesive

ProLink™ utilizes advanced Colloidal Glass Technology to provide a first-class dental adhesive. It is easy to use and works for every restoration.

4.3
out of 5

Dental **PRODUCT SHOPPER**
BEST PRODUCT

SUITABLE FOR:

- Direct Composite or Compomer Restorations.
- Adhesive cementation.
- Composite repair.

ProLink™ delivers outstanding shear bonding values to dentin, enamel & all popular light, dual or auto-cure composites.



Bond Strength to Dentin: > 14 MPa
Bond Strength to Enamel: > 20 MPa

FEATURES & BENEFITS:

- Acetone Solvent.
- HEMA Functional Group.
- 5th Generation.
- Light-Cured.
- Compatible with Light, Dual, and Self Cured Materials.
- Compatible with Both Direct & Indirect Restorations.
- Minimal treatment time.
- Increased adhesion & durability.
- Low sensitivity.

ProLink product line incorporates innovative Chemistry for every restoration challenge and is designed for various needs and preferences.

MORE FROM THE PROLINK ADHESIVES FAMILY:



10 MDP
MONOMER
Adhesion to various substrates

ProLink™ Universal 8th Generation Adhesive

To be used with self or dual-cure composite in any bonding technique.
One adhesive in the clinic for all restorations.



EDITOR'S CHOICE

ProLink™ SE

One Component Self Etching Light-cure Dental Adhesive

One component dentin enamel bonding agent. A single component that combines etch, primer, and Adhesive, for ease of application.



ProLink™ Ethanol

Ethanol Based Single Step Bonding Agent

Direct Composite or Compomer restorations in the total-etch bonding technique.

Description	SKU	Quantity
ProLink	150005	5ml
ProLink Ethanol	150005-ET	5ml
ProLink SE	15SEKT	4ml
ProLink Universal	1505UN	5ml
	1503UN	3ml