

Tray Aways®

Make a Great First Impression

Tray Aways Disposable Impression Trays are
THE original plastic tray!

Updated throughout the years, used not only in dental but also with implants. It has earned its place in your office.

Tray Aways are the disposable impression trays of choice suited for use with alginate or other types of impression materials where detailed reproduction and dimensionally-stable results are vital. Made from rigid, high-impact FDA approved styrene, Tray Aways assist in generating a consistent, unsurpassed impression free of distortion or bite through. Smooth, rounded edges, interior retention slots, and vertical perforation enhance material retention. The best value on the market.

Can be bagged in quantities other than 12.



FEATURES AND BENEFITS

- Made from rigid, high-impact styrene
- Produces dimensionally-stable results
- Creates impression free of distortion & bite through
- Also used with implants

ORDERING INFORMATION

#1 Perforated, Large Upper – 3.06" x 2.41" – 12/bag	#0921884
#2 Perforated, Large Lower – 3.02" x 2.30" – 12/bag	#0921887
#3 Perforated, Medium Upper – 2.91" x 2.21" – 12/bag	#0921885
#4 Perforated, Medium Lower – 2.93" x 2.18" – 12/bag	#0921888
#5 Perforated, Small Upper – 2.66" x 2.00" – 12/bag	#0921886
#6 Perforated Small Lower – 2.78" x 2.06" – 12/bag	#0921889
#7 Solid Partial Upper R, Lower L – 2.24" x 2.29" – 12/bag	#0921890
#8 Solid Partial Upper L, Lower R – 2.15" x 2.33" – 12/bag	#0921891
#9 Solid Lower Anterior – 2.56" x 1.50" – 12/bag	#0921892
#10 Solid Upper Anterior – 2.52" x 1.47" – 12/bag	#0921893
#7s Perforated Partial Upper R, Lower L – 12/bag	#0921894
#8s Perforated Partial Upper L, Lower R – 12/bag	#0921895
#9s Perforated Lower Anterior – 12/bag	#0921896
#10s Perforated Upper Anterior – 12/bag	#0921897

