

AURELIA® Protegé®

Examination Gloves



93995	(XS)	ø 100	10 x 100
93996	(S)	ø 100	10 x 100
93997	(M)	ø 100	10 x 100
93998	(L)	ø 100	10 x 100
93999	(XL)	ø 100	10 x 100



Nitrile



Powder Free



Sky Blue



4.4 +/- 0.2 g



Finger-Textured




5 Years

Premium Quality Sky Blue Nitrile examination gloves that offer greater tensile strength and protection than both Latex and Vinyl.

Nitrile examination gloves offer extra comfort and flexibility.

- 4.0mil Powder-Free Nitrile
- Non-Sterile
- Beaded cuff
- Ambidextrous
- Latex Free
- AQL 1.5

Sizes & Dimensions



XS	70 ± 10 mm	min. 240 mm
S	80 ± 10 mm	min. 240 mm
M	95 ± 10 mm	min. 240 mm
L	110 ± 10 mm	min. 240 mm
XL	≥ 110 mm	min. 240 mm

Single wall thickness



- ① min. 0.12 mm
- ② min. 0.08 mm
- ③ min. 0.05 mm

Physical properties

	Unaged	Aged
Tensile Strength	min. 14.0 MPa	min. 14.0 MPa
Elongation	min. 500 %	min. 400 %
Strength at Break	min. 6N	min. 6N

Quality Signs & Certificates

- 89/686/EEC, 93/42/EWG
- ISO 9001:2008, ISO 13485, CMDCAS 13485 (SGS)
- EN 455 1-4
- EN 420 : 2003 + A 1 : 2009
- EN 374 2-3, EN 388
- EN 1186: Food Contact Approved (SGS)

