

Ketac™ Silver Aplicap™/Maxicap™

Glass Ionomer Filling Material



When patient satisfaction is your #1 priority...use the #1 metal-reinforced glass ionomer solution

- Contains no amalgam
- High sustained fluoride release
- Chemical bond to enamel and dentin provides excellent marginal integrity
- Radiopaque
- Aplicap and Maxicap delivery systems provide for consistent mixing and easy dispensing
- Maxicap provides four times the volume of Aplicap. Designed for larger scale restorations and core build-ups.



Ketac™ Silver Aplicap™/Maxicap Glass Ionomer Filling Material is a proven alternative to amalgam for core build-ups.

Silver particles are sintered into the glass for excellent thermal and electrical insulation—thus reducing post-operative sensitivity.

Moisture tolerance, long-term fluoride release, and adhesion to enamel and dentin give you added confidence of sealed, lasting restorations.



Ketac Silver material is a fluoride-releasing alternative to amalgam for large fillings.

Clinical advantages:

- Ketac Silver is biocompatible—this makes it particularly suitable for core build-ups
- An excellent insulator, which protects the pulp and helps eliminate patient sensitivity
- No need for undercutting—so you conserve tooth structure
- An easy-to-mix hand mix formula and two capsule delivery systems makes it versatile in any application
- Less shrinkage than traditional composites—leads to better marginal integrity

Indications:

- Core build-ups
- Primary teeth fillings
- Small fillings in posterior teeth
- Temporaries
- Base under amalgams

3M ESPE

Glass Ionomer solutions from 3M ESPE produce the results you depend on.

After more than 30 years developing glass ionomer products, 3M ESPE has you covered. Use this handy reference to quickly determine which material is best suited to meet your procedural needs.

Indications									
	Small Class I	Class III	Class V	Minimal Invasive Dentistry (MID)	Primary Teeth	Base Under Composite	Temporary Treatment	Core Build-up	Delivery
RESIN MODIFIED									
Ketac™ Nano	●	●	●	●	●	●	●	●	
Photac™ Fil	●	●	●	●	●		●		
Vitremer™	●	●	●	●	●	●	●	●	
CONVENTIONAL									
Ketac™ Molar Quick	●		●	●	●	●	●	●	
Ketac™ Molar	●		●		●	●	●	●	
Ketac™ Fil Plus	●	●	●	●	●			●	
METAL REINFORCED									
Ketac™ Silver	●		●		●		●	●	



Quick Mix Capsule



Clicker™ Dispenser



Capsule



Powder Liquid

Ordering Information



Ketac™ Silver Aplicap™/Maxicap™ Introductory Kit (37006)



3M™ ESPE™ RotoMix™ Capsule Mixing Unit (76310)

Item Number	Product Information
37006	Ketac™ Silver Aplicap™ Silver GI Cement Restorative Intro Pack Contains: 50 Capsules (37011); 1 Activator (73040); 1 Applier (73050).
37010	Ketac™ Silver Aplicap™ Silver GI Cement Restorative Capsule Refill (50)
56203	Ketac™ Silver Maxicap™ GI Cermet Restorative Intro Kit Contains: 50 Capsules; 1 Activator; 1 Applier
56210	Ketac™ Silver Maxicap™ GI Cermet Restorative Refill (50)
37800	Ketac™ Silver Glass Ionomer Intro Kit: (Handmix) Contains: 1–25g Powder; 1–12mL Liquid; 1 Mixing Pad; 1 Spoon
37810	Ketac™ Silver Powder (25g)
37820	Ketac™ Silver Liquid (12mL)
37470	Ketac™ Conditioner (10mL)
37190	Ketac™ Glaze, 2.5g Bottle
76310	3M™ ESPE™ RotoMix™ Capsule Mixing Unit At the touch of a button, the RotoMix Capsule Mixing Unit provides a void-free, completely homogenous mix which is easier to handle so you do your best work. The unique design of 3M™ ESPE™ Aplicap™ and Maxicap™ capsules prevents “duds” – saving your practice time and money.

Customer Care Center: 1-800-634-2249 www.3MESPE.com

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
U.S.A.
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

Please recycle. Printed in U.S.A.
© 3M 2011. All rights reserved.
70-2009-3715-2

3M, ESPE, Aplicap, Ketac, Maxicap, Photac and Vitremer are trademarks of 3M or 3M ESPE AG. Used under license in Canada.