



Cotton / Non-Woven Sponges

UltraPure® Non-Woven Sponges 4 Ply, Non-Sterile

Made from 100% USA-grown organic pure cotton.

Crosstex introduces a great way to reduce your carbon footprint: eco-friendly UltraPure Non-Woven Sponges. They're soft, lint-free and absorb faster than synthetics, rayon and cotton blends. And, unlike typical synthetic non-woven sponges, they're biodegradable. So you can provide improved treatment for your patients...and protect the environment!

- **Superior, faster absorbency and fluid retention** (40-45 GSM) compared to synthetics, rayon and cotton blends.
- **Improved infection control** – Exclusive automatic manufacturing process reduces cross-contamination risks and provides low bioburden levels.
- **Safe for all** – Hypoallergenic; no toxic chemicals used in manufacturing – all-natural whitening process; no bleach or chlorine; pesticide-free cotton used.
- **Multi-purpose** – Use for cleaning, prepping, placement on incision/extraction sites.

TOP GAUZE 3 YEARS IN A ROW



Ordering Information

| Ref. # | Description | Quantity |
|---------|----------------------|-----------|
| ENC100 | 2" x 2" / 5 x 5 cm | 5000/case |
| ENC4100 | 4" x 4" / 10 x 10 cm | 2000/case |



100% cotton grown in USA



Ultra Gauze® Sponges 4 Ply, Non-Sterile

- **More absorbent than traditional sponges.** Unique multi-layer design with 40 gram basis weight offers faster absorption and higher fluid retention. 70/30 rayon/polyester formed non-woven material replicates the look of cotton woven gauze.
- **Releases easily** from incision/extraction sites.
- **Extremely soft;** produces little or no lint.
- **Multi-purpose** – can be used for cleaning, prepping and placement on incision/extraction sites.


Ordering Information

| Ref. # | Description | Quantity |
|---------|--------------------------|-----------|
| ENCNWU | 2" x 2" / 5 x 5 cm | 5000/case |
| ENCNWUB | 2" x 2" / 5 x 5 cm | 8000/case |
| ENC3NWU | 3" x 3" / 7.62 x 7.62 cm | 4000/case |
| ENC4NWU | 4" x 4" / 10 x 10 cm | 2000/case |



Non-Woven Sponges 4 Ply, Non-Sterile

- **Less adhesion to wounds** for faster healing and virtually lint free.
- **Ideal for sustained patient use**, disinfecting and general cleaning.
- **Excellent absorbency** with unique 70/30 rayon/polyester blend.
- **Choice of three weights** – Premium (40-50 gram weight); Advantage Plus® (35 gram weight); Advantage (30 gram weight).



Ordering Information

| Ref. # | Description | Quantity |
|---|------------------------------------|-----------|
| PREMIUM 4 Ply Non-Sterile - Made in U.S.A. | | |
| ENCNW | 2" x 2" / 5 x 5 cm - 50 GR WT* | 5000/case |
| ENC3NW | 3" x 3" / 7.5 x 7.5 cm - 40 GR WT* | 4000/case |
| ENCNW4K | 2" x 2" / 5 x 5 cm - 50 GR WT* | 4000/case |
| ENC4NW | 4" x 4" / 10 x 10 cm - 40 GR WT* | 2000/case |
| ADVANTAGE PLUS® 4 Ply Non-Sterile - Made in U.S.A. | | |
| ENCNWL | 2" x 2" / 5 x 5 cm - 35 GR WT* | 5000/case |
| ADVANTAGE 4 Ply Non-Sterile | | |
| ENCNWLA | 2" x 2" / 5 x 5 cm - 30 GR WT* | 5000/case |
| ENC3NWLA | 3" x 3" / 7.5 x 7.5 cm - 30 GR WT* | 4000/case |
| ENC4NWLA | 4" x 4" / 10 x 10 cm - 30 GR WT* | 2000/case |

* Gram Weight: Weight per square meter of rayon/poly fabric.



Cotton-Filled Exodontia Sponges 4 Ply, Non-Sterile

- **Pure, virgin white cotton** with highly absorbent 100% cotton-fiber filling - Made from USP* cotton gauze
- **Choice of three materials** – Premium, Sterile and Advantage Plus®

Ordering Information

| Ref. # | Description | Quantity |
|---|-------------------------------------|-----------|
| ADVANTAGE PLUS® Non-Sterile - Made in U.S.A. | | |
| ENACC | 2" x 2" / 5 x 5 cm - 90 GR WT* | 5000/case |
| PREMIUM Non-Sterile - Made in U.S.A. | | |
| ENCC | 2" x 2" / 5 x 5 cm - 125 GR WT* | 5000/case |
| ENCC3 | 3" x 3" / 7.5 x 7.5 cm - 125 GR WT* | 4000/case |
| ENCC4 | 4" x 4" / 10 x 10 cm - 125 GR WT* | 2000/case |
| STERILE | | |
| ESCC | 2" x 2" / 5 x 5 cm - 125 GR WT* | 5000/case |

* Gram Weight: Weight per sleeve.



All-Gauze Sponges Premium, Non-Sterile

- **Highly absorbent** – 100% gauze - Made from USP* 100% gauze
- **Choice of three plies** – 8 ply, 12 ply, 16 ply.

Ordering Information

| Ref. # | Description | Quantity |
|--------|---------------------------------|-----------|
| ENC212 | 2" x 2" / 5 x 5 cm - 12 ply | 8000/case |
| ENC2 | 2" x 2" / 5 x 5 cm - 8 ply | 5000/case |
| ENC312 | 3" x 3" / 7.5 x 7.5 cm - 12 ply | 4000/case |
| ENC38 | 3" x 3" / 7.5 x 7.5 cm - 8 ply | 4000/case |
| ENC416 | 4" x 4" / 10 x 10 cm - 16 ply | 2000/case |
| ENC412 | 4" x 4" / 10 x 10 cm - 12 ply | 2000/case |
| ENC48 | 4" x 4" / 10 x 10 cm - 8 ply | 4000/case |



Manufactured by or for Crosstex - ©2015 Crosstex - CTX2592 - 3/4/2015

= NOT MADE WITH NATURAL RUBBER LATEX

* United States Pharmacopeia (USP) cotton/cotton gauze defines specific standards based on thread count, length, width, weight, absorbency, ignited residue, fatty matter and alcohol-soluble dyes. Manufacturers of this material must conduct ongoing tests that assure the product meets and/or exceeds the standards established by the USP. All Crosstex sponges are manufactured in facilities that have an in-house laboratory that continually monitors for USP compliance.