

INTRODUCING...

IMPRESSION MATERIAL

Alginate

Features and Benefits:



- Faster mixing time, saves you time
- Faster setting time, decreases chances of gag reflex
- Better tear and compression strength, resulting in less tearing or deformation
- Dust free
- More material, 500 grams compare to 450 grams
- 9 days dimensional stability

Tearing or deformation of alginate impressions can result in inaccuracy of gypsum casts that form the initial base for prosthetic appliances. Strength and elasticity of alginate impression materials are important characteristics determining the performance of these products.

The MARK3 Line was created to provide High Quality Products at an Affordable Price

Look at the advantages of MARK3 Alginate

Properties	MARK3 Fast	Kromopan from Lascod	Jeltrate Plus Fast from Dentsply
Mixing time	30 secs.	45 secs.	45 secs.
Setting time in mouth	1 min.	1.45 min	1.45 min.
Compressive strength	1.2 MPa	0.87 MPa	0.68 MPa
Tear strength	1,11 N/mm	0,81 N/mm	0,45 N/mm
Nominal weight package	500g	450g	454g
Dustfree	yes	no	yes