

# Maxso® EP12-NS

## Hygiene Handpieces

Beyes EP12-NS hygiene handpiece is a 4:1 reduction straight nose cone designed to be used with disposable prophylaxis angle for polishing. The EP12-NS offers the highest quality ceramic bearings for long-lasting smooth operation and a slim, lightweight, and compact handpiece design for ultimate user comfort.

- Beyes EP12-NS will fit all standard ISO E-Type air motor or electric micromotor
- 4:1 reduction straight nose
- Ceramic bearings
- Amazingly quiet operation
- Autoclavable
- 1 Year warranty

### Order Info

**Code:** SL1002  
**Model:** EP12-NS



For Disposable  
Prophy Angles



### Key Specification

<b>Transmission ratio</b>	4:1
<b>Connection</b>	ISO E-Type
<b>For use with</b>	Disposable prophylaxis angles
<b>Speed type</b>	Low speed
<b>Maximum speed</b>	10,000 rpm
<b>Optic</b>	Non optic
<b>Overall length</b>	64 mm
<b>Shape</b>	Straight
<b>Weight</b>	57 g
<b>Sterilization</b>	Autoclavable up to 135°C
<b>Warranty</b>	1 year warranty

**Beyes®**

*The Next Standard*

**Beyes Dental Canada Inc.**  
23-595 Middlefield Road  
Toronto, Ontario, M1V 3S2  
Canada  
**Tel:** 1-416-281-1888  
**Email:** info@beyes.ca  
**Website:** www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for more information on EP12-NS



Printed in Canada  
Document # ENB021  
Rev.2/21.08.19

© Beyes Canada Inc. 2019. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under Canadian Medical Devices Conformity Assessment System (CMDCAS).