

# MARK3<sup>®</sup>

Premium Products. Affordable Prices.

## MARK3 UniCore X5 Bond

5th Generation Total Etch Bonding Agent | 5 mL Bottle

MARK3 UniCore X5 Bond is a 5th generation total-etch adhesive for direct and indirect restorations, offering strong bond strength to enamel, dentin, and popular composites. Its one-bottle formula allows fast, single-step application, while Colloidal Glass Technology helps reduce polymerization shrinkage.

### Features & Benefits

- Primer and sealer combined in one bottle
- One application only
- Reduced polymerization shrinkage
- Less sensitivity
- For self-cure, light-cure, or dual-cure modes
- Easy flip-top bottle

### Performance Snapshot

Shear bond strength - Dentin >14 MPa

Shear bond strength - Enamel >20 MPa

Advanced colloidal glass technology helps reduce polymerization shrinkage.



### Intended Use

- Direct Composite or Compomer restorations
- Adhesive cementation
- Composite repair

**COMPARES TO** OptiBond, Prime & Bond NT, Single Bond

**PACKAGING** Contains (1) 5 mL bottle

Premium bonding performance with efficient one-step handling.

MARK3 UniCore X5 Bond