

# SAFETY DATA SHEET

## SECTION 1 - General Information

<b>Product Name:</b> AIR TECHNIQUES FIXER	<b>Catalog #:</b> 40245, 40270
Product Use: Photographic Fixer	HMIS SCALE: Health 1      Flammability 0
Supplier: Air Techniques Inc.	Reactivity 0      PPE C
Address: 1295 Walt Whitman Road, Melville, NY 11747	
Phone No.: (516)433-7676	Emergency Tel.: CHEMTREC 800-424-9300

## SECTION 2 - Hazards Identification

Precautionary Statement.	Classification	WARNING
Oral: Harmful if swallowed or inhaled.	4	
Skin: Causes skin irritation	2	
Eye: Causes eye irritation	2B	

### EFFECTS OF CHRONIC EXPOSURE TO PRODUCT

None known

## SECTION 3 - Product & Hazardous Ingredients Information

HAZARDOUS INGREDIENTS	CAS No.	W/W %	PEL	LD50	LC50
Water	7732-18-5	90-95	not listed	42.6 gm/kg	not available
Ammonium Thiosulfate	7783-18-8	8-12	not listed	2890 mg/kg	not available
Sodium Sulfite	7757-83-7	<2	not listed	2825 mg/kg	not available
Sodium acetate	127-09-03	<2	not listed	3530 mg/kg	30 gm/m <sup>3</sup>

## SECTION 4 - First Aid Measures

### SKIN CONTACT:

Wash exposed skin thoroughly with soap and water. If irritation persists seek medical advice.

### EYE CONTACT:

Flush eyes with water for 15 minutes. Seek medical advice if symptoms persist.

### INGESTION:

DO NOT induce vomiting. Seek medical attention immediately.

### INHALATION:

If inhaled, remove to fresh air.

## SECTION 5 - Fire Fighting Measures

Flammability    YES     NO     IF YES UNDER WHAT CONDITIONS

### MEANS OF EXTINGUISHING

As needed to extinguish adjacent fire source.

### FLASH POINT (C/F)

Upper limit % by volume

Lower limit % by volume

Auto ignition Temperature (C/F)

Not appropriate

Not appropriate

Not appropriate

Not appropriate

### HAZARDOUS COMBUSTION PRODUCTS

Explosion Data    Sensitivity to impact

Sensitivity to Static Discharge

"See Hazardous Decomposition Products"

not susceptible

not susceptible

## SECTION 6 - Accidental Release

Flush to sewer with plenty of water, if permitted. If not, soak up with dry absorbent material. Consult with governmental regulatory agencies for appropriate disposal of this material. Wear appropriate protective equipment.

## SECTION 7 - Handling and Storage

Read and follow the label, safety data sheet, and instructions before using. Avoid any contact with skin and eyes.

Keep container tightly closed. Avoid incompatible substances. Wash thoroughly after handling this product.

**Product Name:** Air Techniques Fixer **Catalog #:** 40245, 40270

**SECTION 8 - Exposure Control**

**VENTILLATION** Good general ventilation  
**SKIN** Latex or Nitril **RESPIRATORY** Not normally required  
**OTHER** Eye wash station **EYE:** Splash goggles

**SECTION 9 - Physical & Chemical Properties**

<b>PHYSICAL STATE</b> Liquid	<b>ODOR &amp; APPEARANCE</b> Clear, colorless with slight vinegar odor.	<b>ODOR THRESHOLD</b> Not available	
<b>VAPOUR PRESSURE (mm-Hg)</b> <17 mm Hg	<b>Vapor Density (Air=1)</b> 0.6 mm Hg	<b>Evaporation Rate (Water=1)</b> Not Available	<b>Boiling Point</b> > 212 F
<b>pH</b> 4.3	<b>Specific Gravity</b> 1.085	<b>Solubility in Water</b> Complete	<b>Freeze Point</b> < 32 F
		<b>Coef H2O/oil</b> N/A	<b>% Volatile (by volume)</b> 50%

**SECTION 10 - Stability & Reactivity Data**

**Chemical Stability** IF NO, UNDER WHAT CONDITIONS  
 YES  NO   
**INCOMPATIBLE WITH OTHER SUBSTANCES** IF YES, WHICH ONES?  
 YES  NO  Strong reducing agents & Strong bases.  
**HAZARDOUS DECOMPOSITION PRODUCTS:** Sulfur dioxide & ammonia  
**REACTIVITY:** Stable under ambient pressure & temperature.

**SECTION 11 - Toxicological Information**

<b>EXPOSURE LIMITS</b> None Known	<b>Irritancy of Product</b> Skin: Irritant	<b>Sensitation to Product</b> None known	<b>CARCINOGENICITY</b> Not listed by IARC, OSHA, or ACGIH
<b>TERATOGENICITY</b> None Known	<b>REPRODUCTIVE TOXIN</b> Not Known	<b>MUTAGENICITY</b> Not Known	<b>SYNERGISTIC PRODUCTS</b> Not Known

**SECTION 12 - Ecological Information**

No Data Available

**SECTION 13 - Disposal Considerations**

Dispose of this product in accordance with federal, state and/or local laws.

**SECTION 14 - Transportation Information**

Not regulated for transportation.

**SECTION 15 - Regulatory Information**

SARA: Reportable quantity (RQ) for Ammonium Thiosulfate 5000 lb.  
 SARA Section 304, CERCLA Hazardous components - Ammonium Thiosulfate  
 SARA Section 302,304,313 Hazardous components (None)  
 Not a controlled product under Canadian WHMIS regulation

**SECTION 16 - Other Information**

The information contained in this SDS is furnished without warranty of any kind. The user should consider this data a supplement to other data gathered and must make an independent determination of suitability from this and other sources to assure proper use and disposal of this material and the health and safety of employees and customers.