

SAFETY DATA SHEET

FS Hemostatic™

1. IDENTIFICATION

Product Identifier: 9007057, 9007058
Product name: FS HEMOSTATIC

Premier® Dental Products Company
1710 Romano Drive
Plymouth Meeting, PA 19462 USA

Phone: 610-239-6000
Fax: 610-239-6171
Emergency Phone: 610-239-6000

Restrictions for Use: Professional hemostatic agent

Recommended Use: For professional use only.

2. HAZARDS IDENTIFICATION

Classification of the chemical

Skin Irritant Cat 2

Eye Irritant Cat 2



Hazard pictogram(s) Exclamation Mark

Signal word: Warning

3. COMPOSITION AND INGREDIENTS

Chemical Name:	CAS #	WT %	Classifications
Ferric Sulfate	15244-10-7	10-30%	Skin Irritant Cat 2 Eye Irritant Cat 2

4. FIRST AID INFORMATION

(SYMPTION)

Skin contact: There may be mild irritation at the site of contact.

Eye contact: There may be irritation and redness.

Ingestion: There may be irritation of the throat.

Inhalation: No symptoms

(ACTION)

Skin contact: Wash immediately with plenty of soap and water.

Eye contact: Bathe the eye with running water for 15 minutes.

Ingestion: Wash out mouth with water.

Inhalation: Consult a doctor.

5. FIRE-FIGHTING MEASURES

Extinguishing Media: Suitable extinguishing media for the surrounding fire should be used.
Use water spray to cool containers.

Exposure Hazards: In combustion, emits toxic fumes.

Protection of Fire-fighters: Wear self-contained breathing apparatus. Wear protective clothing to prevent contact with skin and eyes.

6. ACCIDENTAL RELEASE MEASURES

Personal Precautions: Refer to Section 8 of SDS for personal protection details. Turn leaking containers leak-side up to prevent the escape of liquid.

Environmental Precautions: Do not discharge into drains or rivers. Contain the spillage using bunding.

Clean-up Procedures: Absorb into dry earth or sand. Transfer to a closable, labeled salvage container for disposal by an appropriate method.

7. HANDLING AND STORAGE

Storage Conditions: Store in cool, well ventilated area. Keep container tightly closed.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection: Respiratory protection not required.
Hand Protection: Protective Gloves
Eye Protection: Safety glasses. Ensure eye bath is at hand.
Skin Protection: Protective clothing.

9. PHYSICAL AND CHEMICAL PROPERTIES

State: Liquid

10. STABILITY AND REACTIVITY

Stability: Stable under normal conditions
Conditions to avoid: Heat
Materials to avoid: Strong oxidizing agents. Strong acids.
Haz. Decomp Products: In combustion, emits toxic fumes.

11. TOXICOLOGICAL INFORMATION

Routes of Exposure: Refer to Section 4 of SDS for routes of exposure and corresponding symptoms.

12. ECOLOGICAL INFORMATION

Mobility: Readily absorbed into soil
Persistence & Degradability: Biodegradable
Bioaccumulative Potential: No bioaccumulation potential
Other Adverse Effects: Negligible ecotoxicity

13. DISPOSAL CONSIDERATIONS

NB: The user's attention is drawn to the possible existence of regional or national regulations regarding disposal.

14. TRANSPORT INFORMATION

Product Not classified as dangerous in the meaning of transport regulations

ADR/RID: n/a
UN no: -
IMDG/IMO: n/a
IATA/ICAO: n/a

15. REGULATORY INFORMATION

Risk phrases:
R36/38: Irritating to eyes and skin

Hazard Statements:
H319 Causes serious eye irritation
H315 Causes skin irritation

Precautionary phrases:
P264 Wash thoroughly after handling.
P280 Wear protective gloves, protective clothing, eye protections and face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313 If eye irritation persists: Get medical attention.
P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P332+P313 If skin irritation occurs: Get medical attention.
P362 Take off contaminated clothing and wash before reuse.

Note:The regulatory information given above only indicates the principal regulations specifically applicable to the product described in the safety data sheet. The user's attention is drawn to the possible existence of additional provisions which complete these regulations. Refer to all applicable national, international and local regulations or provisions.

16. OTHER INFORMATION

Premier's revision date: 11/25/2015

Revision number: 4

Supplier's number: 090910

The information contained herein is based on our present knowledge. However, this information shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. Premier Dental Products Company makes no warranties, express or implied with respect to, and assumes no responsibility or liability for, the accuracy or completeness of the information contained herein. Premier Dental Products Company urges persons receiving this information to make their own determination as to the information suitability for their particular application.

111560 Rev4