



NOBLE
DENTAL SUPPLIES

PREMIUM DENTAL SUPPLIES - SIGNIFICANT DISCOUNTS
TRUSTED BY 1,000+ DENTISTS

2018 CATALOG

135 Route 9W, Congers, New York, 10920

Phone 866.333.6825

Web www.Nobledentalsupplies.com



MISSION STATEMENT

We have been serving the global dental community since 1978. Over the last 40 years we have seen the challenge that dentists have faced in finding high quality, affordable dental products. Name-brand manufacturers have become expensive, increasing the cost of their products every year regardless of the income of dental practices.

For these reasons, we have used our extensive knowledge in sourcing to bring dentists high quality products across a wide range of categories that are 30-50% below the expensive name-brand products.

OUR HISTORY

1978

Our corporate office was founded in 1978 by Carol and Gus Markiet out of their home in Pomona New York.

They are well respected in the dental industry, delivering what they promise; honesty, trust and reliability.

1994

Carol and Gus Markiet's 3 children come into the business. Our focus then changed from exporting dental supplies to selling domestically across the United States and Canada.

2004

We launched our retail division, Dental Deals and began selling low priced high quality dental supplies directly to dentists.

2006

We relocated to Congers, New York to a 15,000 square foot facility to accommodate our fast growing dental supply business. In the same year the Dental Deals website was launched and our company was renamed "Dental Deals Online". Our website has brought us thousands of customers. But it is our values, instilled in us by our parents 40 years ago - honesty, trust and reliability that has kept our satisfied customers loyal.

2018

Dental Deals Online has rebranded as **Noble Dental Supplies.**

WHY THE NAME CHANGE?

Dental Deals Online was always a promise to our customers more than a classic brand name: fantastic deals on dental supplies for dentists nationwide.

Online commerce has matured over the last decade, and we think it's time to modernize our brand presence as well. So instead of a descriptive but impersonal name for our business, we have updated our name to Noble Dental Supplies to reflect our family-driven integrity.

TABLE OF CONTENTS

Acrylics & Reline Materials	8
Alloys.....	10
Anesthetics	14
Articulating Paper/Accessories	15
Burs & Diamonds	16
Cements & Liners	20
Composites & Restoratives	25
Core Materials.....	29
Cosmetic Dentistry	30
Crown & Bridge.....	34
Crown Forms.....	36
Disposables.....	38
Endodontic	40
Equipment.....	48
Evacuation	55
Finishing & Polishing	56
Impression Material	58
Infection Control.....	64
Instruments.....	70
Matrix Materials	76
Pins & Posts.....	77
Preventatives	80
Retraction Materials.....	87
Rubber Dam.....	88
Small Equipment.....	89
Sutures.....	94
X-Ray	96
Miscellaneous	102
Index	103

CONTACT INFO



ADDRESS:

135 North Route 9W
Congers, New York 10920



WEB ADDRESS:

www.Nobledentalsupplies.com



COMPANY CONTACTS :

MAIN PHONE: 866-333-6825



EMAIL ORDERS:

info@nobledentalsupplies.com
orders@nobledentalsupplies.com

TEXT ORDERS OR QUESTIONS:

845-499-7369

FAX: 845-510-3191

Want to increase your bottom line?

Check out our Alternative Line of Products

BUY THIS

MARK 3™ PRODUCTS

INSTEAD OF THIS

COMPARABLE PRODUCTS

Bite Registration 2/bx

Defend Bite, Genie Bite, Blu-Mousse

VPS Putty 300ml. each Base & Catalyst

Express Putty, Genie Putty, Defend Putty

VPS Impression Material 50ml. Cartridges 4/bx

Genie, Defend, Take 1 Advanced, Extrude

VPS Tray Adhesive 10ml

Any other VPS adhesive

Disposable Impression Trays 12/pk

Vista Trays, Any other perforated disposable trays

Pressure Indicator Paste 2.25oz. jar

Mizzy PIP, Sultan PIP

K-Files 6/pk

Kerr K-files, Mani Files, Premier, Maillifer Files

Evacuation Traps 144/bx

Dispos-a-Traps

Evacuation Traps Canister Type

Evac-u-traps

Endo Irrigation Kit 27ga. Yellow or 23ga. Blue 3cc Tips 100/bx

Monoject Endo Irrigation Kits

Quick Tip Air Water Syringe Tips White/Metal Core 100 or 1600/pk

Safe Tips EZ

Multicolored Plastic Air Water Syringe Tips 250/pk

Crystal Tips

Topical Anesthetic Gel 1oz. jar

Topex Topical, Hurracaine, Kolorz, Topicale

Plastic Needles 100/bx

Septoject, Monoject, Accuject

Prophy Angles Disposable

Young Prophy Angles, Denticator

Fluoride Trays Dual Arch 100/pk

Topex Fluoride Trays, Any other dual arch trays

APF Thixotropic Fluoride Gel 17oz. bottle

Topex Fluoride Gel, Enamel Pro Gel, Kolorz Gel

Prophy Paste 200/bx

Topex, Gelato, Glitter, D-Lish, Nupro

Polycarbonate Crowns

3m Polycarbonate Crowns

Hemostat 30cc

Hemodent

Questions? Call 866-333-6825 or Email us at info@nobledentalsupplies.com

FREE SAMPLES AVAILABLE

As good or better than name brands at 50% less

BUY THIS



INSTEAD OF THIS



SILMET PRODUCTS

COMPARABLE PRODUCTS

Profil Composite Syringe & Unidose	TPH3, XRV Ultra, Filtek Supreme, Esthet-X HD
Profil Flowable Composite	Revolution, Filtek Flow
Prolink Bond	Optibond, Prime & Bond, Singlebond, Integrabond
Prolink SE Bond	Optibond All In One, Clearfil S3 Bond, Adper Easy Bond
Spherodon-M	Contour, GS-80
Spherodon +	Tytin, Megalloy, Logic+
Nogama	Dispersalloy
Profil-Temp	Cavit
Proglass One Kit	Fuji One Kit, Rely-X Luting cement
Proglass 9 Kit	Ketac-cem, Fuji IX Kit
Proglass 9 Caps	Fuji IX Caps, Ketac-cem Applicap/Maxicap
Proglass Silver	Miracle Mix, Ketac Silver

CORY LABORATORIES PRODUCTS

ProOne Sonic Tablets	Defend Tabs, Vibrakleen, MaxiTabs, Pro-Portion
Evacshock Tabs	VacuKleen, Vacuum Shock
EvacShock 32oz	Vacusol, Purevac, Neutravac

DemeTECH PRODUCTS

Demetech Sutures	Ethicon, Look, Hu Friedy Sutures
------------------	----------------------------------

Premium Products, Significant Discounts.



HANDPIECE REPAIR

Slow Speed Motor Repairs

Bearing Replacement, Rotor Blades, "O" Ring Replacement, Gears, Pins, Drive Rings, Cleaning, Adjustments, Rethread Backends

High Speed Repairs

Cleaning, Bearing Replacement, "O" Ring Replacement, Balancing, Testing, Clearing Water Lines, Fiber Optic Replacement, Repairing Dented Heads

NoseCone & Angle Repairs

Bearing Replacement, Chuck, Gears, Pins, Springs, Cartridges, Latches & Screws, Dog Drives, Spindle, "O" Rings

What We Repair

All High-speed, Low-speed, Electric and Surgical Handpieces
All Manufacturers (Star, KaVo, NSK, Midwest, including generic and house brands)

Repair Time

We offer a 24-48 hour turnaround time on most repairs.

Cost

Repairs are 30-40% less expensive than sending the unit back to the manufacturer for repair.

Warranty on Repair

All repairs come with a 3-6 month warranty depending on the service needed.

Process

Fill out the repair workorder. Then we will pick up your broken handpiece by sending a UPS call tag. Upon receipt, we will provide an estimate with 24 hours. If the work order is confirmed, the handpiece will be repaired immediately and sent back with 48 hours of initial receipt.

Repair & Replacement Parts

We are fully stocked of replacement parts for most handpieces to ensure repairs are complete as soon as possible.

TOO BUSY TO COMPARE PRICES?

Let us do it for you!

Send us your purchase history or invoices and we will do a price comparison for you.

Email to:

anthonyt@nobledentalsupplies.com
info@nobledentalsupplies.com
quotes@nobledentalsupplies.com

Fax to:

845-510-3191

Text to:

845-499-7369

Receive a quote within 48 hours



We will save you between 30-50%

EXAMPLE:



DENTAL PRODUCT PRICE COMPARISON

Quote for: Summit Dental Arts
31 Lowland Road, Stony Point, NY 10980
Attention: Theresa

All products compared are of the same quality - samples available

Buy this	Our price	Instead of this	Retail price	SAVINGS
100-0220 Bite Registration Super Fast Set 2/bx	\$27.95	Blu-Mousse Classic Standard Cart 50mL 2/Pk	\$45.99	\$18.04
100-0404 Intraoral Mixing Tips Yellow 100/pk	\$15.99	Intra-Oral Tips Small Yellow 100/Bx Kerr	\$54.49	\$38.50
100-1200 Pressure Indicator Paste 2.25oz. Jar	\$17.99	Pressure Indicator Paste 2.25oz Mizzy	\$61.99	\$44.00
100-1202 Multicolored Plastic Air WaterTips 250/pk	\$44.95	Crystal Tip Rainbow 250/pk	\$65.99	\$21.04
100-12106 K-Files 21mm #06 6/pk	\$4.99	K-Files 21mm #06 6/Bx Sybron	\$13.49	\$8.50
100-1249 Xray Film Mounts Universal Roll 1000/bx	\$19.99	Universal Film Mounts Double 1000/Bx HS		\$64.00
100-1402 Headrest Covers Regular 11-1/4" x 10" 250/bx	\$9.99	Headrest Cover 11x9.5" Clear 250/Bx F		\$5.00
100-1411 Quick Tip Air Water Syr. Tips w/Metal Core 1600/bx	\$278.99	Safe-Tips EZ Bulk 1600/Pk Pinnacle		\$118.00
100-1602 Topical Anesthetic Gel Mint 1oz. jar	\$4.99	Topex Topical Gel Mint 1oz/Jr Sulta		\$5.50
100-1612 APF Thixotropic Fluoride Gel Mint 17oz. Bottle	\$10.49	Topex APF Thixo Gel Mint 16oz/Bt		\$16.50
100-1623 Prophy Paste Coarse Mint 200/bx	\$18.49	Topex Prophy Paste Coarse Mint 2		\$21.00
100-16306 Plastic Needles 30ga. Short 100/bx	\$7.99	Septoject Needles 30Ga Short Blue		\$6.00
100-1701 Endo Irrigation Kit 27ga. 3cc Yellow Tips 100/bx	\$22.49	Irrigating Needles Closed End 27Ga		\$47.50
100-2101 Barrier Film Blue 4" x 6" Roll of 1,200/box	\$13.49	Allrap Cover Film 4x6 Clear 1200/RI P		\$13.50
100-2290 PreBent Needle Tips Black 19ga. 100/pk	\$10.19	Needle Tips Pre-Bent 25ga 100/Bx Vista	\$268.99	\$5.10
100-247-144 Prophy Angles Disposable Soft 144/bx	\$40.49	Emerald Prophy Angle LF 500/Bx Denticator	\$40.99	\$228.50
100-3015 VPS Impression Material Heavy FS 50ml. Cart 4/bx	\$29.99	Defend VPS Heavy Body Fast Set 4/Pk Defend	\$165.99	\$11.00
100-385227 Carbide Burs FG 557 100/pk. - Cargus	\$89.00	Midwest Carbide Burs FG 557 100/pk	\$199.49	\$76.99
100-4001 Deluxe Crown & Bridge A2 72 ml. Cartridge	\$37.99	Luxatemp Plus Refill A2 DMG	\$73.99	\$109.50
100-6200 Evacuation Traps Universal 144/pk.	\$29.99	Dispos-A-Trap #6200C 144/Bx	\$55.99	\$36.00
100-7095 Bite Blocks Disposable Fits Panoramic Corp100/bx	\$3.69	Dispos-A-Bite #7095 100/Bx Pinnacle	\$18.99	\$26.00
100-P300 Polycarbonate Crowns #300 5/pk.	\$2.99	Ion Polycarbonate Crowns #300 5/Bx 3M	\$14.79	\$15.30
101-200114 Chlorhexidine Gluconate 0.12% Rinse 16oz.	\$23.99	Chlorhexidine Gluconate Rinse 0.12% 16oz/Bt	\$41.99	\$11.80
	\$23.99	Suture Surg Plain Gut Bge DS18 3-0 18" 12/Bx Look	\$41.99	\$18.00
	\$23.99	Contour Caps Reg Set 2 Spill 500/Jr Kerr	\$1,176.99	\$641.00
	\$23.99		\$1,347.99	\$648.00
	\$23.99			\$135.94

OVER \$2500 SAVED

MARK3™

THE SMART ALTERNATIVE

The MARK3 Line was created to provide
Premium Products
at an **Affordable Price**



Look for the MARK3 line of products throughout the catalog.

PRODUCT EVALUATION



My Dentist was done in 2 minutes!

PROVENEER™


The Ultimate System For Creating LifeLike Veneers

Why Dentists choose ProVeneer:

- Greater acceptance by patient
- No Cad Cam
- No preformed veneers
- No lab bill
- Fun, stress free, impactful

Why patients choose ProVeneer:

- Non-invasive, no filing down of teeth
- Beautiful smiles created in one visit
- Affordable
- Wear and stain resistant
- Looks natural



Each ProVeneer Kit makes at least 10 veneers

Accelerate your revenue with minimum chair time.

The only system with a two minute smile preview veneer.

To learn more visit: WWW.PRO-VENEER.COM

Assorted Kit **\$325.00**

Kit Includes: 10x RBFM Single-use Packs (3x A1, 2x A2, 2x B1, 3x BW), 1 x ProVeneer Bonding Agent, 1 x ProFil Flow Cement, 1 x ProEtch and Accessories

FREE Marketing Materials come with your kit first order:
1 Poster, 50 Brochures, 50 Dental History Forms and Digital Flipbook

Resin Based Foundation Blister Refills **\$225.00** 10/pk

Shades Available: A1, A2, B1, BW (Bleach White)

See page 33 for more details



Featured on DENTALPRODUCTSREPORT.com



ACRYLICS & RELINE MATERIALS

Acrylic Materials

Trim II

Bosworth Co.

230-921090\$98.95

Complete Kit Includes:

6, 1.5 oz. Bottles, 1 each of #59, #61, #62, #65, #69, #77 and 4 oz. Liquid

230-921091\$25.49 Liquid - 4 oz.

230-921092Powder, #59.

230-921093Powder, #61

230-921094Powder, #62

230-921095Powder, #65

230-921096Powder, #69

230-921097Powder, #77

\$20.95 1.5 oz.



Alike

GC America

Liquid

250-340591\$27.49 4 oz./btl

A quick-setting, self-cured temporary crown and bridge resin

Powder

250-340512#62 - Shade A1

250-340513#65 - Shade A2

250-340514#67 - Shade A3

250-340516#81 - Shade A3.5

250-340511#59 - Shade B1

250-340515#77 - Shade B3

\$27.49 45g/btl



Cargus

Glass Dappen Dish

100-2615.....Light Blue

100-2616.....Clear

\$3.95

Durable, deep dish can hold acrylic liquid, powder and other liquids



MARK3™

Disposable Dappen Dishes

100-0036.....Blue

100-0037.....White

100-0038.....Yellow

100-0039.....Pink

\$13.95 250/pk

100-0040.....Assorted \$35.99 1000/bx

- Disposable plastic dappen dish is easy-to-use
- Provides convenience when applying polish
- Comes in assorted colors



Jet Tooth Shade Acrylic

Lang

Liquid

455-1403.....\$17.95
4 oz./btl

455-1404.....\$29.95
8 oz./btl

Powder

455-1410-59Shade B1

455-1410-60Shade A1

455-1410-61Shade B2

455-1410-62Shade A2

455-1410-65Shade D3

455-1410-66Shade A3

455-1410-67Shade B3

455-1410-69Shade C2

455-1410-77Shade B4

455-1410-81Shade A3.5

\$17.95 2 oz.



With orders over \$300
Excludes Canada. Certain restrictions apply.

Tissue Conditioner

Coe Comfort

GC America

250-341001\$134.95

Professional Pack

Contains:
6 oz. Powder, 6 oz. Liquid, Lubricant,
and Accessories



Reline Materials

Truliner Original

Bosworth Co.

230-0921940\$96.95

Standard Pink - Hard Reline

Contains:
8 oz. Powder, 6 oz. Liquid
and Accessories



New Truliner

Bosworth Co.

230-0921970\$96.95

Standard Pink - Hard Reline

Contains:
1 - 8 oz. Powder, 1 - 6 oz. Liquid,
1.33 oz. Bonding Liquid,
20 Applicator Brushes, 2 Mixing Cups,
1 Powder Nozzle, 1 Powder Measure
and 1 Liquid Measure



Softone

Bosworth Co.

230-0921775\$74.95

Standard White Kit - Soft Reline

Contains:
3 oz. Powder, 4 oz. Liquid



Coe Pak

GC America

250-135001Regular

250-135301Hard & Fast

\$83.99

Each Package Contains:
90g Tube each of Base
and Catalyst



Coe Pak Automix NDS

GC America

250-135003\$95.99

Two-Pack Contains:
2 - 50ml Cartridges,
and 12 Mixing Tips



Coe Soft

GC America

250-344001\$137.50

Professional Pack - Soft Reline

Contains:
6 oz. Powder, 6 oz. Liquid, Lubricant,
and Accessories



Kooliner

GC America

250-345001\$99.95

Professional Pack - Hard Reline

Contains:
3 oz. Powder, 2 oz. Liquid and
Accessories



CAN'T FIND IT?

We're always adding new products.
If you don't see it here, visit our
website, give us a call or email us at:
info@nobledentalsupplies.com

Amalgam



Spherodon-M® Admix Alloy 45.5% Ag (Silver), High Copper, Non Gamma II

A spherical lathe cut blended amalgam alloy (type III) that provides excellent stability, enhanced corrosion resistance and reduced ditching.

FEATURES

- Outstanding dimensional stability
- Improved marginal integrity
- More body & easy condensation versatility
- Reduced risk of fracture
- High corrosion resistance

PHYSICAL PROPERTIES

- Composition Ag 45.5%, Sn 31.5%, Cu 23%
- Particle shape 40% Lathe cut 60% spherical
- Average Creep (7 days) 0.35%
- Average Dimensional Change +0.08%
- Average Comprehensive Strength (1 hr) 178 MPa
- Average Comprehensive Strength (24 hr) 460 MPa



500 per bottle

- | | | |
|-------------------|---------------------|-----------------|
| 104-120491 | 1 Spill Fast Set | \$425.00 |
| 104-120492 | 2 Spill Fast Set | \$595.00 |
| 104-120493 | 3 Spill Fast Set | \$699.00 |
| 104-121491 | 1 Spill Regular Set | \$425.00 |
| 104-121492 | 2 Spill Regular Set | \$595.00 |
| 104-121493 | 3 Spill Regular Set | \$699.00 |

50 per bottle

- | | | |
|-------------------|---------------------|----------------|
| 104-120471 | 1 Spill Fast Set | \$46.95 |
| 104-120472 | 2 Spill Fast Set | \$64.95 |
| 104-120473 | 3 Spill Fast Set | \$74.50 |
| 104-121471 | 1 Spill Regular Set | \$46.95 |
| 104-121472 | 2 Spill Regular Set | \$64.95 |
| 104-121473 | 3 Spill Regular Set | \$74.50 |

Spherodon+® Spherical Alloy 60% Ag (Silver), High Copper, Non Gamma II

A single composition, ultra-fine spherical ternary. This formulation utilizes modern resin technology that offers unique dimensional stability and resistance to corrosion. Exhibits high comprehensive strength due to its unique formulation. This is achieved by spraying inert nitrogen under high pressure to produce a fine spherical powder.

FEATURES

- 100% spherical particles
- Decreases micro leakage
- Superb compressive & tensile strength
- Easy handling & smooth carving
- Extra low mercury alloy ratio

PHYSICAL PROPERTIES

- Composition Ag 60%, Sn 28%, Cu 12%
- Particle shape 100% spherical
- Average Creep (7 days) 0.62%
- Average Dimensional Change +0.01%
- Average Comprehensive Strength (1 hr) 171 MPa
- Average Comprehensive Strength (24 hr) 485 MPa



500 per bottle

- | | | |
|-------------------|---------------------|------------------|
| 104-130491 | 1 Spill Fast Set | \$575.00 |
| 104-130492 | 2 Spill Fast Set | \$825.00 |
| 104-130493 | 3 Spill Fast Set | \$1050.00 |
| 104-131491 | 1 Spill Regular Set | \$575.00 |
| 104-131492 | 2 Spill Regular Set | \$825.00 |
| 104-131493 | 3 Spill Regular Set | \$1050.00 |

50 per bottle

- | | | |
|-------------------|---------------------|-----------------|
| 104-130471 | 1 Spill Fast Set | \$65.95 |
| 104-130472 | 2 Spill Fast Set | \$94.20 |
| 104-130473 | 3 Spill Fast Set | \$119.50 |
| 104-131471 | 1 Spill Regular Set | \$65.95 |
| 104-131472 | 2 Spill Regular Set | \$94.20 |
| 104-131473 | 3 Spill Regular Set | \$119.50 |

Nogama® Dispersed Phase Alloy 69.5% Ag (Silver), High Copper, Non Gamma II

An admix of spherical, copper-silver eutectic particles and lathe-cut conventional alloy. It has a high compressive strength that minimizes the risk of restoration fracture.

FEATURES

- Provides a firm pack & smooth carving
- Increased strength of restored teeth
- Supreme marginal adaptability
- Corrosion resistant
- Excellent durability

PHYSICAL PROPERTIES

- Composition Ag 69.5%, Sn 19.5%, Cu 10.5%, Zn 0.5%
- Particle shape 80% Lathe cut 20% spherical
- Average Creep (7 days) 0.35%
- Average Dimensional Change +0.08%
- Average Comprehensive Strength (1 hr) 164 MPa
- Average Comprehensive Strength (24 hr) 470 MPa



500 per bottle

- | | | |
|-------------------|---------------------|------------------|
| 104-110491 | 1 Spill Fast Set | \$599.00 |
| 104-110492 | 2 Spill Fast Set | \$849.00 |
| 104-110493 | 3 Spill Fast Set | \$1095.00 |
| 104-111491 | 1 Spill Regular Set | \$599.00 |
| 104-111492 | 2 Spill Regular Set | \$849.00 |
| 104-111493 | 3 Spill Regular Set | \$1095.00 |

50 per bottle


- | | | |
|-------------------|---------------------|-----------------|
| 104-110471 | 1 Spill Fast Set | \$62.95 |
| 104-110472 | 2 Spill Fast Set | \$89.95 |
| 104-110473 | 3 Spill Fast Set | \$114.95 |
| 104-111471 | 1 Spill Regular Set | \$62.95 |
| 104-111472 | 2 Spill Regular Set | \$89.95 |
| 104-111473 | 3 Spill Regular Set | \$114.95 |

Amalgam Separators

**Amalsed Conversion Kit
Medentex**

910-A0US0002\$499.00

- Perfect fit for pre-existing Air Techniques Guardian systems
- ISO 11143 certified at over 99% efficiency



**Amalsed Direct Dry Vacuum
Amalgam Separator
Medentex**

910-A0US1803Small \$799.00

For use in offices with four operatories/ 1-2 Dentists. Includes one amalsed green cartridge.



910-A0US1804Large \$899.00

For use in offices with eight operatories/3-4 Dentists. Includes two amalsed green cartridges

- The amalsed® direct wet kit is a perfect fit for all dental offices
- Wet ring vacuum system
- ISO 11143 certified at over 99% efficiency

**Amalsed Direct Wet Vacuum
Amalgam Separator
Medentex**

910-A0US1801\$849.00

For use in offices with four operatories/ 1-2 Dentists. Includes one amalsed green cartridge.



910-A0US1802\$949.00

For use in offices with eight operatories/3-4 Dentists. Includes two amalsed green cartridges

- The amalsed® direct dry kit is a perfect fit for all dental offices
- Dry vacuum system
- ISO 11143 certified at over 99% efficiency

Amalgam Separator Waste Containers

**Medentotainer Amalgam
Waste Container
Medentex**

**910-B0US1901Small - 1.5 Gallon
\$149.00**

**910-B0US1902Large - 4.0 Gallon
\$249.00**

**910-B0US1904X-Large - 6.5 Gallon
\$349.00**

Includes:
Container, Packaging and Prepaid Return Shipping Label

UN/DOT approved compliant containers for the disposal of amalgam waste from chair-side traps, filters, leftover or spent amalgam, capsules and amalgam containing teeth.



**Medentotainer Cap Guard
Medentex**

910-B0US19051.5 Gallon \$189.00

Includes:
Container, Packaging and Prepaid Return Shipping Label

UN/DOT approved compliant waste container for empty, left-over and expired amalgam capsules. Includes container, packaging and prepaid return shipping label.



**Medentotainer Lead Foils
Medentex**

910-B0US19031.5 Gallon \$119.00

Includes:
Container, Packaging and Prepaid Return Shipping Label

UN/DOT approved compliant waste container for x-ray lead foils.



**Medentotainer Tooth Box
Medentex**

910-B0US19000.35 Gallon \$79.00

Includes:
Container, Packaging and Prepaid Return Shipping Label

UN/DOT approved compliant waste container for removed teeth containing amalgam.



Manage your costs
without compromising
quality with a leading
imported brand.

Silmet Ltd. specializes in the development and production of high quality dental restorative materials. With years of experience and successful clinical trials, our products are used by dentists worldwide. Silmet is known as a source for private label and OEM products for many dental manufacturers.

Silmet Ltd. is ISO certified. All products meet the European Medical Device Directive and are 510K registered with the FDA.



Profil nano hybrid is more advanced than XRV and Z-250

ProFil

UNIVERSAL MICROHYBRID COMPOSITE with
ADVANCED NANO TECHNOLOGY

4g Syringe **\$39.95**

Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D3, OA2

Compare to Filtek Supreme Syringes

Compules .315g ea. - 20/pk **\$53.95**

Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, BW, C2, C3, D3, OA2

Compare to Filtek Supreme Unidose

ProFil Flow

FLOWABLE MICROHYBRID COMPOSITE

4, 1g Syringes plus 10 Tips **\$41.95**

Shades: A1, A2, A3, A3.5, B1, B2, B3, C2,

Compare to Revolution Formula2



FIFTH GENERATION



ProLink

SINGLE STEP BONDING AGENT

5 ml Bottle **\$46.95**

Compare to Prime & Bond NT Refill

ProLink Self Etch

SELF ETCH BONDING AGENT

4ml Bottle, 25 Disposable Applicators, and
1 Dispensing Dish **\$82.95**

Compare to Optibond All-in-One

Dental Products Report
Issue 10/2010, Volume 44



PRODUCT REFERENCE GUIDE

Compare these Silmet products to the more expensive name brands

Description	Compare To:				
	Dentsply	Kerr	3M	GC	SDI
ProFil Composite Syringe & Unidose	Esthet-X HD	XRV Ultra	Filtek Supreme		
Profil Flowable Composite		Revolution	Filtek Flow		
Prolink Bond	Prime & Bond	Optibond	Singlebond		
Prolink SE Bond		Optibond All In One	Adper Easy Bond		
Spherodon-M		Contour			GS-80
Spherodon +	Megalloy	Tytin			Logic+
Nogama	Dispersalloy				
Profil-Temp			Cavit		
Proglass One Kit				Fuji One Kit	
Proglass Two Light Cure Caps				Fuji II	
Proglass 9 Kit			Ketac-cem	Fuji IX Kit	
Proglass 9 Caps			Ketac-cem Caps	Fuji IX Caps	
Proglass Silver Kit			Ketac Silver	Miracle Mix	
Proglass Plus Kit				Fuji Plus	
Proglass Plus Caps				Fuji Plus	



Anesthetic Cartridges

Lidocaine HCL 2% with Epinephrine

Novocol

101-0100 **\$28.95** 50/bx
1:100,000



**ONLY
\$28.95**

- This versatile, widely applicable injection delivers local anesthesia of intermediate duration (pulpal anesthesia of at least 60 minutes and soft tissue anesthetic of approximately 2.5 hours)
- Now packaged in easily stackable boxes
- Convenient blister packs protect ampoules from light and air

Lignospan (Lidocaine) Forte 2% with Epi Septodont

290-01-A1050 **\$37.95** 50/bx
2% w/ Epinephrine 1:500,000



Vivacaine HCL 0.5% with Epinephrine Septodont

290-01-A1600 **\$59.95** 50/bx
1:200,000



Mepivacaine HCL 3% Plain Novocol

101-MP3 **\$39.95** 50/bx



Septocaine HCL 4% with Epinephrine Septodont

290-01-A1400 **\$46.95** 50/bx
1:100,000




290-01-A1200 **\$46.95** 50/bx
1:200,000

Aspirating Syringes

J&J Instruments

6444 Type A **\$16.75**

6333 Type C-W **\$16.75**



Anesthetic Cartridge Warmer & Dispenser

Premier Dental

580-1048021 **\$64.99**



Needles

Monoject Needles

Covidien/Kendall

- 520-40003325G Long Plastic
 - 520-40005827G Long Plastic
 - 520-40006627G Short Plastic
 - 520-40007430G Short Plastic
 - 520-40017330G X-Short Plastic
 - 520-40103125G Long Metal
 - 520-40104925G Short Metal
 - 520-40105627G Long Metal
 - 520-40106427G Short Metal
 - 520-40107230G Short Metal
 - 520-40117130G X-Short Metal
- \$12.99** 100/bx



MARK 3™

BEST VALUE

Plastic Hub Needles

- 100-1630327G Short
 - 100-1630427G Long
 - 100-1630530G X-Short
 - 100-1630630G Short
- \$9.95** 100/bx



- Clean puncture into the septum without punching out rubber particles due to the accurate grinding
- Bevel mark for optimal injection position
- High quality disposable needles

Septoject Needles (Plastic Hub)

Septodont

- 290-01-N127227G Long
- 290-01-N130030G X-Short
- 290-01-N127127G Short
- 290-01-N130130G Short

\$13.99 100/bx



Topicals

MARK3™

BEST VALUE

Anesthetic Gel

- 100-1601.....Cherry
- 100-1602.....Mint
- 100-1603.....Strawberry
- 100-1604.....Pina Colada
- 100-1605.....Bubble Gum
- 100-1606.....Grape



\$4.99 1 oz./jar

COMPARE TO
**Sultan
Topex**



- 20% Benzocaine
- Fast acting Topical Anesthetic
- Can also be used to comfort patients with ulcers, wounds and other minor mouth irritations
- Eliminates the gag reflex when performing radiographs, periodontal, and prosthetic impressions

Topicale Gel Pump

Premier Dental

- 580-9007140**\$39.99** Assorted 3/pk
Assorted Package Contains:
3 - 43g pumps (1 each) of Cherry, Raspberry, Piña Colada

- 580-9007141Cherry
- 580-9007142Raspberry
- 580-9007143Strawberry
- 580-9007144Piña Colada



\$13.50 43g - 1/btl

Topicale Gel Patch

Premier Dental

- 580-9007149**\$52.50** 25/bx



Topicale Xtra

Premier Dental

- 580-9007152Very Berry Cherry
- 580-9007153Mint Freeze
- 580-9007154Jammin' Strawberry



\$7.49 - 1 oz. jar

- 580-9007155**\$20.95** Assorted 3/pk
Contains: 1 oz. each Cherry, Mint, Strawberry

Articulating Paper



Bausch

Bausch

200 Microns .008"

- 610-BK-01 - Blue - Pre Cut Strips.....**\$17.95** 300/bx
- 610-BK-02 - Red - Pre Cut Strips.....**\$17.95** 300/bx
- 610-BK-03 - Blue - Horseshoe.....**\$15.25** 50/pk
- 610-BK-05 - Blue - Booklet**\$15.49** 300/pk
- 610-BK-28 - Black/Red - Arti-Fol**\$16.95** 20m/roll

40 Microns .0016" Micron Thin

- 610-BK-09 - Blue - Strips**\$14.95** 200/pk
- 610-BK-10 - Red - Strips**\$14.95** 200/pk
- 610-BK-61 - Blue - Pre Cut Strips.....**\$9.25** 200/bx
- 610-BK-62 - Red - Pre Cut Strips.....**\$9.25** 200/bx

100 Microns Progress 100

- 610-BK-53 - Blue - Horseshoe**\$15.49** 50/pk
- 610-BK-54 - Red - Horseshoe**\$15.49** 50/pk

Occlusal & Pressure Indicators

MARK3™

Pressure Indicator Paste

- 100-1200.....**\$22.95**
2.25 oz./jar



Hydent Indicator Spray

Pascal

- 560-05-100.....**\$16.95** 30g/can



Occlude Aerosol Powder

Pascal


- 560-05-300.....Green
 - 560-05-320.....Red
- \$17.49** 23g/can




100 per pack **\$99.95**

Cargus CARBIDE BURS


Round

Item#	Bur Size	Shank Type	
100-386201	1/4	FG	
100-386202	1/2	FG	
100-385204	2	FG	
100-386268	3	FG	
100-385205	4	FG	
100-385206	6	FG	
100-385207	8	FG	
100-386101	1/4	RA	
100-386102	1/2	RA	
100-386103	1	RA	
100-386104	2	RA	
100-386169	5	RA	
100-385105	4	RA	
100-385106	6	RA	
100-386107	8	RA	
100-386304	2	HP	
100-386305	4	HP	
100-386306	6	HP	
100-386307	8	HP	


Taper Fissure

Item#	Bur Size	Shank Type	
100-386221	169	FG	
100-386223	171	FG	


Taper Fissure Crosscut

Item#	Bur Size	Shank Type	
100-386232	700	FG	
100-385233	701	FG	
100-386234	702	FG	


Straight Dome

Item#	Bur Size	Shank Type	
100-385255	1156	FG	
100-385256	1157	FG	
100-385257	1158	FG	


Straight Dome Crosscut

Item#	Bur Size	Shank Type	
100-386259	1557	FG	
100-386456	1557SS	FG	
100-385260	1558	FG	


Amalgam Prep

Item#	Bur Size	Shank Type	
100-385371	245	FG	
100-386471	245SS	FG	


Inverted Cone

Item#	Bur Size	Shank Type	
100-385209	33 1/2	FG	
100-386210	34	FG	
100-385211	35	FG	


Pear Shape

Item#	Bur Size	Shank Type	
100-385261	330	FG	
100-385461	330SS	FG	
100-386262	331	FG	
100-385264	331L	FG	
100-386263	332	FG	

Straight Fissure

Item#	Bur Size	Shank Type	
100-385216	56	FG	
100-386217	57	FG	

Straight Fissure Crosscut

Item#	Bur Size	Shank Type	
100-385226	556	FG	
100-386426	556SS	FG	
100-385227	557	FG	
100-385427	557SS	FG	
100-385228	558	FG	

Surgical

Item#	Bur Size	Shank Type	100 per pack \$179.95
100-385605	4	FGOS	
100-385606	6	FGOS	
100-385627	557	FGOS	
100-385628	558	FGOS	
100-399633	701	FGOS	
100-399634	702	FGOS	

Trimming & Finishing

Item#	Bur Size	Shank Type	100 per pack \$499.99
100-386522	7404	FG	
100-386510	7901	FG	
100-385523	7406	FG	
100-385524	7408	FG	
100-385509	7803	FG	
100-385511	7902	FG	
100-383512	7903	FG	

Compare to Midwest

10 per pack

Cargus CARBIDE BURS

Round

Item#	Bur Size	Shank Type	Price/10pk
100-389104	2	RA	\$12.95
100-389105	4	RA	\$12.95
100-389106	6	RA	\$12.95
100-389107	8	RA	\$12.95
100-389201	1/4	FG	\$12.95
100-389202	1/2	FG	\$12.95
100-389204	2	FG	\$12.95
100-389205	4	FG	\$12.95
100-389206	6	FG	\$12.95
100-389207	8	FG	\$12.95
100-389301	1/4	HP	\$12.95
100-389302	1/2	HP	\$12.95
100-389304	2	HP	\$12.95
100-389305	4	HP	\$12.95
100-389306	6	HP	\$12.95
100-389355	5	HP	\$12.95

Taper Fissure Crosscut

Item#	Bur Size	Shank Type	Price/10pk
10-389233	701	FG	\$12.95

Straight Dome

Item#	Bur Size	Shank Type	Price/10pk
100-389256	1157	FG	\$12.95
100-389257	1158	FG	\$12.95

Straight Dome Crosscut

Item#	Bur Size	Shank Type	Price/10pk
100-389260	1558	FG	\$12.95
100-389360	1558	HP	\$12.95

Amalgam Prep

Item#	Bur Size	Shank Type	Price/10pk
100-389371	245	FG	\$12.95

Inverted Cone

Item#	Bur Size	Shank Type	Price/10pk
100-389110	34	RA	\$12.95
100-389209	33 1/2	FG	\$12.95
100-389210	34	FG	\$12.95
100-389211	35	FG	\$12.95

Pear Shape

Item#	Bur Size	Shank Type	Price/10pk
100-389261	330	FG	\$12.95

Straight Fissure Crosscut

Item#	Bur Size	Shank Type	Price/10pk
100-389226	556	FG	\$12.95
100-389227	557	FG	\$12.95
100-389427	557SS	FG	\$12.95
100-389228	558	FG	\$12.95

Surgical

Item#	Bur Size	Shank Type	Price/10pk
100-399602	2	FGOS	\$18.49
100-399604	4	FGOS	\$18.49
100-399606	6	FGOS	\$18.49
100-399607	8	FGOS	\$18.49
100-399627	557	FGOS	\$18.49
100-389628	558	FGOS	\$18.49
100-389633	701	FGOS	\$18.49

Trimming & Finishing

Item#	Bur Size	Shank Type	Price/10pk
100-389510	7901	FG	\$49.99
100-389511	7902	FG	\$49.99
100-389513	7104	FG	\$49.99
100-389522	7404	FG	\$49.99
100-389523	7406	FG	\$49.99
100-389524	7408	FG	\$49.99

GALIL DIAMOND BURS - Multi Use

10 per pack **\$31.49**

	BALL				INVERTED CONE				TAPER ROUND END			FRICTION GRIP								
Shaded spaces below illustrations are the grits available																				
Catalog Number	A-1	A-2	A-3	A-4	B-1	B-2	B-3	B-4	C-1	C-2	C-3	K-1	K-2	L-1	M-1	M-2	M-3	M-4	M-5	M-6
ISO Number	801-012	801-014	801-018	801-022	805-012	805-018	807-012	807-018	855-012	805-016	807-025	368-016	368-023	431-020	813-016	813-018	855-018	811-033	813-012	813-014
Extra Fine Grit (White)																				
Fine Grit (Yellow)																				
Medium Grit (Blue)																				
Coarse Grit (Green)																				

	FLAME								CYLINDER				TAPER FLAT END ROUND END				CYLINDER			
Shaded spaces below illustrations are the grits available																				
Catalog Number	E-1	E-2	E-3	E-4	E-5	E-6	E-7	E-8	D-1	D-2	D-3	D-4	F-1	F-2	F1R	F2R	G-1	G-2	G-3	G-4
ISO Number	852-012	858-014	859-016	860-014	862-014	862-016	885-012	860-016	835-012	835-016	836-014	837-014	848-016	848-018	850-016	850-018	855-018	811-033	813-012	813-014
Extra Fine Grit (White)																				
Fine Grit (Yellow)																				
Medium Grit (Blue)																				
Coarse Grit (Green)																				

	GINGIVAL CURETTAGE SET											
These shapes are available in medium grit only												
Catalog Number	S-1	S-2	S-3	S-4	S-5	T-1	T-2	T-3	T-4	T-5	T-6	
ISO Number	879k-014	878k-012	879k-012	877k-012	878k-014	878k-020	878k-022	878k-018	877k-018	877k-014	879k-018	

	FOR GROSS REDUCTION					
These shapes are available in super-coarse grit only						
Catalog Number	FP2	CP2	EP5	EP6	KP2	MP3
ISO Number	848-018	805-016	862-014	862-016	368-023	855-018

Accessories

MARK 3™

Bur Blocks

100-19996.....Magnetic - 24 holes

\$8.99 1/pk

100-19997.....FG - 15 holes

\$12.59 1/pk

100-19999.....FG & RA - 30 holes

\$15.99 1/pk



Clean-A-Diamond

Premier Dental

580-2014001.....Stone - **\$16.95**

580-2014015.....Mini Squares - **\$23.95** 12/bx



DID YOU KNOW?

GALIL diamond burs were tested at the most prestigious dental testing facility in the United States, CRA.

They were proven to be better than most well known, widely used brand name diamond burs. (Brasseler, NTI, Two Striper)

In lab testing, using "macor," a substance stronger than glass, Galil had the fastest cut and the best longevity.



DIAMOND BURS

5 per pack FROM **\$49.50**



Please Check Our Website
For many more Diamond types and sizes plus sizing charts

Round End Taper

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015459	770.10	coarse	1.1/1.8	856L-018
580-2015457	770.8	coarse	1.3/1.8	856-018
580-2015467	780.9	coarse	1.4/2.1	856L-021
580-2015455	767.9	coarse	0.9/1.8	36006-018
580-2015453	767.7	coarse	1.1/1.8	855-018
580-2015475	799.11	coarse	1.3/1.8	856-018
580-2015473	799.6	coarse	0.8/1.5	869-014
580-2015469	785.7	coarse	1.5/2.4	KS6

Flat End Taper

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015345	703.8	coarse	1.2/1.8	847-018
580-2015343	702.8	coarse	0.9/1.4	847-014
580-2015341	701.9	coarse	0.9/1.8	848-018
580-2015339	701.7	coarse	1.1/1.8	846-018
580-2015347	703.10	coarse	1.1/1.8	848-018
580-2015335	700.9	coarse	0.6/1.4	848-014
580-2015337	700.11	coarse	0.7/1.8	848-018

Bevel End Cylinder

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015067	250.8	coarse	0.4/1.2	885-010
580-2015069	251.8	coarse	0.4/1.2	885-012
580-2015087	250.8	fine	0.4/1.2	885-010

Flat End Cylinder

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015277	515.7	coarse	1.2	837-012
580-2015275	515.5	coarse	1.2	836-012

Football

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015143	285.5	coarse	0.7/2.3	368-023
580-2015173	285.5	fine	0.7/2.3	368-023
580-2015203	285.5	very fine	0.7/2.3	368-023
580-2015142	283.4	coarse	2.5	368-025
580-2015145	292.3	coarse	0.4/1.6	368-016

Crown Cut

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015575	SC8	coarse	1.7	
580-2015583	ST8	coarse	1.4/1.8	

Round

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015007	135	coarse	2.6	801-025
580-2015008	130	coarse	2.3	801-021

Flame

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015123	260.8	coarse	0.4/1.4	862-014
580-2015131	265.8	fine	0.4/1.4	862-014
580-2015153	260.8	fine	0.4/1.4	862-014
580-2015125	260.10	coarse	0.5/1.4	863-014
580-2015183	260.8	very fine	0.4/1.4	862-014
580-2015185	260.10	very fine	0.5/1.4	863-014
580-2015155	260.10	fine	0.5/1.4	863-014
580-2015139	270.9	coarse	0.5/1.9	862-021
580-2015161	265.8	fine	0.4/1.4	30011-54
580-2015124	261.8	coarse	0.4/1.1	862-012
580-2015071	253.10	coarse	1.0/1.6	879K-016

Barrel

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015063	230	coarse	1.3/3.4	811-033

Interproximal Trimmer

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015043	201.3	fine	0.3/0.8	889-010

Wheel

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015545	863	coarse	4.7	909-050

Gold Plated Curettage

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015607	252S	medium	0.4/1.8	879K-014
580-2015605	252SB	medium	0.4/1.4	878K-014
580-2015611	253SB	medium	0.4/1.6	878K-016
580-2015617	254SB	medium	0.4/1.8	878K-018

Cone

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015037	209.10	coarse	0.5/2.3	859-021

Amalgam Removers

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2015599	1250	medium	0.9/1.3	806-014

Amalgam Removers

Item#	Type	Grit	Maj/Min Diam	Ref. #
580-2013281	OP 557	medium	0.9/1.0	H31-010

CEMENTS & LINERS

Cavity Liners & Bases

Fuji Lining LC Paste Pak

GC America

250-001887\$119.95

Refill Pack

Contains:

1 - 7g Cartridge and 1 Mixing Pad



Cavity Liner Liquid

Pulpdent

590-PCL\$23.49 15cc/btl



Lime-Light

Pulpdent

590-LIME\$41.95

Cavity Liner Syringe Kit

Contains: 4 - 1.2ml Syringes and 8 Pre-bent Tips

590-LIME-3\$24.95 Refill Syringe 3ml



Multi-Cal

Pulpdent

590-MULTI\$42.49 4/pk

Calcium Hydroxide Syringes

Contains: 4 - 1.2ml Syringes and Applicator Tips



Vitrebond

3M/ESPE

760-7510\$101.00

Intro Kit Contains:

9g Powder, 5.5ml Liquid and Mixing Pad



760-7512L\$69.99 Liquid - 5.5ml/btl

Vitrebond Plus Clicker

3M/ESPE

760-7572\$259.95

Double Pack Contains:

2 - 20g Liner/Base Clickers



SILMET

ProBase Liner/Base

104-170001.....\$59.95

Powder/16g and
Liquid/12ml

Based on its excellent adhesion qualities, ProBase™ can be used under any kind of restorative material. ProBase™ exhibits high radiopacity, high strength and snap set.



COMPARE TO
Vitrebond
Intro Kit

Desensitizers

Dentin Desensitizer

Pulpdent

590-DES\$39.99 12ml/btl



MARK3™

BEST VALUE

Dentin Desensitizer

100-4500\$34.99 10ml/btl

- Contains benzalkonium chloride, HEMA and fluoride as an added source of fluoride ion
- Does NOT contain glutaraldehyde, a toxic irritant
- Superior relief from sensitivity



COMPARE TO
Gluma

Glass Ionomer Cements

Fuji Cem 2 Automix Paste Pak Slide & Lock System

GC America

E250-004603\$179.99

Refill Pack Contains:

2 - 13.3g Cartridges and Mixing Pad

No Tips



FujiCem 2 Mixing Tips

GC America

250-431578White Slide & Lock

\$31.99 20/bg



Fuji II LC Capsules

GC America

250-425001Shade A1

250-425002Shade A2

250-425003Shade A3

\$224.95 48/bx



Fuji IX GP Caps Fast Set

GC America

250-425060Shade A1

250-425061Shade A2

250-425062Shade A3

\$224.95 48/bx



Fuji I

GC America

E250-002503**\$69.95** 1:1 Kit

Contains:

35g Powder, 20g Liquid, Spatula and Mixing Pad

250-425095**\$154.99** Capsules - 48/bx



Fuji Plus Capsules

GC America

250-425030**\$149.95**

A3 Yellow

48/bx



Ketac Cem Aplicap

3M/ESPE

E760-56060**\$184.95**

Refill Pack

50/bx



SILMET

BEST VALUE

ProGlass One Kit

104-171001.....**\$34.99**

Powder/17g and
Liquid/7ml

Chemically bonds
to tooth structure
& metal providing
excellent strength
& marginal integrity
for long term restorations.



COMPARE TO
Fuji I

ProGlass Plus Kit

104-17P001**\$79.50**

Powder/15g and
Liquid/9.6ml

Higher bond strength
than conventional
glass ionomer
cements.



COMPARE TO
Fuji Plus

ProGlass Plus Capsules

104-17PA00R**\$179.95**

50/pk

Simple placement
technique with low
coefficient of
thermal expansion.

Excellent biocompatibility
to tooth structure and
soft tissues.



COMPARE TO
Fuji Plus

ProGlass Nine Capsules

104-179A02R**\$149.95**

Shade A2

50/pk

High translucency
ensures use in
anterior restorations
without compromising
aesthetics.

*The ONLY brand
with removable tips!*



COMPARE TO
Fuji IX Caps

Ketac Cem Radiopaque

3M/ESPE

E760-37200.....\$57.95

Intro Package Contains:
12ml Liquid and 33g Powder

E760-37230.....\$149.95

Triple Pack Contains:
3 - 12ml Liquid, 3 - 33g Powder,

E760-37220.....\$33.99

Liquid Only 12ml/btl

E760-37210.....\$69.99

Powder Only 33g/btl



RelyX Luting Cement

3M/ESPE

E760-3515.....\$86.95

Complete Kit Contains:
16g Powder, 9ml Liquid,
2 Mixing Pads and Scoop



RelyX Luting Plus

3M/ESPE

E760-3525.....\$145.95

Refill Pack Contains:
2 - 11.5g Clicker Dispensers



Permanent Cements

Flexi-Flow Auto E

Essential Dental Systems

235-880-00.....\$122.95

Complete Kit Contains:
2 Automix Syringes (each contains
1 - 4.5g Base and 1 - 4.5g Catalyst),
and 20 Intra Oral Tips (Self Cure)



Flexi-Flow Cem

Essential Dental Systems

235-850-00.....\$107.95

Contains:
9g of Base Gel, 9g of Catalyst Gel,
Mixing Pad and Plastic Spatula



Flexi-Flow Natural

Essential Dental Systems

235-860-00.....\$107.95

Contains:
9g of Base, 9g of Catalyst, Mixing Pad
and 100 Plastic Spatulas



MARK 3™

BEST VALUE

Permanent Resin Cement

Self Adhesive - Self Curing

100-5380.....\$49.99

7ml Automix Syringe and 6 Mixing Tips

- Versatile solution to many luting indications
- Biocompatibility minimizes post-operative discomfort
- Insoluble in oral fluids
- High bond strength
- Non-Eugenol



COMPARE TO
Maxcem
Elite

Maxcem Elite Clear

Kerr

430-33872.....\$119.95

Clear Syringe Refills
Contains:
2 - 5g Dual Syringes,
20 - 4:1 Automix Tips
(10 regular & 10 wide),
10 Intraoral Tips and
10 Post-canal Tips



IntegraCem

Premier Dental

580-3001410\$112.99

Contains:
4g of each (Catalyst and Base)



Implant Cement

Premier Dental

580-3001450\$**59.95**

Contents:

1 - 5ml Automix Dual-Barrel Syringe,
10 Mixing Tips, 1 Mixing Pad



Embrace Wetbond Resin Cement

Pulpdent

Syringe Kit\$**98.95**

590-EMCA.....Low Viscosity

590-EMCM.....Medium Viscosity

Contains: 7g Automix Syringe, 20 Automix Tips,
1.2ml Seal-n-Shine and 1.2ml OxOut

Syringe Refill\$**72.49**

590-EMCARLow Viscosity

590-EMCMRMedium Viscosity

Contains: 7g Automix Syringe and 20 Automix Tips



Multilink

Vivadent

Primer Bottle

820-613626Refill A - \$**46.99** 3g

820-613627Refill B - \$**46.99** 3g

820-576825.....Refill A&B - \$**79.95** 3g each

Automix NG Refill

820-615216Translucent

820-615217.....Yellow

820-615218Opaque

\$**159.95** 3g each



Durelon

3M/ESPE

E760-38236\$**49.95**

Powder - 60g jar

E760-38216.....\$**49.95**

Liquid - 40ml



RelyX Unicem Aplicap

3M/ESPE

E760-56818A2 Universal

\$**299.95** 50/bx



Temporary Cements

MARK 3™

BEST VALUE

Temporary Cement NE (Non-Eugenol)

100-4125.....\$**49.95**

Contains: 2 - 5ml (8.6gm) Automix Syringe
plus 20 Mixing Tips

- Easy and comfortable to use
- Suitable for cementing temporary crowns and bridges as well as trial cementing of permanent restorations
- Excellent adhesion to temporary crown and bridge materials
- Easily removed from prepared tooth stumps
- Doesn't influence setting of resin based filling and luting materials
- Based on Zinc oxide cement organic acids

COMPARE TO
**Tempbond
NE**



TempBond

Kerr

Clear Automix Syringe

E430-33351Clear - \$**57.99**

Contains:

6g Syringe and 10 Mixing Tips



NexTemp

Premier Dental

580-3001470**\$69.95** Bulk Package

Contains: 4 - 5ml Automix Dual Barrel Syringes and 10 Mixing Tips



580-3001471**\$32.49** Standard Package

Contains: 1 - 5ml Automix Dual Barrel Syringe, 10 Mixing Tips

Temporary Filling Material

SILMET

ProFil Temp Grey

Silmet

104-19T028**\$11.99** 28g jar



MARK 3™

BEST VALUE

Reinforced Z.O.E. Cement Intermediate Restorative Material

100-5370**\$29.99**

Contains: 1 - 50gm Powder, 1 - 15ml Liquid, 1 Scoop and 1 Dropper

- Good for temporary fillings or as a base
- Sedative qualities on hypersensitive pulp
- Abrasion resistant with good sealing qualities
- 37 MPa compressive strength



COMPARE TO Dentsply IRM

Cavit G Grey

3M/ESPE

E760-44313**\$16.95** 28g/jar



Cavit Pink

3M/ESPE

E760-44030**\$15.95**

Hard 28g/jar



MARK 3

THE SMART ALTERNATIVE

NEW FROM MARK 3 Premium Products at an Affordable Price

Self Adhesive - Self Curing Permanent Resin Cement



See details on page 22

Compare to Maxcem, Multilink Auto and Rely-X Unicem Automix

Non-Eugenol - Automix Syringe Temporary Cement



MADE IN GERMANY

See details on page 23

Compare to TempBond NE, TempGrip, NexTemp, and RelyX Temporary

Composite Accessories

MARK3™

Intra Oral Mixing Tips

- 100-0059Clear
- 100-0058Pink
- 100-0060Yellow

\$16.95 100/pk

- Innovative color-coding system
- Allows impression materials to flow easily while keeping the materials bubble free.



Pre-Bent Needle Tips

BEST VALUE

- 100-2270Pink Tips - 18 ga.
- 100-2280Yellow Tips - 20 ga.
- 100-2290Black Tips - 19 ga.
- 100-2300Gray Tips - 22 ga.
- 100-2310Blue Tips - 25 ga.

\$9.95 100/pk

- Disposable, blunt-end for use with etchants, resins and flowable composites.



Infusor Tips

- 100-8025**\$39.99**
100/pk

- Pre-bent needle with brush-padded end



Mixing Wells

- 100-0002**\$16.95** 2 Well
- 100-0004**\$17.95** 4 Well
200/pk



MARK3™

Mixing Pads

- 100-85843x3 - **\$2.99** 100 sheets/pad
- 10-85876x6 - **\$4.49** 100 sheets/pad



Brush Tips

- 100-2613**\$8.99** 100/bg



Composite Gun/Compute Dispenser

- 100-0790**\$34.95**

- Fits Most Unidose Computes - Including ProFil, TPH3 XRV and Z250



Intraoral Tips XXF

Premier Dental

- 580-3001468**\$24.95** 50/pk
SyringeMix Extra-Fine



Pic-n-Stic

Pulpdent

- 590-PIC**\$18.95** 60/bx



Curing Lights

Ledex WL-090+ Dental Curing Light & Caries Detector

Dentmate

- 725-WL-090+ **\$499.99**



Composites & Restoratives

Ledex Dental Curing Light WL070

Dentmate

- 725-LEDEX-BLLight Blue
 - 725-LEDEX-BLK ..Black
 - 725-LEDEX-GR....Apple Green
 - 725-LEDEX-OR....Orange
 - 725-LEDEX-PKPink
 - 725-LEDEX-WH...White
 - 725-LEDEX-YEL...Yellow
- \$259.99**



Optical Fiber Light Guide Rod

- 725-DC030**\$34.50**
- 8mm/Black Replacement - Universal



Curing Light Sleeves

- 725-DC051**\$13.49** 100/pk



Flowable Composites

Clearfil Majesty Flow

Kuraray America

- 450-2610KAShade A1
- 450-2611KAShade A2
- 450-2612KAShade A3
- 450-2613KAShade A3.5
- 450-2615KAShade B1
- 450-2616KAShade B2

\$58.95

Contains: 1 - 3.2g Syringe and 15 Needle Tips



With orders over \$300
Excludes Canada. Certain restrictions apply.

SILMET

BEST VALUE

ProFil Flow

- 104-160003Shade A1
- 104-160004Shade A2
- 104-160005Shade A3
- 104-160006Shade A3.5
- 104-160007Shade B1
- 104-160008Shade B2
- 104-160009Shade B3
- 104-160010Shade C2

\$41.95

Contains: 4 - 1g Syringe and 10 Tips

The excellent flowable consistency characteristics makes ProFil Flow an ideal composite for filling cavities.



SAMPLES
FREE
available



COMPARE TO
Revolution
Formula 2

Heliomolar Flow Refill

Vivadent

- 820-557030Shade 110/A1
- 820-557031Shade 140/A2
- 820-557032Shade 210/A3
- 820-557033Shade 410/D3

\$29.49

Contains: 1 - 1.6g Syringe and 5 Application Tips



Tetric EvoFlow Syringe

Vivadent

- 820-595953Shade A1
- 820-595954Shade A2
- 820-595955Shade A3
- 820-595956Shade A3.5
- 820-595978Shade A4
- 820-641358Shade B1
- 820-641361Shade B2
- 820-595957Shade B3

2g Syringe

\$40.95



Tetric Evoflow Cavifil

Vivadent

- 820-595987Shade A1
- 820-595988Shade A2
- 820-595989Shade A3
- 820-595990Shade A3.5
- 820-595992Shade A4
- .2g ea. - 20/bx **\$87.95**



Universal Composites

SILMET

ProFil Bulk Fill

104-185A02 \$79.99

1 - 5ml Automix Syringe
Shade A2

Dual cure restorative composite designed for direct restorations. Suitable for bulk fill technique in the posterior region and can be applied in an unlimited layer thickness.

- Unlimited Depth of Cure
- Compatible with all methacrylate-based bonding agents
- Dual Curing Mechanism (can self-cure)
- Faster than the conventional increment technique
- Easy-to-use automix syringe



COMPARE TO
SonicFill

100% MONEY BACK
GUARANTEE
on every product we sell

See Terms and Conditions for details.

SILMET

BEST VALUE

ProFil Syringe Hybrid Composite

- 104-184A01Shade A1
- 104-184A02Shade A2
- 104-184A02Shade A2 Opaque
- 104-184A03Shade A3
- 104-184A35Shade A3.5
- 104-184A04Shade A4
- 104-184B01Shade B1
- 104-184B02Shade B2
- 104-184B03Shade B3
- 104-184C02Shade C2
- 104-184C03Shade C3
- 104-184D03Shade D3
- 4g Syringe **\$39.95**

THE DENTAL
ADVISOR

++++



SAMPLES
FREE
available

COMPARE TO
**Filtek
Supreme**

ProFil Unidose Hybrid Composite

- 104-18UA01Shade A1
- 104-18UA02Shade A2
- 104-18UA02Shade A2 Opaque
- 104-18UA03Shade A3
- 104-18UA35Shade A3.5
- 104-18UA04Shade A4
- 104-18UB01Shade B1
- 104-18UB02Shade B2
- 104-18UB03Shade B3
- 104-18UBWShade BW
- 104-18UC02Shade C2
- 104-18UC03Shade C3
- 104-18UD03Shade D3
- .315 g ea. - 20/pk **\$53.95**

COMPARE TO
**Filtek
Supreme
Unidose**

THE DENTAL
ADVISOR

++++



SAMPLES
FREE
available

ProFil Features:

- Advanced Nano Technology
- Micro hybrid composite resin
- Light cured, radiopaque & color stable
- High polishability
- Non sticky formula
- Low polymerization shrinkage
- Fast placement technique, easy to apply & sculpt
- High strength & wear resistance

Tetric EvoCeram Cavifil

Vivadent

820-590331Shade A1	820-590337Shade B3
820-590332Shade A2	820-590339Shade C2
820-590333Shade A3	820-590340Shade C3
820-590334Shade A3.5	820-637355Shade D2
820-590335Shade A4	.2g ea. - 20/bx \$87.95
820-637354Shade B1	
820-590336Shade B2	



Tetric EvoCeram Syringe

Vivadent

820-590312Shade A1	820-590312Shade C3
820-590313Shade A2	820-637548Shade D2
820-590314Shade A3	820-590322Shade D3
820-590315Shade A3.5	3g Syringe \$63.95
820-590316Shade A4	
820-637547Shade B1	
820-590317Shade B2	
820-590319Shade C1	



Heliomolar Cavifil

Vivadent

820-541501Shade 110/A1
820-541502Shade 140/A2
820-541503Shade 210/A3
820-541505Shade 340/A4
820-541504Shade 310/B3
820-541508Shade 510/C3
.25g ea. - 20/bx \$98.49



Heliomolar Syringe

Vivadent

820-532686Shade 110/A1
820-532687Shade 140/A2
820-532685Shade 210/A3
820-532688Shade 340/A4
820-532689Shade 310/B3
820-532691Shade 410/D3
820-532690Shade 440/D4
3g Syringe \$56.95



Filtek Z-100 Syringes

3M/ESPE

E760-8004A1Shade A1
E760-8004A2Shade A2
E760-8004A3Shade A3
E760-8004A3.5 ..Shade A3.5
4g Syringe \$42.95



Filtek Z-100 Restorative Capsules

3M/ESPE

760-5905A1Shade A1
760-5905A2Shade A2
760-5905A3Shade A3
760-5905A3.5Shade A3.5
.20g - 18/bx \$79.95



Filtek Z-250 Caps

3M/ESPE

760-6021A1Shade A1
760-6021A2Shade A2
760-6021A3Shade A3
760-6021A3.5Shade A3.5
760-6021A4Shade A4
760-6021B1Shade B1
760-6021B2Shade B2
.2g - 20/bx \$79.95



Filtek P-60 Syringe

3M/ESPE

E760-8100A3Shade A3
E760-8100B2Shade B2
4g Refill \$49.95



MARK 3™

BEST VALUE

Core Build Up

- 100-5351.....Shade A3
- 100-5352.....Shade White
- 100-5353.....Shade Blue

\$84.99

Contains: 2 - Automix 5ml Syringes 1:1, 10 Mixing Tips and 10 Intra-Oral Tips

- High Radiopacity
- High Compressive strength
- Contains Nano-zirconium dioxide
- Excellent flow properties and consistency



COMPARE TO
Clearfil
DC Core

**Ketac Silver Aplicap
3M/ESPE**

E760-37010.....**\$229.95** 50/bx

Standard Pack

A silver glass ionomer in a hand-mix format



**Miracle Mix
GC America**

250-000239**\$99.99**

Complete Kit Includes:
15g Powder, 10g Liquid
and 17gm Alloy Powder



250-452100**\$199.95** Capsules, 48/bx

**CompCore
Premier Dental**

Economy Kit

580-3001401Shade White
580-3001404Shade Natural
\$183.95

Contains: 25g Base, 25g Catalyst, Mixing Pads and Spatulas

Single Shade Kit

580-3001402Shade White
580-3001405Shade Natural
\$109.95

Contains: 14g Base, 25g Catalyst, Mixing Pads and Spatulas



**Compcore AF Cartridge STACK:
Dual Cure/Automix**

Premier Dental

Refill Pack

580-3001432Shade White
580-3001433Shade A3
\$148.95

Contains: 1 - 50g Cartridge, 25 Mixing Tips, 25 Intraoral Tips and 10 ACCOR Matrices



Mixing Tips Package

580-3001417**\$50.95**

Contains: 35 Mixing Tips and 35 Intraoral Tips

**CompCore AF SyringeMix STACK:
Dual Cure/Automix**

Premier Dental

Dual-Cure/Automix

580-3001512Shade White
580-3001513Shade A3
580-3001514Shade Blue
\$99.95 per kit

Kit Contains: 2 - 9g Automix Syringes, 15 Mixing Tips, 15 Standard and 10 Extra-Fine Intraoral Tips



**CompCore AF Twist STACK:
Auto Cure/Syringe Hand Mix**

Premier Dental

580-3001490Shade White
580-3001491Shade White A3
\$68.95 per kit

Kit Contains: 6g Syringe Catalyst, 6g Syringe Base, Mixing Pads and 20 Spatulas



SILMET

COMPARE TO
Ketac Silver

ProGlass Silver

104-170SIL.....**\$79.95**

Powder/20g and Liquid/10ml

- Mercury free fluoride releasing silver restorative
- Designed to offer strong abrasion resistance
- Unique blend of spherical silver-tin alloy with glass ionomer, providing an ultimate balance between working and setting time
- Ideal for lasting core build-ups, block-outs and repair, can also be used where radiopacity is required



Bonding Systems

SILMET

BEST VALUE

ProLink Self Etch

104-15SEKT.....\$82.95

Kit Contains:

4ml Bottle,
25 Disposable Applicators
and 1 Dispensing Dish



7th Generation

No Etch Needed

COMPARE TO
Opti Bond
All-in-One

ProLink Universal Adhesive

104-150005.....\$46.95 5ml/btl

5th Generation

Compare to:

Dentsply's Prime & Bond,
Kerr's Optibond and
3M's Singlebond

SAMPLES
FREE
available

Dental PRODUCT SHOPPER
BEST PRODUCT
2014



G-Bond Kit

GC America

250-002277.....\$134.95

Kit Contains:

1 - 5 ml Bottle, 50 Micro-Tips,
1 Micro Tip Handle, 1 Dispensing Dish
and Technique Card



7th Generation Bonding

Optibond Solo Plus

Kerr

E430-29669.....\$199.95

Unidose Refill 100/pk



Clearfil SE Bond

Kuraray

450-1970KA.....\$124.95

Complete Kit Contains:

1- 5ml Bond and 1- 6ml Primer,
Dispenser Case, 2 Brush Holders,
2 - 50 Piece Brush Tips, Mixing Dish
and Light Blocking Plate



Clearfil SE Protect

Kuraray

450-2870KA.....\$194.50

Introductory Kit Contains:

1 - 6ml Primer, 1 - 5ml Bond, 2-50 Disposable Brush Tips,
2 Brush Tip Handles (black/white), 1 Mixing Dish,
1 Light Blocking Cover, and 1 Outer Dispensing Case



IntegraBond

Premier Dental

580-3001411.....\$94.95 7ml/btl
Refill Bottle

580-3001412.....\$94.95 7ml/btl
Auto-Cure Activator



Silane Bond Enhancer

Pulpdent

590-SIL.....\$37.49

Complete Kit Contains:

4 - 1.2 ml Syringe and 8 Dropper Tips

590-SIL-3.....\$21.75 3ml/syringe



DentTASTIC Uno-Duo

Pulpdent

590-UNDO.....\$122.95

Kit Contains:

6ml Bottle of Uno, 3ml Bottle DUO,
5ml Bottle of EtchRite Etching Gel, 20 Applicator Tips,
50 Brush Tips and Brush Handle



590-UNO.....\$122.95

Kit Contains:

2 - 6ml Bottles of Uno, 5ml Bottle of EtchRite Etching Gel,
20 Applicator Tips, 50 Brush Tips and Brush Handle

Embrace WetBond Seal-N-Shine

Pulpdent

590-EMSNY.....\$57.95 Syringe Kit

Kit Contains:

2 - 1.2ml Syringes and 40 Flocked Applicator Tips



590-EMSNS\$67.95 Bottle Kit

Kit Contains:

1 - 6ml Syringes, 1 Brush Handle and 100 Brush Tips

Embrace WetBond Repair Kit

Pulpdent

590-EMPFM.....\$109.95

Intra Oral Repair Kit Contains:

5 - 1.2 ml Syringes (1 of each): First-Coat, Seal-n-Shine, Opaquer, Porcelain Etch Gel and Kool Dam



Monobond Plus

Vivadent

820-626221\$89.99 5g Bottle



Excite F Refill

Vivadent

820-630375\$99.95 5g Bottle



Heliobond Liquid

Vivadent

820-532907\$72.95 6g Bottle



Adper Prompt L-Pop

3M/ESPE

E760-41926\$199.99 100/bx



Adper Single Bond Plus

3M/ESPE

E76051202\$74.95

Vial Refill 6g



37% Phosphoric Acid

Etchant Gel

BEST VALUE

FEATURES & BENEFITS

- Thick, creamy consistency
- Dark blue color creates a good contrast in the mouth
- Rinses completely and quickly



100-9090Four Pack \$12.99 ea.

Contains 4 - 1.2ml Etch Syringes, 20 Pre-Bent Applicator Tips and Instructions

100-9091Bulk Pack \$39.95 ea.

Contains 25 - 1.2ml Etch Syringes, 50 Pre-Bent Applicator Tips and instructions

100-9092Jumbo Kit \$30.95 ea.

Contains 1 - 50ml Etch Syringe, 4 - 1.2ml Empty Syringes, 50 Pre-Bent Applicator Tips and instructions

Compare to
Etch-Rite

Etching Products

SILMET

ProEtch

4 Pack

104-190100\$14.95

37% Ph. Acid Etch Gel -
1.2ml Syringe - 4/pk
Includes Delivery Tips

- High quality, non drying, medium viscosity gel stays put and does not run.
- Excellent water solubility. Washes off quickly & easily.
- Large & comfortable to use syringes



Premier Etch

Premier Dental

580-3001422\$24.95

37% Etch

Kit Contains: 1 - 12 g jumbo Luer-Lock Syringe, 10 empty 1.2 ml Syringes and 50 Disposable Tips



38% Etch Microdose

580-Premier Dental

580-4009051\$32.99

Refill Kit Contains: 50 - 0.2cc Pipettes



Etch Rite Etching Gel

Puldent

590-ETCH\$16.75

Etch Gel Syringe - 4/bx

590-ET-24\$66.95

Bulk Pack 1.2ml - 24/pk

590-ET-50\$44.95

Jumbo Syringe Kit - 2 - 25ml

590-ET-50R\$34.95

Jumbo Refill - 2 - 25ml

590-25B50\$10.15

Applicator Tips - 50/bg



Etch Royale

Puldent

590-ER24\$64.95

Bulk Pack: 24 - 1.2ml Syringes



Porcelain Etch Gel

Puldent

590-PEG\$41.49

Kit Contains: 4 x 1.2ml Syringes and 8 Pre-bent Applicators



Porcelain Prep Kit

Puldent

590-PPK\$43.50

Kit Contains: 1 - 1.2 ml Syringe Porcelain Etch Gel, 2 - 1.2 ml Syringes Dry-Rite Drying Agent, 1 - 1.2 ml Syringe Silane Bond Enhancer and 8 Dispenser Tips.



Micro Applicators

Microbrush Disposable Applicators

Microbrush

470-MRA400Regular Assorted

470-MRB400Regular Blue

470-MRG400Regular Green

470-MRS400Regular Peach

470-MRP400Regular Purple

470-MSF400Regular Super Fine White

470-MFA400Fine Assorted

470-MFY400Fine Yellow

\$34.95 400/pk



470-PF400Plus Fine Assorted

470-PR400Plus Regular Assorted

470-PSF400Plus Superfine White

\$35.95 400/pk



Micro Applicator Brushes

Unipak

790-UMA-9012-FFine Assorted Blue and Yellow

790-UMA-9013-RRegular Assorted Green & Orange

\$17.95 400/bx

Tooth Whitening

iBrite LED Whitening System

PacDent

108-IB-888\$1339.95

Unit Contains:

1 Stand, 1 Control Box, 1 Base, 5 Wheels, 1 Headlamp and 1 Mouthpiece



iBrite Automix Gel Type Whitening System

PacDent

Single Patient Pack

108-IB-801\$38.49

Single Patient Pack Contains:
1 - 3.5ml 30% HP Syringe and accessories



Five Patient Pack

108-IB-803\$153.29

5 Patient Pack Contains: 2 - 6.4ml 30% HP Syringes and accessories

Perfecta REV Tooth Whitening

Premier Dental

Patient Pack

580-4000142\$35.95

Contains: 2 - 3cc Syringes, 1 - 4oz. Finishing Rinse and Doctor and Patient Instructions



Refresher Pack

580-4000141\$15.95

Contains: 1 - 3cc Syringe, 1 - 4 oz. Finishing Rinse and Patient Instructions



Perfecta Bravo Tooth Whitening

Premier Dental

Patient Pack

580-4000092\$45.99

Package Contains: 4 - 3cc Syringes, 4 Dispensing Tips, Patient Instructions with Shade Guide, Doctor's Instructions, Cosmetic Patient Pouch with Mirror, Tray Material and Tray Case



Value Pack

580-4000093\$181.95

Package Contains: 24 - 3cc Syringes, 24 Dispensing Tips, Patient Instructions with Shade Guide, Doctor's Instructions, Cosmetic Pouches - Enough Material for 6 Patients

50 Pack Bulk4000094

580-4000094\$329.95

Package Contains: 50 - 3cc Syringes, 50 Dispensing Tips, Patient Instructions with Shade Guide, Doctor's Instructions Enough Material for 12 Patients

Veneer System

SILMET

ProVeneer

In Office Veneer System

Assorted Kit

104-18PVAST\$325.00

Kit Includes: 10 x RBFM Single-use Packs (3x A1, 2x A2, 2x B1, 3x BW), 1 x ProVeneer Bonding Agent, 1 x ProFil Flow Cement, 1 x ProEtch and Accessories



Kit makes 10 Veneers

FREE Marketing Materials come with your kit first order:

1 Poster, 50 Brochures, 50 Dental History Forms and Flipbook

Resin Based Foundation Blister Refills

104-18PV10-A1.....Shade A1

104-18PV10-A2.....Shade A2

104-18PV10-B1.....Shade B1

104-18PV10-BW.....Shade BW (Bleach White)

\$225.00 10/pk

ProVeneer Sculpting Instrument

104-GFCCI.....Double Ended \$24.99

- New approach for preparing in office veneers
- Create a single veneer in 15 minutes or less
- Natural looking results
- Affordable price

In just 2 minutes, ProVeneers' patent pending Smile Preview Veneer System pre-qualifies patients, without preformed veneers, and assures gorgeous results.



Core Materials

LuxaCore Dual

DMG America

Cartridge Refill

- 190-212011A3
- 190-212012Blue
- 190-212126White

Contains: **\$184.95**

1 - 50g Cartridge, 35 Mixing Tips and 35 Intraoral Tips



Smartmix Handheld Syringe

- 190-212039A3
- 190-212038Blue
- 190-212064White

Contains: **\$139.95**

2 - 9g Syringes, 10 Mixing Tips and 10 Intraoral Tips



Ti-Core Natural

Essential Dental Systems

- 235-810-00.....Regular Set
- 235-815-00.....Fast Set

Contains: **\$143.95** Kit

15g Catalyst, 14g Base, 1 Mixing Pad and 100 Plastic Spatulas



Ti-Core Gray

Essential Dental Systems

- 235-800-00.....Regular Set
- 235-805-00.....Fast Set

Contains: **\$143.95** Kit

15g Catalyst, 14g Base, 1 Mixing Pad and 100 Plastic Spatulas



CompCore

Premier Dental

Economy Kit

- 580-3001401Type White
- 580-3001404Type Natural

Contains: **\$183.95**

25g Base, 25g Catalyst, Mixing Pads and Spatulas



Single Shade Kit

- 580-3001402Type White
- 580-3001405Type Natural

Contains: **\$109.95**

14g Base, 25g Catalyst, Mixing Pads and Spatulas

HardCore

Pulpdent

590-HC.....**\$95.95**

Standard Kit Contains:

2 - 10g Jars of each (Base & Catalyst),
20 Assorted Core Forms, Mixing Pad
and Spatulas



590-HC-R**\$46.95**

Refill Package Contains: 10g Jar of each (Base & Catalyst)

Temporary Crown & Bridge

Luxatemp Automix Plus

DMG America

Refill Package

- 190-110401Shade A1
- 190-110402Shade A2
- 190-110403Shade A3.5
- 190-110404Shade B1

\$186.95

Contains: 1 - 76g Cartridge and 15 Mixing Tips



Mixing Tips for Luxatemp

100-110351**\$21.95** 25/bx

Luxatemp Ultra

DMG America

- 190-110906Shade A1
- 190-110907Shade A2
- 190-110908Shade A3
- 190-110909Shade A3.5
- 190-110910Shade B1

\$189.99

Refill Package Contains: 1 - 76g Cartridge and 15 Mixing Tips



MARK3™
THE SMART ALTERNATIVE



10:1 Automix
Temporary Crown & Bridge

FEATURES & BENEFITS

- Biocompatible
- Easy handling and polishing
- Improved precision of application
- Lifelike restoration with natural translucency
- Excellent resistance to fracture

COMPARE TO
Luxatemp

Technical Data

- Dosage and application time: Approx. 45 seconds
- Setting time in mouth: Approx. 2-3 minutes (hard elastic phase, from beginning of mixing)
- Setting time: Approx. 6-7 minutes (until setting time is completed)
- Mixing ratio: 10:1

- 100-4101.....Shade A1
- 100-4102.....Shade A2
- 100-4103.....Shade A3

BEST VALUE

\$99.49 1 Cartridge of 50ml (75 gm) plus 10 Mixing Tips

LuxaFlow Ultra

DMG America

- 190-224001Shade A1
- 190-224002Shade A2
- 190-224003Shade A3
- 190-224004Shade A3.5
- 190 224005Shade B1



\$44.49

Refill Package Contains: 2 - 1.5g Syringes and 10 Tips

Protemp 4 Garant

3M/ESPE

- E760-46954Shade A1
- E760-46956Shade A2
- E760-46957.....Shade A3



\$124.95

Standard Pack Contains:
1 - 50ml Cartridge and 16 Garant Mixing Tips

MARK3™

Dispensing Gun 10:1/4:1 - Type 50

100-110411**\$65.50**

- Dispensing Gun Luxatemp Integrity 10:1 by MARK3 is an all-purpose tool used for dispensing most crown and bridge materials.
- It can also be used for Automix Crown and Bridge materials; Auto Mix Core Build-Up; and many more.
- This dispenser fits all standard cartridge sizes with a 10:1 or 4:1 ratio.
- Dispensing Gun Luxatemp Integrity 10:1 by MARK3 has Superior flexural strength for durable temporaries.



Crown & Bridge Products

CROWN FORMS

Polycarbonate Crowns

MARK3™

MARK3 Crowns

Will not crack, split, shatter or chip. Easily scissor trimmed, contoured, crimped and polished. Thin design prevents gingival irritation. May be cemented with any plastic, zinc oxide or orthophosphate cement.



C180\$69.95

Crown Kit - 180 pcs (3 of each size)

Replacement Crowns \$4.95 5/bx

UPPER CENTRAL Right

Size	Item#
10	100-P10
11	100-P11
12	100-P12
13	100-P13
14	100-P14
100	100-P100
101	100-P101

UPPER CENTRAL Left

Size	Item#
15	100-P15
16	100-P16
17	100-P17
18	100-P18
19	100-P19
102	100-P102
103	100-P103

UPPER LATERAL Right

Size	Item#
2	100-P2
20	100-P20
21	100-P21
22	100-P22
23	100-P23
24	100-P24

UPPER LATERAL Left

Size	Item#
25	100-P25
26	100-P26
27	100-P27
28	100-P28
29	100-P29
200	100-P200

RIGHT CUSPID Upper & Lower

Size	Item#
30	100-P30
31	100-P31
32	100-P32
33	100-P33
34	100-P34
300	100-P300
301	100-P301

LEFT CUSPID Upper & Lower

Size	Item#
35	100-P35
36	100-P36
37	100-P37
38	100-P38
39	100-P39
302	100-P302
303	100-P303

1st BICUSPID

Size	Item#
40	100-P40
41	100-P41
42	100-P42
43	100-P43
44	100-P44

2nd BICUSPID

Size	Item#
50	100-P50
51	100-P51
52	100-P52
53	100-P53
54	100-P54

LOWER ANTERIOR Long

Size	Item#
60	100-P60
61	100-P61
62	100-P62
63	100-P63
64	100-P64

LOWER ANTERIOR Short

Size	Item#
65	100-P65
66	100-P66
67	100-P67
68	100-P68
69	100-P69

Stainless Steel Crowns

MARK3™

BEST VALUE

Stainless Steel 1st & 2nd Primary Molar Pedo Crowns

Surgical grade steel that is Pre-formed and Pre-Crimped for easy placement. Shallow occlusal anatomy to preserve tooth structure. 100% compatible with other leading name brands.

50-60% below expensive name brands.



Replacement Crowns \$24.95 5/bx

1st Primary Molar

UPPER RIGHT	UPPER LEFT	LOWER RIGHT	LOWER LEFT
100-DUR2	100-DUL2	100-DLR2	100-DLL2
100-DUR3	100-DUL3	100-DLR3	100-DLL3
100-DUR4	100-DUL4	100-DLR4	100-DLL4
100-DUR5	100-DUL5	100-DLR5	100-DLL5
100-DUR6	100-DUL6	100-DLR6	100-DLL6
100-DUR7	100-DUL7	100-DLR7	100-DLL7

2nd Primary Molar

UPPER RIGHT	UPPER LEFT	LOWER RIGHT	LOWER LEFT
100-EUR2	100-EUL2	100-ELR2	100-ELL2
100-EUR3	100-EUL3	100-ELR3	100-ELL3
100-EUR4	100-EUL4	100-ELR4	100-ELL4
100-EUR5	100-EUL5	100-ELR5	100-ELL5
100-EUR6	100-EUL6	100-ELR6	100-ELL6
100-EUR7	100-EUL7	100-ELR7	100-ELL7

COMPARE TO
**3M
Crowns**



Primary Molar Micro Crowns

3M/ESPE

Replacement Crowns **\$34.95** 5/bx



1st UPPER LEFT

Size	Item#
D-UL-2	E760-DUL2
D-UL-3	E760-DUL3
D-UL-4	E760-DUL4
D-UL-5	E760-DUL5
D-UL-6	E760-DUL6
D-UL-7	E760-DUL7

1st UPPER RIGHT

Size	Item#
D-UR-2	E760-DUR2
D-UR-3	E760-DUR3
D-UR-4	E760-DUR4
D-UR-5	E760-DUR5
D-UR-6	E760-DUR6
D-UR-7	E760-DUR7

1st LOWER LEFT

Size	Item#
D-LL-2	E760-DLL2
D-LL-3	E760-DLL3
D-LL-4	E760-DLL4
D-LL-5	E760-DLL5
D-LL-6	E760-DLL6
D-LL-7	E760-DLL7

1st LOWER RIGHT

Size	Item#
D-LR-2	E760-DLR2
D-LR-3	E760-DLR3
D-LR-4	E760-DLR4
D-LR-5	E760-DLR5
D-LR-6	E760-DLR6
D-LR-7	E760-DLR7

2nd UPPER LEFT

Size	Item#
E-UL-2	E760-EUL2
E-UL-3	E760-EUL3
E-UL-4	E760-EUL4
E-UL-5	E760-EUL5
E-UL-6	E760-EUL6
E-UL-7	E760-EUL7

2nd UPPER RIGHT

Size	Item#
E-UR-2	E760-EUR2
E-UR-3	E760-EUR3
E-UR-4	E760-EUR4
E-UR-5	E760-EUR5
E-UR-6	E760-EUR6
E-UR-7	E760-EUR7

2nd LOWER LEFT

Size	Item#
E-LL-2	E760-ELL2
E-LL-3	E760-ELL3
E-LL-4	E760-ELL4
E-LL-5	E760-ELL5
E-LL-6	E760-ELL6
E-LL-7	E760-ELL7

2nd LOWER RIGHT

Size	Item#
E-LR-2	E760-ELR2
E-LR-3	E760-ELR3
E-LR-4	E760-ELR4
E-LR-5	E760-ELR5
E-LR-6	E760-ELR6
E-LR-7	E760-ELR7

Temporary Crowns

Iso-Form Bicuspid Crowns

3M/ESPE

Replacement Crowns **\$24.95** 5/bx



1st

UPPER RIGHT	UPPER LEFT	LOWER RIGHT	LOWER LEFT
760-U40	760-U41	760-L40	760-L41
760-U42	760-U43	760-L42	760-L43
760-U44	760-U45	760-L44	760-L45
760-U46	760-U47	760-L46	760-L47
760-U48	760-U49	760-L48	760-L49

2nd

UPPER RIGHT	UPPER LEFT	LOWER RIGHT	LOWER LEFT
760-U50	760-U51	760-L50	760-L51
760-U52	760-U53	760-L52	760-L53
760-U54	760-U55	760-L54	760-L55
760-U56	760-U57	760-L56	760-L57
760-U58	760-U59	760-L58	760-L49

Iso-Form Molar Crowns

3M/ESPE

Replacement Crowns **\$24.95** 5/bx

1st

UPPER RIGHT	UPPER LEFT	LOWER RIGHT	LOWER LEFT
760-U60	760-U61	760-L60	760-L61
760-U62	760-U63	760-L62	760-L63
760-U64	760-U65	760-L64	760-L65
760-U66	760-U67	760-L66	760-L67
760-U68	760-U69	760-L68	760-L69

2nd

UPPER RIGHT	UPPER LEFT	LOWER RIGHT	LOWER LEFT
760-U70	760-U71	760-L70	760-L71
760-U72	760-U73	760-L72	760-L73
760-U74	760-U75	760-L74	760-L75
760-U76	760-U77	760-L76	760-L77
760-U78	760-U79	760-L78	760-L79

Applicators

Cotton Tipped Applicators

Safe-Dent

- 102-6030.....**\$5.49** 3"/1000/bx
 102-6060.....**\$6.65** 6"/1000/bx



Consumables

MARK 3™

Cotton Rolls

- 100-2000..... **\$14.95** 2000/bx
 1.5" x .375" Non Sterile #2
 100% Pure White Cotton Fiber



Dental Bibs 2-ply

- 100-1319B Blue
 100-1319W White
 100-1319PI Pink
 100-1319LV Lavender
 100-1319G Green
 100-1319PE Peach
 100-1319Y Yellow



2-ply paper, 1 ply poly - 13"x 18" **\$19.99** 500/bx

- Maximum protection with horizontal embossed pattern and a unique water repellent edge.
- Increased strength and tear resistance with plastic backing for added clothing protection.

Bib Holder - Autoclavable Silicone

- 100-090B..... Blue
 100-090W..... White
 100-090Y..... Yellow
 100-090R Red

\$2.99



- Non-slip design
- Fully autoclavable up to 275°F
- Medical grade silicone
- Stainless steel alligator clips

C-Fold Towels

Safe-Dent

- 102-200 **\$32.99** 2400/cs
 New and improved stronger and thicker towel



Non-Woven N/S Gauze Sponges

UniPak

- 790-UNW-3522 **\$29.95** 5000/bx
 2"x2" 4-Ply
 790-UNW-3044 **\$39.95** 2000/bx
 4"x4" 4-Ply



- No chemicals added to avoid allergic reactions
- 100% more absorbent

Plastic 5 oz. Cups

UniPak

- 790-UBC-6201Blue
 790-UBC-6202White
 790-UBC-6203Pink
 790-UBC-6204Lavender
 790-UBC-6205Green
 790-UBC-6206Yellow



\$21.95 5 oz. - 1000/bx

Tray Covers Ritter "B"

UniPak

- 790-UBC-80141Blue
 790-UBC-80142White
 790-UBC-80144Lavender
 790-UBC-80145Green
 790-UBC-80146Yellow
 790-UBC-80149Mauve



8.5" x 12.25" **\$19.95** 1000/bx

- Tray covers improved moisture resistance helps prevent cross-contamination and unwanted water seepage

MARK 3™

Tray Sleeves

- 100-2105.....**\$17.49** 500/bx
 Regular - 10.5" x 14"



Tray Sleeves - Extra Large

UniPak

- 790-UBC-8011A **\$19.99** 500/bx
 Large - 11.625" x 14.5"

Monoject 412 Syringe

Covidien/Kendall

- 520-412012**\$27.95** 50/bx
 12cc Disposable Syringe Curved Tip




MARK 3™

BEST VALUE

Curved 412 Syringe

100-5212 **\$18.95** 12cc/50/bx

- Disposable, non-sterile syringes
- Designed for surgical site and post-surgical home irrigation
- Can also be used with silicone and rubber based impression materials for intra-oral delivery



Surface Barriers

MARK 3™

Pinnacle Quality Barriers

Chair Covers

100-2124 **\$25.95** 225/bx
Half Chair - 27.5" x 24"

100-2125 **\$33.95** 125/bx
Full Chair - 29" x 80"



T-Style Sleeves

100-2127 **\$14.95** 500/bx

- Clear plastic sleeve for T-Style Handles
- 4" x 5.75"



Headrest Covers

Plastic - Latex Free

100-1402 **\$9.49** 250/pk
Regular - 9.5" x 11"

100-1403 **\$11.95** 250/pk
Large - 9.5" x 14"




Headrest Covers

UniPak

Paper - White

790-UBC-80262 **\$22.95** 500/pk
Regular - 10" x 10"

790-UBC-80272 **\$28.95** 500/pk
Large - 10" x 13"



Syringe Tip Covers

Palmero Health Care

565-1950 **\$15.95** 100/bg



Air/Water Syringe Tips


BEST VALUE

MARK 3™

Multi-Colored Air/Water Syringe Tips

100-1202 **\$35.99** 250/bg

100-1203 **\$159.99** 1500/bg



COMPARE TO Crystal Tips

- Easy to use - no adapter needed
- Superior air/water spray
- Completely beveled ends-no sharp edges
- Made from FDA approved plastics
- Unique quad aircore design ensures dry air when you need it and eliminates water leaks

Clear Tips Plastic Air/Water Syringe Tips

100-1410 **\$238.99** 1500/pk

100-1407 **\$49.99** 250/pk

100-1408 **\$23.99** 100/pk



COMPARE TO Sani-Tip

- Produces dry air for optimal bonding results
- Eliminates possibility of cross contamination
- Clear tip to verify sanitary condition
- Easy to change between patients

Quick Tips Metal Air/Water Syringe Tips

100-1411 **\$278.99** 1600/pk

100-1409 **\$31.95** 100/pk



COMPARE TO Safe-Tips EZ

- Interchangeable with most metal tips
- No conversion kits or adapters necessary
- Separate stainless steel water tube prevents air-water crossover
- Interchangeable with A-dec and DCI style syringes

LockTight Air/Water Syringe Tips

COMPARE TO Seal Tight

100-1205 **\$35.95** 200/pk

100-1206 **\$259.95** 1500/pk



- Eliminate the need for tightening and loosening an adaptor
- Push in place for a perfect seal every time

Drying/Obturation

Meta Paper Points

Regular

106-PPCB-#10#10	106-PPCB-XXFXX Fine
106-PPCB-#15#15	106-PPCB-XFX Fine
106-PPCB-#15/40 ..#15/40	106-PPCB-FFine
106-PPCB-#20#20	106-PPCB-MMedium
106-PPCB-#25#25	106-PPCB-CCoarse
106-PPCB-#30#30	106-PPCB-XCExtra Coarse
106-PPCB-#35#35	
106-PPCB-#40#40	
106-PPCB-#45#45	
106-PPCB-#45/80 ..#45/80	



\$5.99
200/pk

Cell Pack

106-PP200-#10#10	106-PP200-XXFXX Fine
106-PP200-#15#15	106-PP200-XFX Fine
106-PP200-#15/40 ..#15/40	106-PP200-FFine
106-PP200-#20#20	106-PP200-MMedium
106-PP200-#25#25	106-PP200-CCoarse
106-PP200-#30#30	106-PP200-XCExtra Coarse
106-PP200-#35#35	
106-PP200-#40#40	
106-PP200-#45#45	
106-PP200-#45/80 #45/80	



\$5.99
200/pk

.04 Tapered

106-PP04-#10#10	106-PP04-#40#40
106-PP04-#15#15	106-PP04-#45#45
106-PP04-#15/40 ..#15/40	106-PP04-#45/80 ..#45/80
106-PP04-#20#20	
106-PP04-#25#25	
106-PP04-#30#30	
106-PP04-#35#35	



\$6.95
60/pk

.06 Tapered

106-PP06-#10#10	106-PP06-#40#40
106-PP06-#15#15	106-PP06-#45#45
106-PP06-#15/40 ..#15/40	106-PP06-#45/80 ..#45/80
106-PP06-#20#20	
106-PP06-#25#25	
106-PP06-#30#30	
106-PP06-#35#35	



\$6.95
60/pk


Pac-Dent Paper Points

- ISO sizes and color-coded
- Available in .02, .04, and .06 tapers
- Laser mark line at 16, 18, 19, 20, 22 and 24mm

SIZE	.02 TAPER	.04 TAPER	.06 TAPER
15	108-PP-15-02	108-PP-15-04	108-PP-15-06
20	108-PP-20-02	108-PP-20-04	108-PP-20-06
25	108-PP-25-02	108-PP-25-04	108-PP-25-06
30	108-PP-30-02	108-PP-30-04	108-PP-30-06
35	108-PP-35-02	108-PP-35-04	108-PP-35-06
40	108-PP-40-02	108-PP-40-04	108-PP-40-06
45	108-PP-45-02	108-PP-45-04	
50	108-PP-50-02	108-PP-50-04	
55	108-PP-55-02	108-PP-55-04	
60	108-PP-60-02	108-PP-60-04	
70	108-PP-70-02	108-PP-70-04	
80	108-PP-80-02	108-PP-80-04	
15-40	108-PP-AS15-02	108-PP-AS15-04	108-PP-AS15-06
45-80	108-PP-AS45-04	108-PP-AS45-04	

\$6.99

Sold in a box of 120
(6 vials containing
20 points each)



CAN'T FIND IT?

We're always adding new products. If you don't see it here, visit our website, give us a call or email us at: info@nobledentalsupplies.com

Meta Gutta Percha Points

Regular

106-GPCB-#10#10	106-GPCB-FFFine Fine
106-GPCB-#15#15	106-GPCB-MFMedium Fine
106-GPCB-#15/40 ..#15/40	106-GPCB-FFine
106-GPCB-#20#20	106-GPCB-FMFine Medium
106-GPCB-#25#25	106-GPCB-MMedium
106-GPCB-#30#30	106-GPCB-MLMedium Large
106-GPCB-#35#35	106-GPCB-LLarge
106-GPCB-#40#40	106-GPCB-XLX Large
106-GPCB-#45#45	
106-GPCB-#45/80 ..#45/80	
106-GPCB-#50#50	
106-GPCB-XFX Fine	



\$5.99
120/pk

Vials

106-GPCV-#10#10	106-GPCV-FFFine Fine
106-GPCV-#15#15	106-GPCV-MFMedium Fine
106-GPCV-#15/40 ..#15/40	106-GPCV-FFine
106-GPCV-#20#20	106-GPCV-FMFine Medium
106-GPCV-#25#25	106-GPCV-MMedium
106-GPCV-#30#30	106-GPCV-MLMedium Large
106-GPCV-#35#35	106-GPCV-LLarge
106-GPCV-#40#40	106-GPCV-XLX Large
106-GPCV-#45#45	
106-GPCV-#45/80 ..#45/80	
106-GPCV-#50#50	
106-GPCV-XFX Fine	



\$5.99
120/pk

.04 Tapered

106-GP04-#10#10	106-GP04-#40#40
106-GP04-#15#15	106-GP04-#45#45
106-GP04-#15/40 ..#15/40	106-GP04-#45/80 ..#45/80
106-GP04-#20#20	106-GP04-#50#50
106-GP04-#25#25	
106-GP04-#30#30	
106-GP04-#35#35	



\$8.25
60/pk

.06 Tapered

106-GP06-#10#10	106-GP06-#40#40
106-GP06-#15#15	106-GP06-#45#45
106-GP06-#15/40 ..#15/40	106-GP06-#45/80 ..#45/80
106-GP06-#20#20	106-GP06-#50#50
106-GP06-#25#25	
106-GP06-#30#30	
106-GP06-#35#35	



\$8.25
60/pk


Pac-Dent Gutta Percha Points

- ISO sizes and color-coded
- Available in .02, .04, and .06 tapers
- Laser mark line at 16, 18, 19, 20, 22 and 24mm

SIZE	.02 TAPER	.04 TAPER	.06 TAPER
15	108-GP-15-02	108-GP-15-04	108-GP-15-06
20	108-GP-20-02	108-GP-20-04	108-GP-20-06
25	108-GP-25-02	108-GP-25-04	108-GP-25-06
30	108-GP-30-02	108-GP-30-04	108-GP-30-06
35	108-GP-35-02	108-GP-35-04	108-GP-35-06
40	108-GP-40-02	108-GP-40-04	108-GP-40-06
45	108-GP-45-02	108-GP-45-04	
50	108-GP-50-02	108-GP-50-04	
55	108-GP-55-02	108-GP-55-04	
60	108-GP-60-02	108-GP-60-04	
70	108-GP-70-02	108-GP-70-04	
80	108-GP-80-02	108-GP-80-04	
15-40	108-GP-AS15-02	108-GP-AS15-04	108-GP-AS15-06
45-80	108-GP-AS45-04	108-GP-AS45-04	

\$6.99

Sold in a box of 120
(6 vials containing
20 points each)



Please Check Our Website
For many more Point types and sizes

PacFil™ Thermal Obturators

Pac-Dent

- Calibration rings facilitates setting working length
- Works with your current obturator oven at the same settings
- Easy to get great 3D fills for all canal types
- ISO-sized system for compatibility

PacFil™ GPC (Gutta Percha Carrier)

Size	6-Pack	20-Pack
20	108-GOB-20-06	108-GOB-20-20
25	108-GOB-25-06	108-GOB-25-20
30	108-GOB-30-06	108-GOB-30-20
35	108-GOB-35-06	108-GOB-35-20
40	108-GOB-40-06	108-GOB-40-20
45	108-GOB-45-06	108-GOB-45-20
50	108-GOB-50-06	108-GOB-50-20

\$45.99 6/pk
\$139.99 20/pk

PacFil™ PC (Gutta Plastic Core Carrier)

Size	6-Pack	20-Pack
20	108-POB-20-06	108-POB-20-20
25	108-POB-25-06	108-POB-25-20
30	108-POB-30-06	108-POB-30-20
35	108-POB-35-06	108-POB-35-20
40	108-POB-40-06	108-POB-40-20
45	108-POB-45-06	108-POB-45-20
50	108-POB-50-06	108-POB-50-20

\$31.99 6/pk
\$99.99 20/pk



PacFil™ Verifiers

Pac-Dent

- Match obturators to confirm accurate sizing
- Calibration rings for setting accurate working length
- Made of polycarbonate

Size	6-Pack
20	108-PV-20
25	108-PV-25
30	108-PV-30
35	108-PV-35
40	108-PV-40
45	108-PV-45
50	108-PV-50

\$14.99 6/pk



Hand Files

MARK 3™



K-Files - Stainless Steel

100-1210621mm #06	100-1253525mm #35
100-1210821mm #08	100-1254025mm #40
100-1211021mm #10	100-1254525mm #45
100-1211521mm #15	100-1254580 ...25mm #45/80
100-1211540 ...21mm #15/40	100-1255525mm #55
100-1212021mm #20	100-1256025mm #60
100-1212521mm #25	100-1257025mm #70
100-1213021mm #30	100-1310631mm #06
100-1213521mm #35	100-1310831mm #08
100-1214021mm #40	100-1311031mm #10
100-1214521mm #45	100-1311531mm #15
100-1215021mm #50	100-1311540 ...31mm #15/40
100-1215521mm #55	100-1312031mm #20
100-1216021mm #60	100-1312531mm #25
100-1217021mm #70	100-1313031mm #30
100-1218021mm #80	100-1313531mm #35
100-1250625mm #06	100-1314031mm #40
100-1250825mm #08	100-1314531mm #45
100-1251025mm #10	100-1314580 ...31mm #45/80
100-1251525mm #15	100-1315031mm #50
100-1251540 ...25mm #15/40	100-1315531mm #55
100-1252025mm #20	100-1316031mm #60
100-1252525mm #25	100-1317031mm #70
100-1253025mm #30	\$5.99 6/bx

COMPARE TO
Sybron Endo

BEST VALUE

Hedstrom Files

Pac-Dent

NiTi Hedstrom Files (Hand)

Size/ Length	21mm	25mm	28mm	31mm
#40	108-NHH-40-21	108-NHH-40-25	108-NHH-40-28	108-NHH-40-31
#35	108-NHH-35-21	108-NHH-35-25	108-NHH-35-28	108-NHH-35-31
#30	108-NHH-30-21	108-NHH-30-25	108-NHH-30-28	108-NHH-30-31
#25	108-NHH-25-21	108-NHH-25-25	108-NHH-25-28	108-NHH-25-31
#20	108-NHH-20-21	108-NHH-20-25	108-NHH-20-28	108-NHH-20-31
#15	108-NHH-40-21	108-NHH-15-25	108-NHH-15-28	108-NHH-15-31
#15-40	108-NHH-15/40-21	108-NHH-15/40-25	108-NHH-15/40-28	108-NHH-15/40-31



\$8.99 6/pk

Stainless Steel Hedstrom Files (Hand)

Size/ Length	21mm	25mm	28mm	31mm
#40	108-SHH-40-21	108-SHH-40-25	108-SHH-40-28	108-SHH-40-31
#35	108-SHH-35-21	108-SHH-35-25	108-SHH-35-28	108-SHH-35-31
#30	108-SHH-30-21	108-SHH-30-25	108-SHH-30-28	108-SHH-30-31
#25	108-SHH-25-21	108-SHH-25-25	108-SHH-25-28	108-SHH-25-31
#20	108-SHH-20-21	108-SHH-20-25	108-SHH-20-28	108-SHH-20-31
#15	108-SHH-40-21	108-SHH-15-25	108-SHH-15-28	108-SHH-15-31
#15-40	108-SHH-15/40-21	108-SHH-15/40-25	108-SHH-15/40-28	108-SHH-15/40-31

\$5.99 6/pk

Reamers

Pac-Dent

NiTi Reamers (Hand)

Size/ Length	21mm	25mm	28mm	31mm
#40	108-NRH-40-21	108-NRH-40-25	108-NRH-40-28	108-NRH-40-31
#35	108-NRH-35-21	108-NRH-35-25	108-NRH-35-28	108-NRH-35-31
#30	108-NRH-30-21	108-NRH-30-25	108-NRH-30-28	108-NRH-30-31
#25	108-NRH-25-21	108-NRH-25-25	108-NRH-25-28	108-NRH-25-31
#20	108-NRH-20-21	108-NRH-20-25	108-NRH-20-28	108-NRH-20-31
#15	108-NRH-40-21	108-NRH-15-25	108-NRH-15-28	108-NRH-15-31
#15-40	108-NRH-15/40-21	108-NRH-15/40-25	108-NRH-15/40-28	108-NRH-15/40-31



\$8.99 6/pk

Stainless Steel Reamers (Hand)

Size/ Length	21mm	25mm	28mm	31mm
#40	108-SRH-40-21	108-SRH-40-25	108-SRH-40-28	108-SRH-40-31
#35	108-SRH-35-21	108-SRH-35-25	108-SRH-35-28	108-SRH-35-31
#30	108-SRH-30-21	108-SRH-30-25	108-SRH-30-28	108-SRH-30-31
#25	108-SRH-25-21	108-SRH-25-25	108-SRH-25-28	108-SRH-25-31
#20	108-SRH-20-21	108-SRH-20-25	108-SRH-20-28	108-SRH-20-31
#15	108-SRH-40-21	108-SRH-15-25	108-SRH-15-28	108-SRH-15-31
#15-40	108-SRH-15/40-21	108-SRH-15/40-25	108-SRH-15/40-28	108-SRH-15/40-31

\$5.75 6/pk

Irrigation

PacEndo™ Pre-Filled Irrigation Kit

Pac-Dent

- Packed in a convenient ready-to-use tear-away tray
- Includes all necessary solutions for proper canal disinfection
- Assures fresh effective solutions are used
- Prevents spills and accidents
- Number and color coding limits sequence mistakes



Canal Aspirator fits any HVE handpieces

108-EC-1A PacEndo™ Single Canal kit

PacEndo™ Single Canal kit includes:

- 1 X 10 ml pre-filled Sodium Hypochlorite syringe
- 1 X 3 ml pre-filled EDTA syringe
- 1 X 3 ml pre-filled 95% Ethyl Alcohol syringe
- 1 X 3 ml pre-filled Chlorhexidine syringe
- 4 X 30 ga. irrigation needle tips
- 1 X Canal aspirator with tip

\$5.99

108-EC-1C PacEndo™ Multiple Canals kit

PacEndo™ Multiple Canals kit includes:


- 2 X 15 ml pre-filled Sodium Hypochlorite syringe
- 1 X 9 ml pre-filled EDTA syringe
- 1 X 9 ml pre-filled 95% Ethyl Alcohol syringe
- 1 X 9 ml pre-filled Chlorhexidine syringe
- 5 X 30 ga. irrigation needle tips
- 1 X Canal aspirator with tip

\$9.99

MARK 3™

Luer-Lock Endo Irrigation Syringes

- 100-5260 \$9.99** 3cc/100/pk
- 100-5261 \$10.49** 6cc/50/pk
- 100-5262 \$16.99** 12cc/100/pk
- 100-5263 \$17.95** 20cc/50/pk



- Ideal for surgical site irrigation
- Compatible with all luer style needles and hubs
- Disposable, non-sterile

ENDODONTIC PRODUCTS

Medicaments & Canal Prep Solutions

Endo Ice Spray

Coltene/Whaledent

830-H05032\$24.95 6 oz./can



MARK 3™

BEST VALUE

Pulp Vitality Endo Refrigerant Spray

100-5970\$16.49

Contains: 6 oz can of Spray plus 1 Dispensing Tube

- Quick and effective way to check vitality
- Directional nozzle for easy placement
- Colder than ethyl chloride or ice



COMPARE TO
Endo Ice

RC Prep

Premier Dental

580-9007131\$30.49

18g jar

580-9007133\$249.95

227g jar

580-9007127\$13.95

Plastic Syringe Tip Refill (50)



580-9007139\$17.95

Metal Syringe Tip Refill (50) - (use with 3cc Syringe)

580-9007135\$38.95

9g Syringe 2/pk

580-9007129\$36.49

Syringe Kit 3cc

Contains: 5 - 3cc Syringes and 50 Tips

Prep-Rite RC

Puldent

590-PRC.....\$30.49

5g Syringes, 4/pk



EDTA Solution

Puldent

590-EDTA-30.....\$10.75 30ml

590-EDTA-60.....\$16.75 60ml

590-EDTA-120\$21.95 120ml



Multi-Cal

Puldent

590-MULTI\$42.49

Contains: 4 - 1.2ml Syringes and 8 Applicator Tips



Pulp Capping Paste

Puldent

590-PSYK.....\$40.95

Complete Kit - 1 - 3ml and 24 Needles

590-PSYN.....\$14.75

Refill Sterile Needles - 24/bx

590-PSY\$25.49

Refill - 1 - 3ml Syringe



Canal Plus

Septodont

290-E0275\$32.49

Contains: 1 - 5g Syringe, 1 Stopper and 5 Tips

Needles

MARK 3™

Endo Irrigation Combo Needles


100-1700.....23 Gauge, 3cc - Blue

100-1701.....27 Gauge, 3cc - Yellow

\$28.95

Kit Contains:

100 - 3cc luer-lock endo irrigation syringes and 100 - 1" bendable monoject style (notched) needle tips



Rotary Files

- Made of premium Super-Elastic USA
- Constant Taper design available in 04 and 06 taper and ISO tip sizes 15 - 50
- 3 fluted file design
- Ideal for simple straight canals, calcified canals and removing gutta-percha
- Calibration rings for setting accurate working length

Pac-File Flex Rotary Files

Pac-Dent

Taper	Tip	Pac-File™ Conform™		Pac-File™ NiTi™	
		21mm	25mm	21mm	25mm
.04	15	108-FFC210415	108-FFC250415	108-FNT210415	108-FNT250415
.04	20	108-FFC210420	108-FFC250420	108-FNT210420	108-FNT250420
.04	25	108-FFC210425	108-FFC250425	108-FNT210425	108-FNT250425
.04	30	108-FFC210430	108-FFC250430	108-FNT210430	108-FNT250430
.04	35	108-FFC210435	108-FFC250435	108-FNT210435	108-FNT250435
.04	40	108-FFC210440	108-FFC250440	108-FNT210440	108-FNT250440
.04	45	108-FFC210445	108-FFC250445	108-FNT210445	108-FNT250445
.04	50	108-FFC210450	108-FFC250450	108-FNT210450	108-FNT250450
.04	55	108-FFC210455	108-FFC250455	108-FNT210455	108-FNT250455
.04	60	108-FFC210460	108-FFC250460	108-FNT210460	108-FNT250460
.04	15-40 assorted	108-FFC2104AS	108-FFC2504AS	108-FNT2104AS	108-FNT2504AS
.06	15	108-FFC210615	108-FFC250615	108-FNT210615	108-FNT250615
.06	20	108-FFC210620	108-FFC250620	108-FNT210620	108-FNT250620
.06	25	108-FFC210625	108-FFC250625	108-FNT210625	108-FNT250625
.06	30	108-FFC210630	108-FFC250630	108-FNT210630	108-FNT250630
.06	35	108-FFC210635	108-FFC250635	108-FNT210635	108-FNT250635
.06	40	108-FFC210640	108-FFC250640	108-FNT210640	108-FNT250640
.06	45	108-FFC210645	108-FFC250645	108-FNT210645	108-FNT250645
.06	50	108-FFC210650	108-FFC250650	108-FNT210650	108-FNT250650
.06	15-40 assorted	108-FFC2106AS	108-FFC2506AS	108-FNT2106AS	108-FNT2506AS



Pac-File Conform™

21mm \$38.50 6/pk

25mm \$38.50 6/pk

Pac-File NiTi™

21mm \$36.99 6/pk

25mm \$36.99 6/pk

Compare NiTi to K3™, Vortex®, HyFlex®NT™, Endosequence®


Compare Conform to TF® Twisted Files, Vortex Blue®, HyFlex®CM™

Features of Pac-File Conform™ Rotary Files



Can be pre-bent

- No rebound
- Superior resistance against breakage caused by cyclical fatigue
- Excellent canal tracking




Pac-Taper Flex Rotary Files

Pac-Dent

- Made of premium Super-Elastic nickel-titanium
- Ideal for simple straight canals, calcified canals and removing gutta-percha
- 3 fluted file design
- Calibration rings for setting accurate working length

Pac-Taper™		PROTAPER™ UNIVERSAL		Pac-Taper Conform™		Pac-Taper Flex™	
				21mm	25mm	21mm	25mm
CX	SX	108-TFC21CX	108-TFC25CX	108-TNT21CX	108-TNT25CX		
C1	S1	108-TFC21C1	108-TFC25C1	108-TNT21C1	108-TNT25C1		
C2	S2	108-TFC21C2	108-TFC25C2	108-TNT21C2	108-TNT25C2		
A1	S1	108-TFC21A1	108-TFC25A1	108-TNT21A1	108-TNT25A1		
A2	S2	108-TFC21A2	108-TFC25A2	108-TNT21A2	108-TNT25A2		
A3	S3	108-TFC21A3	108-TFC25A3	108-TNT21A3	108-TNT25A3		
A4	S4	108-TFC21A4	108-TFC25A4	108-TNT21A4	108-TNT25A4		
A5	S5	108-TFC21A5	108-TFC25A5	108-TNT21A5	108-TNT25A5		
CX, C1, C2, A1, A2, A3 (assorted)	SX, S1, S2, S1, S2, S3 (assorted)	108-TFC21AS	108-TFC25AS	108-TNT21AS	108-TNT25AS		



Pac-Taper Conform™

21mm \$38.50 6/pk

25mm \$38.50 6/pk

Compare to Pro Taper Gold®

Pac-Taper Flex™

21mm \$36.99 6/pk

25mm \$36.99 6/pk

Compare to Pro Taper®



Features of Pac-Taper Conform™ Rotary Files

Can be pre-bent

- No rebound
- Superior resistance against breakage caused by cyclical fatigue
- Excellent canal tracking



Sealers

Sealapex

Kerr

E430-18432\$54.95

Contains:

12g Base, 18g Catalyst and Mixing Pad



Tubliseal

Kerr

E430-00340\$52.95

Contains: 15g Tube Base,

5g Tube Accelerator and Mixing Pad



Apex Locators

RootPro™ Wireless Endo Motor with Apex Locator

Pac-Dent

Canal Preparation and Apex Locator System

- Automatic calibration ensures accuracy & eliminates the effect of changes in temperature and moisture
- Prepares root canals while displaying accurate & precise measurements of the root canal
- Many automatic safety control options:
 - Auto Motor Start/Stop
 - Auto Torque Reverse
 - Auto Apical Reverse
- Adjustable between 100 and 650 RPM
- Torque can be set between 0.6 to 5.2 N.cm



Wireless Foot Switch

Complete RootPro Wireless Endo Motor Unit Includes:

108-RP01M	1 x Complete RootPro™ Wireless Endo Motor Unit
108-FS01	1 x Wireless Foot Switch
108-RP-HP01	1 x Low Speed Handpiece Replacement
108-RAC240	1 x AC Adapter (110V-240, 50/60Hz)
108-LCH250	Low Sped Contra-Angle Handpiece Sleeves, 250/bx

\$1099.99

BEST VALUE

Root ZX II Apex Locator

J. Morita

860-5336392..... **\$979.99**

Contains: 1 Unit, 1 Probe, 3 File Holders, 5 Contrary Electrodes, 3 AA Batteries and 1 Tester

Slow speed handpiece sold separately



Accessories

Endo Aid Kit & Refills

Pac-Dent

- Made of autoclavable materials
- One-finger-held organizer with ruler ensures convenience
- Effectively organizes endodontic files



108-802	Endo Aid Kit: 1 x Endo Organizer, 50 Teal Foam Refills	\$36.99
108-803	50 - Triangle Foam Refills, Gray	\$15.49
108-803T	50- Triangle Foam Refills, Teal	\$15.49
108-803L	50- Triangle Foam Refills, Lavender	\$15.49
108-803R	Endo Organizer Only	\$19.49
108-804	50 - Regular Foam Refills, Gray	\$15.49

PacEndo™ Air/Water Irrigator

Pac-Dent

- Multi-purpose
- Works with most luer-lock tips
- Compatible with most air/water dental syringes
- Provides excellent vision when locating hidden canals
- Ideal for rinsing and drying during microsurgical procedures
- Use to air-dry the canals prior to obturation for optimal seal
- Remove debris when retrieving broken instruments



108-PEI-1	1 x PacEndo™ Air/Water Irrigator, Regular Size, 2"	\$76.99
108-PEI-2	1 x PacEndo™ Air/Water Irrigator, Extended Size, 3.5"	\$76.99

FREE SHIPPING NATIONWIDE

With orders over \$300
Excludes Canada. Certain restrictions apply.



Essence Complete Dental Package

Custom Colors and Designs Available

Standard Features:

- Factory new with 5 year warranty
- US made and assembled with all US parts
- 2 HP (3 HP option) delivery system
- Assistant's arm HVE/SE - Optional cuspidor
- HVE/SE Triple Syringe with autoclavable tip

Options:

- Cuspidor
- Touch Pad
- UltraLeather
- HVE
- Extra Triple Syringe
- Doctor Stool
- Assistant Stool

Call for more information:
866.333.6825



MADE IN THE USA



Allure Complete Dental Package

Custom Colors and Designs Available

Standard Features:

Factory new with 5 year warranty

US made and assembled with all US parts

2 HP (3 HP option) delivery system

Assistant's arm HVE/SE - Optional cuspidor

HVE/SE Triple Syringe with autoclavable tip

Options:

Cuspidor

Stools

Call for more
information:
866.333.6825



MADE IN THE USA





II Island Unit

Custom Colors and Designs Available

Standard Features:

- Pass-thru upper storage
- 3 drawers on Assistant's side
- 1 drawer on Doctor's side
- Handle pull doors and drawers
- X-ray mounting column
- 100" w x 82" h x 22" d
(add 3" for countertop)

Options:

- Trash drop
- Granite countertop
- Drawers/shelves option
- Round sink
- Glass sink

Call for more
information:
866.333.6825



MADE IN THE USA





II Sterilization Center

Custom Colors and Designs Available

Standard Features:

- Red/Blue lights
- Stainless steel sink
- Two counter waste drop (Regular & Hazardous waste)
- Slide out shelf for sterilizer
- Slide out shelf for statim
- 106" w x 80" h x 27" d
(add 3" for countertop)
- Shipped in 3 pieces

Options:

- Recessed cut-out for ultrasonic
- Tapmaster foot pedal
- Automatic soap dispenser
- Lights shown in photo
- Shelves

Call for more
information:
866.333.6825



MADE IN THE USA





12 O'Clock Rear Cabinet

100% Plywood (No particle or laminate wood used)

Standard Features:

Signature Series 5069 Radius delivery system with Assistant package

3 HP delivery with HVE/SE/Triple Syringe

Built-in Signature Series control panel

Built-in Belmed or Porter nitrous panel

Solid surface countertop

Solid surface round workstation

Swival tray installed

Ultra sliding computer bracket

Sliding keyboard for
implant motor

Call for more
information:
866.333.6825



MADE IN THE USA



MADE IN THE USA

5
YEAR
WARRANTY

Alabama Dental Mobile Rolling Assistant's Cart

Custom Colors and Designs Available

Standard Features:

21 1/2" w x 18 3/4" d x 32" h

Specify requirements for a front-to-rear or side-to-side slide top

Snow white standard color
other colors available

No assembly required

Call for more information: 866.333.6825


Buyer pays shipping

Options:

Choice of accent colors

Choice of countertops

Custom sizes and
designs available



Assistant's North Carolina Dental Mobile Rolling Cart

Custom Colors and Designs Available

Standard Features:

21 1/2" w x 18 3/4" d x 32" h

Specify left or right sliding side-shelf

Snow white standard color
other colors available

No assembly required

Call for more information: 866.333.6825

Buyer pays shipping

Options:

Choice of accent colors

Choice of countertops

Custom sizes and
designs available



Mobile Dental Portable Cabinet

Custom Colors and Designs Available

Standard Features:

21 1/2" w x 18 3/4" d x 32" h

Snow white standard color
other colors available

Seven drawers

No assembly required

Call for more information: 866.333.6825

Buyer pays shipping

Options:

Choice of accent colors

Choice of countertops

Custom sizes and
designs available



DENTAL EQUIPMENT SOLUTIONS



Bring your vision to life

Whether you're opening a new office or updating an existing one, we can help make all your dreams come true by delivering custom solutions and designs.

The Signature Series equipment line was developed by practicing dentists and dental engineers to ensure you are getting the most innovative, unique and reliable products on the market.

**Chairs • Delivery Systems • Operator Lights • Cabinetry
Mobile Carts • Compressors • Vacuums and More**

Questions? Email us: info@nobledentalsupplies.com




Aspirator Tips

MARK 3™

Surgical Tips

- 100-2599Blue - Small
- 100-2601White - Regular
- 100-2600Green - Large



\$5.95 25/bg

- Surgical Aspirator tips for suction
- Angled 30 degree tip

COMPARE TO
Surg-O-Vac

Evacuation Traps

Filters for Central Vacuum Pump

- 100-12200 2.75" x 3.625" - 12/pk
- 100-12300 3.5" x 4.375" - 8/pk
- 100-12350 3.5" x 4.375" - 8/pk Fine Screen

\$49.95 pk

Each kit is complete with:

A Waste Canister and Lid, Filter Bowl, Filter Screen, Replacement Rubber Gasket and Hazardous Warning Sticker



COMPARE TO
Evac-u-trap
&
Dispos-a-trap

In-Line Filters

- | | |
|-----------------------|-----------------------|
| 100-5500 1.25" | 100-5509 1.5" |
| 100-5501 2.125" | 100-5512 2.75" |
| 100-5503 1.875" | \$36.95 144/pk |
| 100-5505 3.75" | |
- \$36.95** 144/pk

- 100-6100 Universal - **\$47.95** 144/pk
- 100-6200 Universal - **\$47.95** 144/pk
- 100-6500 Universal - **\$36.95** 144/pk

Each kit is complete with:

A Waste Canister and Lid, Filter Bowl, Filter Screen, Replacement Rubber Gasket and Hazardous Warning Sticker



Evacuation System Cleaners

NeutraVAC

Biotrol

- 901-NV800 **\$53.95** 32 oz./btl
- Metered Dose Bottle



Vacusol Ultra

Biotrol

- 901-ED900 **\$49.75** 32 oz./btl
- Solution - Metered Dose Bottle



Evac Shock Evacuation Cleaner

Cory Labs

- 205-EVST64 **\$39.95** Tablets - 64/bx
- 205-5233 **\$32.95** Liquid - 32 oz.

- High level of concentrated formula - 64 uses
- Safe for all type of traps that require neutral ph 6.5 - 7.5
- Cleans cuspidors, pipe lines (particles, soils)
- Removes sludge build up
- Increases system efficiency and life of pump
- Compact, fast and easy to use
- Fresh lemon scent



**FAST
AND
EASY**

BEST VALUE

Vac Attak Evacuation System Cleaner

Premier Dental

- 580-9011101 **\$54.95** 800g Jar




Evacuator Tips

High Volume Vented Evacuator Tips

Unipak

- 706-UST-90712 .. **\$2.99** 100/bg



Saliva Ejectors

Unipak

- 791-BBS15 White
- 791-2900 Clear **\$2.99** 100/bg



Abrasives & Polishers

MARK 3™

Finishing & Polishing System

100-6340Intro Kit - 10x Point, 10x Disc, 10x cup

100-6341Cups

100-6342Disc

100-6343Point

\$49.99 30/bx

- Disposable - One Step
- Unique silicone rubber impregnated with aluminum oxide particles is 100% latex free
- Delivers perfect control to reduce surface roughness and bacterial retention
- Reinforced mandrel allows for optimal concentric torque



Poli-Pro Disks Composite Polishing System

Premier Dental

580-2019070**\$89.95** System Kit

Contains:

180 disks (20 of each size) and 4 RA mandrels

580-2019099**\$21.95** Starter Kit

Contains:

35 disks (5 each of 7 disk types) and 1 RA mandrel

Refills - 40 per Box **\$18.49**

Specify:

Coarse - Black

580-20190711/2" Grit Up

580-20190723/8" Grit Up

Medium - Blue

580-20190731/2" Grit Up

580-20190741/2" Grit Down

580-20190753/8" Grit Up

Fine - Red

580-20190761/2" Grit Double Sided

580-20190773/8" Grit Double Sided

Xtra Fine - White

580-20190781/2" Grit Double Sided

580-20190793/8" Grit Double Sided



Compo Strips

Premier Dental

580-2013127**\$30.95** Assorted Package

Contains:

2 Trimming (60 Micron Grit),
2 Finishing (45 Micron Grit),
and 2 Ultrafinishing (20 Micron Grit)



2.5mm Width **\$29.95** 6/bx

Specify:

580-2013107100T Blue (60 Micron Grit)

580-2013111100F Red (45 Micron Grit)

580-2013115100UF Yellow (20 Micron Grit)

3.75mm Width **\$29.95** 6/bx

Specify:


580-2013109150T Blue (60 Micron Grit)

580-2013113150F Red (45 Micron Grit)

580-2013117150UF Yellow (20 Micron Grit)

Sof-Lex Discs

3M/ESPE



1/2" and 3/8" Regular

1/2" and 3/8" Extra Thin

E760-1981SF Superfine 3/8"

E760-2381SF Superfine XT 3/8"

E760-1981F Fine 3/8"

E760-2381F Fine XT 3/8"

E760-1981M Medium 3/8"

E760-2381M Medium XT 3/8"

E760-1981C Coarse 3/8"

E760-2381C Coarse XT 3/8"

E760-1982SF Superfine 1/2"

E760-2382SF Superfine XT 1/2"

E760-1982F Fine 1/2"

E760-2382F Fine XT 1/2"

E760-1982M Medium 1/2"

E760-2382M Medium XT 1/2"

E760-1982C Coarse 1/2"

E760-2382C Coarse XT 1/2"

\$41.95 85/pk

MARK 3™

Pop-On Mandrels

100-1983RA**\$14.49** 3/pk

Right Angle



Composite Finishing

CARGUS

White Stones

100-1105.....FG CN1 - Conical

100-1155.....FG FL2 - Flame

100-9955.....FG RD1 - Round

\$16.75 12/bx

- White Arkansas stone (Aluminum Oxide) are specially designed for polishing composite and porcelain
- Very useful for removing excess of composite and polishing teeth

COMPARE TO
Shofu



Dura-White Stones

Shofu Dental Corp.

640-0241.....FG CN1 - Conical


640-0244.....FG FL2 - Flame

640-0247.....FG RD1 - Round

640-0248.....FG RD2 - Round

640-0249.....FG TC1 - Tapered Cone

\$24.49 12/bx



640-0201.....HP CN1 - Conical

640-0202.....HP FL1 - Flame

640-0203.....HP FL2 - Flame

640-0204.....HP FL3 - Flame

640-0205.....HP IC4 - Inverted Cone


640-0207.....HP RD1 - Round

640-0209.....HP RE1 - Round Wheel

640-0210.....HP TC1 - Tapered Cone

640-0211.....HP TC4 - Tapered Cone

\$24.49 12/bx



640-0221.....RA CN1 - Conical

640-0223.....RA FL2 - Flame

640-0224.....RA FL3 - Flame

640-0227.....RA RD1 - Round

\$24.49 12/bx

OneGloss

Shofu Dental Corp.

PS Kit Assortment

640-0175.....**\$94.95** 50/bx

PS Refill Cup

640-0176.....**\$94.95** 50/bx

PS Refill Midi

640-0177.....**\$94.95** 50/bx

PS Refill IC

640-0178.....**\$94.95** 50/bx



Polishing Pastes

Luminescence

Premier Dental

580-2014202.....**\$43.95** Gel - 3g Syringe

580-2014206.....**\$34.95** Applicator Tips - 50/pk



Diamond Twist SCL Extra-Oral Polishing Kit

Premier Dental

580-2019000.....**\$124.95**

Complete Kit

Contains:

- 1 - 6g Syringe Diamond Twist Polishing Paste, 2 HP Mandrels, 4 Prepolishing Fibra Shapes and 4 Felt Paste Applicators



580-2019001.....**\$57.95** Refill

6g Diamond Twist SCL Polishing Paste (Unflavored)


Strips

Lightening Strips - Metal

Miltex

540-620-81090.....**\$14.95** 10/tube

Medium Width - Medium Grit



IMPRESSION MATERIAL

Alginates

MARK3™

BEST VALUE

Alginate Impression Material

100-9020Fast Set

Dust Free - Pink **\$9.95** 1.1 lb/bg

- Patient friendly setting times
Normal Set 1½ minute - Fast Set 1 minute
- Firm consistency
- More rigid impression
- Good tear and deformation resistance
- Detail reproduction at least 50µm
- Peppermint aroma
- Compatible with all types of gypsum, smooth gypsum surface



COMPARE TO
Jeltrate Plus
Fast Set

Alginate Alternative

100-9010Fast Set

\$57.95 8, 50ml Cartridges/bx

- Patient friendly setting times
Work time 1 minute - Intra Oral set time 1 minute
- Pre-mixed cartridges.
No mess - Saves time
- Excellent dimensional stability
- Pleasant mint flavor
- Cost effective alternative



COMPARE TO
Alginit

Supergel

Bosworth

230-0921825Regular Set

230-0921826Fast Set

\$14.95 1 lb/can



Supergel Fresh

Bosworth

230-0921822Regular Set

230-0921822Fast Set

\$15.49 1 lb/can



Kromopan Alginate

Lascod

609-KRM302**\$10.95** 1 lb/bg



MARK3™
THE SMART ALTERNATIVE



Vinyl Polysiloxane

Bite Registration Material

FEATURES & BENEFITS

- Mousse-like consistency allows exceptional accuracy and greater patient comfort
- Easier to extrude - less force required
- Super Fast, Fast and Regular Set working times
- Plaster-like hardness
- Made in the USA

BEST VALUE



Standard Unflavored \$22.95 ea.

100-0218Regular Set

100-0219Fast Set

100-0220Super Fast Set

Hard Bite Unflavored \$25.95 ea.

100-0221Fast Set

Clear Bite Unflavored \$28.95 ea.

100-0222Clear Bite Set

Mint Chocolate Chip \$22.95 ea.

100-0224Fast Set

100-0225Super Fast Set

COMPARE TO
Blu-Mousse

Contains 2, 50ml cartridges plus 6 mixing tips per box

Bite Registration Trays


MARK 3™

Premium Bite Registration Trays

- 100-0330**.....Anterior- 35/pk
 - 100-0331**.....Posterior - 50/pk
 - 100-0332**.....Sideless Posterior - 50/pk
 - 100-0333**.....Wide Body Posterior - 50/pk
 - 100-0334**.....Quadrant -35/pk
- \$19.95**

- Ergonomic handles for easy seating and removal
- Anatomical design for maximum comfort
- Consistent, high quality impressions
- Thin, non-tear mesh for superior strength and maximization of intercuspation of teeth ensuring accurate bite

COMPARE TO
**Premier
Triple Trays**



Bite Registration Trays

- 100-0201**.....Anterior- 35/bx
 - 100-0202**.....Quadrant - 35/bx
 - 100-0203**.....Posterior - 50/bx
 - 100-0204**.....Sideless Posterior - 50/bx
 - 100-0205**.....Full Arch -30/bg
- \$15.95**

- Disposable impression trays for counter impressions and bite registration
- Rigid plastic frame
- Thin non-tear mesh
- Ergonomically designed handles

SAMPLES
FREE
available



Bite Registration

Exabite II NDS

GC America

250-132412**\$41.95** 2/pk Refill
Contains:
2 - 50ml Cartridges and 12 Mixing Tips

250-132418**\$149.95** 8/pk Refill
Contains:
8 - 50ml Cartridges and 48 Mixing Tips



Triple Tray

Premier Dental

- 580-1006213**Anterior - 35/bx
- 580-1006203**Posterior - 48/bx
- 580-1006223**Quadrant - 40/bx
- 580-1006233**Assorted - 35/bx
- 580-1006243**Sideless Posterior - 48/bx
- 580-1006253**Extended Quadrant - 35/bx
- 580-1006273**3/4 Arch - 35/bx
- 580-1006283**Posterior Wide Body - 48/bx

\$43.95

T-Loc Triple Tray

- 580-1006204**T-Loc Anterior - 35/bx
- 580-1006214**T-Loc Posterior - 28/bx

\$38.95



Alpha Triple Trays

Premier Dental

- 580-1006313**Small - (Lavender)
- 580-1006323**Medium - White
- 580-1006333**Large - Green
- 580-1006343**Assorted

\$34.49 24/bx

Assorted Contains: 5 Small, 14 Medium and 5 Large Trays



Neo Trays

Premier Dental

- 580-2201NT**Posterior - 48/bx
- 580-2202NT**Sideless - 48/bx
- 580-2203NT**Quadrant - 40/bx
- 580-2204NT**Anterior - 30/bx
- 580-2205NT**Full Arch - 30/bx
- 580-2207NT**Extended Quadrant - 35/bx
- 580-2208NT**Posterior Wide - 48/bx
- 580-2206NT**Introductory Kit - 40/bx

\$45.95

Introductory Kit Contains: 10 Posterior, 10 Sidless, 10 Quadrant, 5 Anterior and 5 Full Arch



Impression Trays

MARK 3™

Disposable Impression Trays

- 100-1001.....#1 Upper Large
 - 100-1002.....#2 Lower Large
 - 100-1003.....#3 Upper Medium
 - 100-1004.....#4 Lower Medium
 - 100-1005.....#5 Upper Small
 - 100-1006.....#6 Lower Small
 - 100-1007.....#7 UL-LR
 - 100-1008.....#8 UR -LL
 - 100-1009.....#9 Anterior Upper
 - 100-1010.....#10 Anterior Lower
- \$5.99** 12/bg



SAMPLES
FREE
available

COMPARE TO
Vista Trays

- Strong, rigid, distortion free construction
- Anatomically designed
- Smooth edges for comfort
- Perforations are designed to enhance consistent flow and retention
- 100% latex free
- Highly impact resistant
- Compatible with a wide variety of impression materials

Tray-Aways

Bosworth

Perforated/Disposable

- 230-0921884.....#1 Large Upper
 - 230-0921887.....#2 Large Lower
 - 230-0921885.....#3 Medium Upper
 - 230-0921888.....#4 Medium Lower
 - 230-0921886.....#5 Small Upper
 - 230-0921889.....#6 Small Lower
 - 230-0921894.....#7S Upper Right/Lower Left
 - 230-0921895.....#8S Upper Left/Lower Right
 - 230-0921896.....#9S Lower Anterior
 - 230-0921897.....#10S Upper Anterior
- \$9.49** 12/bg



Tra-Ten Trays

Water Pik, Inc

Perforated/Disposable

- 710-013451.....#1 Upper - Small
 - 710-013452.....#2 Upper - Medium
 - 710-013453.....#3 Upper - Large
 - 710-013454.....#4 Lower - Small
 - 710-013455.....#5 Lower - Medium
 - 710-013456.....#6 Lower - Large
 - 710-013457.....#7 Anterior
 - 710-013458.....#8 Bi-Lateral - UL/LR
 - 710-013459.....#9 Bi-Lateral - UR/LL
- \$10.24** 12/bg



High Performance Mixing Tips

MARK 3™

Short HP Mixing Tips

- 100-1430.....Green 6.5 mm
 - 100-1432.....Yellow 4.2 mm
- \$28.95** 48/pk



Compcore AF Mixing Tips

Premier Dental

580-3001417.....**\$51.95**

Contains: 35 Mixing Tips & 35 Intraoral Tips

Syringe Mix Mixing Tips

Premier Dental

580-3001418.....**\$53.49** 70/pk

Contains: 35 Mixing Tips & 35 Intraoral Tips

MARK 3™

Dynamic (Penta) Mixer Red

100-1414.....\$31.95 50/pk

- Dynamic mixer compatible with all Penta machines
- Comparable to 3M's and Penta™ Mixing tips
- Compatible with most dynamic machines
- Red tips are designed for Penta™ Red



COMPARE TO
Penta Red

Dynamic Yellow Mixing Tips

100-1419.....\$31.95 50/pk

- Fit on any automatic dynamic mixing units for use with 380ml VPS Impression Material cartridges



Bayonet Replacement Ring

100-1420.....\$6.99 2/pk

- For Automatic Mixing Cartridges. Used for 380ml VPS Cartridges



HP Mixing Tips

100-1418.....Blue 1:1 \$23.95 25/pk

100-1412.....Yellow 4.2 mm

100-1413.....Green 6.5 mm

100-1416.....Purple 7.5 mm

100-1417.....Pink 5.0 mm

\$27.95 48/pk



Intra Oral Mixing Tips

100-0059.....Clear

100-0060.....Yellow

\$16.95 100/pk

- Innovative color-coding system
- Allows impression materials to flow easily while keeping the materials bubble free.



Temporary Crown & Bridge Mixing Tips

100-110351.....\$21.95 25/bg

- Has clear tip with an orange mixing device inside
- Has a light blue nozzle and is 3" in length

Works with MARK3 Crown & Bridge For 10:1 Cartridges



Mixing Tools

MARK 3™

Mixing Bowls

100-1520.....Extra Small, 160ml - \$3.99

100-1521.....Small, 320ml - \$4.25

100-1522.....Medium, 500ml - \$4.75

100-1523.....Large, 700ml - \$5.25

- Flexible silicone mixing bowls
- Autoclavable up to 275°F



Mixing Pads

100-8584.....3"x 3" - \$2.99 100/pad


100-8587.....6"x 6" - \$4.49 100/pad



Mixing Wells

100-0002.....2 well - \$16.95 200/pk

100-0004.....4 well - \$17.95 200/pk



Polyether

Impregum F

3M/ESPE

E760-31473.....\$28.95

Catalyst - 15ml/tube

E760-31710.....\$94.95

Double Package

Contains: 4 tubes (2 each of 120ml base and 15ml catalyst) and mixing pad



Impregum Penta

3M/ESPE

E760-31644R.....Medium Body

E760-31730R.....Soft Body

\$229.95

Contains:

2-300ml base paste and 2-60ml catalyst paste



Permadyne Garant 2:1

3M/ESPE

E760-30412.....\$239.95

Refill Pack Contains:

4 - 50ml Cartridges and 20 Garant Mixing Tips



IMPRESSION MATERIAL

Tray Adhesives

TAC Adhesive Spray

Bosworth

230-0921880**\$25.49** 3.5 oz./can



Universal VPS Adhesive

GC America

250-138912**\$13.49** 7ml

MARK 3™

COMPARE TO
Caulk

VPS Tray Adhesive

100-3000**\$13.99**
10ml/btl

- Easy to apply
- Applicator brush top included
- Works with any VPS material
- Forms a strong bond between VPS impression material and tray



Hold Adhesive

Water-Pik

710-011460**\$20.95** 2/pk
.5 oz. Liquid

710-011461**\$27.95** ea.
2.1 oz. Aerosol Spray



Tray Cleaner

Coe Tray Cleaner

GC America

250-250010**\$32.95** 1.27 lb/btl



Vinyl Polysiloxane

Examix NDS

GC America

Refill Package with Tips - **\$30.99**

Contains: 2 - 48ml Cartridges and 6 Mixing Tips

250-137016Regular - Blue

250-137106Injection - Pink

250-137206Heavy - Yellow

250-137406Monophase - Purple



Refill Package without Tips - **\$26.99**

Contains: 2 - 48ml Cartridges

250-137032Regular - Blue

250-137132Injection - Pink

250-137232Heavy - Yellow

250-137432Monophase - Purple

Exafast NDS

GC America

Refill Package with Tips - **\$29.95**

Contains: 2 - 48ml Cartridges and 6 Mixing Tips

250-137007Regular

250-137107Injection

250-137207Heavy

250-137407Monophase

250-137307Putty Package **\$87.99**



Genie VPS

Sultan Healthcare

Standard Package **\$31.95**

700-77610Rapid Set Light Body

700-77615Standard Set Light Body

700-77620Rapid Set Regular Body

700-77625Standard Set Regular Body

700-77630Rapid Set Heavy Body

700-77635Standard Set Heavy Body

Each Contains: 2 - 50ml Cartridges and 6 mixing tips



MARK 3™ THE SMART ALTERNATIVE

Dental **PRODUCT SHOPPER**
RECOMMENDED
2015



Vinyl Polysiloxane • Super Hydrophilic VPS Impression Material

FEATURES & BENEFITS

- High performance cartridge system - easier to extrude
- Excellent dimensional stability decreases distortions
- Fast and Regular Set working times
- No odor or irritation allowing greater patient comfort
- Pleasant mint flavor
- Made in the USA



BEST VALUE

- 100-3010Light Body - Regular Set
- 100-3011Light Body - Fast Set
- 100-3012Medium Body - Regular Set
- 100-3013Medium Body - Fast Set
- 100-3014Heavy Body - Regular Set
- 100-3015Heavy Body - Fast Set
- 100-3016Monophase - Regular set
- 100-3017Monophase - Fast Set

Contains 4, 50ml cartridges per box **\$31.95**

Available in Bulk: 100 x 50ml Cartridges **\$695.00**

COMPARE TO
**Take 1
and
Examix**

MARK 3™

VPS 380 Impression Material

100-3025.....Heavy Body

100-3026.....Monophase

Work Time - 2:30 minute

Set Time - 1 minute

\$139.95 2 x 380ml cartridges



- Same great VPS impression material in a new delivery system that provides a more economic way to take accurate impressions
- For use in any dynamic mixing machines
- Excellent dimensional stability decreases distortions
- No odor or irritation allowing greater patient comfort
- Pleasant mint flavor
- Made in the USA

COMPARE TO
**Genie
Magic Mix**

VPS Putty

BEST VALUE

100-3004.....Regular Set

Work Time - 2 minutes Set Time - 4 minutes

100-3005.....Fast Set

Work Time - 1 minute Set Time - 3 minutes

\$49.95

Contains: 300ml each Base and Catalyst,
2 Color-coded Scoops and Instructions

COMPARE TO
**7312
Express
Putty**

- Precise detail reproduction
- Good dimensional stability
- Fast and Regular Set working times
- Excellent elastic recovery
- Made in the USA



Express STD Putty 3M/ESPE

760-7312.....**\$89.95**

Putty Kit - Firm Set 610ml



Enzymatic Cleaners



PRO one Sonictab™
The ULTRA Bioenzymatic Tablet

Will degrade **all** complex organics and soils that traditional protease enzymatic formulations cannot.

- Biodegradable, non-chlorine, anti-rust, dual enzyme formula
- Solution digests soils that may be resistant to other detergents
- Eliminates the handling and storage of heavy liquids
- Non corrosive, neutral pH
- Two tablets per gallon of water
- Pleasant fresh mint fragrance

BEST VALUE

Cory Labs

205-POS80 **\$47.95** Contains 80 Tablets

Purity Waterline Cleaning Tablets

Cory Labs

205-PRT60..... **\$32.95**

750ml Tablets 60/bx



- Safe and effective procedure maintains clean dental unit waterlines. Natural preservatives allow continuous cleaning
- Eliminates the need to empty bottles, purge and air dry lines each night
- Effective, easy to use, orderless, tasteless, non-irritating and has no effect on bonding. Compatible with all dental unit systems, or self contained water systems
- Environmental friendly, no harsh chemicals

Stain Out Tarter & Stain Remover

Cory Labs

205-TSR40..... **\$32.99** 40/bx



- Stain Out is an ultrasonic cleaning tablet formulated for removal of tartar, stain, tobacco, calculus from dentures, bridges and all your orthodontic appliances
- Fast and Easy (only 5 minutes to clean)
- Highly concentrated table, Single dose tablet of Stain Out can be mixed with 5-6 oz of water(ultrasonic beaker) for an effective cleaning procedure
- Neutral pH formulation safe for Instruments

MARK3™

Ultrasonic Bioenzymatic Cleaner

100-7640 **\$38.99** 64/bx

- Easy to use microbial bio-enzymatic tablets create a powerful, fast and effective ultrasonic cleaning solution
- USA/CANADA classified as biological safety level 1
- Non foaming, non-chlorine, anti-rust formula
- Digests and removes blood, bone, tissue, food particles and other soils
- May also be used as an instrument pre-soak or evacuation system cleaner
- Fresh mint fragrance



COMPARE TO
**Defend
 Tablets &
 Vibrakleen**

Eyewear

Econo Wrap - Pro-Vision Safety Glasses

Palmero Health Care

565-3601C Clear Lens

565-3601G Grey Lens

\$6.49 1 pair



Face Masks

Procedural Earloop Masks

Crosstex

150-GCPPK..... Pink

150-GCPBL..... Blue

\$11.95 50/pk



Isofluid Earloop Masks

Crosstex

150-GCIBL Blue

150-GCIPK Pink

150-GCILV Lavender

\$10.95 50/pk



Ultra Sensitive Earloop Masks

Crosstex

150-GCFCXS White
\$14.95 50/bx



Premium Extra Soft Earloop Masks

UniPak

790-ULM-6181..... ASTM Level 1 Blue

790-ULM-6184..... ASTM Level 1 Lavender

790-ULM-6189..... ASTM Level 1 Pink
\$3.95 50/bx



790-ULM-6381..... ASTM Level 3 Blue
\$7.95 50/bx

- PFE Exceeds 99%
- High Breathability
- Extra Soft & Comfortable
- Latex & Fiberglass Free
- Fluid Resistant
- 3-Ply Pleated

Extra-Safe Earloop Masks

Valuemax

570-5330E-PGO Multicolor

570-5330E-LP Pink

570-5330E-PP Purple

570-5330E-SB Sky Blue

570-5330E-TE..... Teal
\$9.95 50/bx



Extra-Safe Sensitive Earloop Masks

Valuemax

570-5430E-PGO Multicolor

570-5430E-LP Pink

570-5430E-PP Purple

570-5430E-SB Sky Blue

570-5430E-TE..... Teal
\$9.95 50/bx



Gloves

Perform Nitrile Gloves

Aurelia

705-92395 X-Small

705-92396 Small

705-92397 Medium

705-92398 Large

Powder Free/Latex Free

\$12.49 200/bx



Protege Nitrile Exam Gloves

Aurelia

705-93995 X-Small

705-93996 Small

705-93997 Medium

705-93998 Large

705-93999 X-Large

Powder Free/Latex Free

\$7.99 100/bx



Vibrant Micro-Textured Latex Exam Gloves

Aurelia

705-98225 X-Small

705-98226 Small

705-98227 Medium

705-98228 Large

705-98229 X-Large

Powder Free **\$7.49** 100/bx



Distinct Latex Gloves

Aurelia

705-29225 X-Small

705-29226 Small

705-29227 Medium

705-29228 Large

Honeycomb textured Polymer-Coated

Powder Free **\$8.99** 100/bx



Touch 200 Nitrile Exam Gloves - 3.5 mil.

Emerald

706-5400..... X-Small

706-5401..... Small

706-5402..... Medium

706-5403..... Large

Latex Free/Powder Free

3.5 mil.

\$8.79 200/bx



BEST VALUE

Latex Exam Gloves - 5 mil.

Emerald

706-5990..... X-Small

706-5991..... Small

706-5992..... Medium

706-5993..... Large

Powder Free - 5 mil.

\$5.99 100/bx



BEST VALUE

Nitrile Powder Free Gloves

Emerald

706-5799..... X-Small

706-5800..... Small

706-5801..... Medium

706-5802..... Large

Powder Free - 3.5 mil.

\$6.49 100/bx



Jackets & Gowns

Jacket X-Safe - Knee Length

Valumax

570-3660BBL	Blueberry/Large	570-3660RBL	Raspberry/Large
570-3660BBM	Blueberry/Medium	570-3660RBM	Raspberry/Medium
570-3660BBS	Blueberry/Small	570-3660RBS	Raspberry/Small
570-3660BBXL	Blueberry/X-Large	570-3660RBL	Raspberry/X-Large
570-3660CBL	Ceil Blue/Large	570-3660SBL	Sky Blue/Large
570-3660CBM	Ceil Blue/Medium	570-3660SBM	Sky Blue/Medium
570-3660CBS	Ceil Blue/Small	570-3660SBS	Sky Blue/Small
570-3660CBXL	Ceil Blue/X-Large	570-3660SBL	Sky Blue/X-Large
570-3660LPL	Light Pink/Large	570-3660TEL	Teal/Large
570-3660LPM	Light Pink/Medium	570-3660TEM	Teal/Medium
570-3660LPS	Light Pink/Small	570-3660TES	Teal/Small
570-3660LPXL	Light Pink/X-Large	570-3660TEXL	Teal/X-Large
570-3660PPL	Purple/Large	570-3660WHL	White/Large
570-3660PPM	Purple/Medium	570-3660WHM	White/Medium
570-3660PPS	Purple/Small	570-3660WHS	White/Small
570-3660PPXL	Purple/X-Large	570-3660WHXL	White/X-Large

\$32.99 10/pk

- Soft, protective 3-layer SMS fabrics
- Fluid-resistant and breathable
- Latex-free & fire retardant
- Autoclavable & washable up to 5 times
- Meets OSHA requirements
- Anti-static treated - No clinging
- Knee length
- 3 Pockets (2 hip pockets; 1 breast pocket)
- Snap front
- Knit cuffs and collar



Jacket X-Safe - Hip Length

Valumax

570-3630WHXS White/X-Small	570-3630PPS Purple/Small
570-3630WHS White/Small	570-3630PPM Purple/Medium
570-3630WHM White/Medium	570-3630PPL Purple/Large
570-3630WHL White/Large	570-3630PPXL Purple/X-Large
570-3630WHXL White/X-Large	570-3630RBS Raspberry/Small
570-3630BBXS Blueberry/X-Small	570-3630RBM Raspberry/Medium
570-3630BBS Blueberry/Small	570-3630RBL Raspberry/Large
570-3630BBM Blueberry/Medium	570-3630RBXL Raspberry/X-Large
570-3630BBL Blueberry/Large	570-3630RB2XL Raspberry/2X-Large
570-3630BBXL Blueberry/X-Large	570-3630SBS Sky Blue/Small
570-3630CBS Ceil Blue/Small	570-3630SBM Sky Blue/Medium
570-3630CBM Ceil Blue/Medium	570-3630SBL Sky Blue/Large
570-3630CBL Ceil Blue/Large	570-3630SBXL Sky Blue/X-Large
570-3630CBXL Ceil Blue/X-Large	570-3630TES Teal/Small
570-3630CB2XL Ceil Blue/2X-Large	570-3630TEM Teal/Medium
570-3630BB3XL Ceil Blue/3X-Large	570-3630TEL Teal/Large
570-3630CYS Cranberry/Small	570-3630TEXL Teal/X-Large
570-3630CYM Cranberry/Medium	570-3630MBM SF Medical Blue/Med
570-3630CYL Cranberry/Large	570-3630PGOM Multi Color/Medium
570-3630CYXL Cranberry/X-Large	570-3630PGOL Multi Color/Large
570-3630LPXS Light Pink/X-Small	\$32.99 10/pk
570-3630LPS Light Pink/Small	
570-3630LPM Light Pink/Medium	
570-3630LPL Light Pink/Large	
570-3630LPXL Light Pink/X-Large	

- Soft, protective 3-layer SMS fabrics
- Fluid-resistant and breathable
- Autoclavable & washable up to 5 times
- Latex-free & fire retardant
- Meets OSHA requirements
- Anti-static – No clinging
- Hip length
- 3 Pockets
(2 hip pockets; 1 breast pocket)
- Snap front
- Knit cuffs and collar



Isolation Gown Elastic Cuffs

Valumax

570-3230B Blue/Regular
570-3230BLX Blue/X-Large
570-3230G Green/Regular
570-3230P Pink/Regular
570-3230W White/Regular
570-3230Y Yellow/Regular
\$49.95 50/cs

- 1-layer non-woven polypropylene fabric
- Light-weight
- Breathable
- Splash Resistant
- Latex-Free
- Tie At Back
- Elastic Cuff



Isolation Gown Knit Cuffs

Valumax

570-3260B Blue/Regular
570-3260BLX Blue/X-Large
570-3260G Green/Regular
570-3260P Pink/Regular
570-3260W White/Regular
570-3260Y Yellow/Regular
\$49.95 50/cs

- 1-layer non-woven fabric
- Breathable
- Splash Resistant
- Latex-Free
- Tie At Back
- Knit Cuff



Patient Bibs

MARK 3™

Dental Bibs 2-ply

- 100-1319B Blue
- 100-1319W White
- 100-1319PI Pink
- 100-1319LV Lavender
- 100-1319G Green
- 100-1319PE Peach
- 100-1319Y Yellow



2-ply, 13"x 18" **\$19.99** 500/bx

- Maximum protection with horizontal embossed pattern and a unique water repellent edge.
- Increased strength and tear resistance with plastic backing for added clothing protection.
- 2 ply paper plus 1 ply poly.

Sterilization Pouches

Self Seal Sterilization Pouches

UniPak

- 790-21072.25" x 5"**\$7.49**
- 790-21082.75" x 10"**\$7.39**
- 790-21092" x 9"**\$7.49**
- 790-21103.5" x 6.5"**\$8.25**
- 790-21113.5" x 10"**\$6.19**
- 790-21123.25" x 12"**\$9.75**
- 790-21134.25" x 12"**\$11.49**
- 790-21145.25" x 7.5"**\$9.49**
- 790-21155.25" x 11"**\$12.75**
- 790-21167.5" x 13"**\$24.95**
- 790-211710" x 15"**\$43.95**
- 790-211812" x 18"**\$59.95**

Each 200/bx



Surface Barriers

MARK 3™

Air Water Syringe Sleeves

100-2104**\$7.09** 500/pk
2.5" x 10" - Blue

100-2103**\$7.09** 500/bx
2.5" x 10" - Clear

- Prevents spreading of infection caused by the splashing and splattering of excess blood and saliva
- Provides patient with a clean and relaxing rinse and dry process
- Disposable
- Prevents cross contamination




Barrier Film

100-21014" x 6" - Blue

100-21024" x 6" - Clear
1200 perforated sheets/bx **\$14.49**

COMPARE TO
Allrap

- Extra Tacky Adhesive
- No Residue
- Disposable
- Seconds to Apply
- Non Stick Edges



Digital X-Ray Sensor Sleeves

100-7798Size 1 - Small/Medium 1 3/8" x 8"

100-7799Size 2 - Large 1 5/8" x 8"
Plastic Sleeves **\$12.95** 500/bx



Bite Block Covers

Palmero

565-1850**\$21.49** 1000/bx



Surface Disinfectants

Birex SE Disinfectant

Biotrol

901-BC144\$519.95

Clinic Pack - .125 oz. - 144/pk

901-BI012\$179.95

Super Refill Packets - .5 oz. - 12/pk

901-BI013\$189.95

Super Pack Kit -

Contains: 12.5 oz. Packets,
Gallon Mixing Container,
4 - 16 oz. Spray Bottles,
Funnel and Squirt Cap

901-BI025\$103.95

Intro Kit - 24 Pack

Contains: 24 - 0.125 oz. Packets,
Spray Bottle, Squirt Cap and Sprayer

901-BI048\$184.95

Operatory Pack - .125 oz. - 48/pk

901-SC103\$13.05

Spray Bottle - 32 oz. - ea.



Caviwipes Towelettes

Metrex

500-13-1100\$9.95

Large - 160/canister

500-13-1150\$12.49

Extra Large - 65/canister



Cavicide

Metrex

500-13-1000\$29.95 Gallon

500-13-1024\$12.95 24 oz./btl

500-13-1025\$62.95 2.5 Gallon



DisCide Ultra Disinfectant

Palermo Health Care

565-3565Q\$14.25 32 oz./pump spray

565-3565G\$32.49 Gallon



DisCide Ultra Towelettes

Palermo Health Care

565-60DIS\$10.75 160/canister
Regular Size - 6" x 6.75"

565-10DIS\$11.99 60/canister
Large Size - 10.5" x 10.5"



**WE'RE ALWAYS ADDING
NEW PRODUCTS**

Please Go To

www.Nobledentalsupplies.com

And find the latest updates.

or Email Us At:

info@nobledentalsupplies.com

J&J Instruments

DIAGNOSTIC

Explorers

854-04-811Explorer Double End #11/12 **\$6.99** ea.



854-04-820Explorer Double End #2 **\$6.99** ea.



854-04-850Explorer Double End #5 (23/17) **\$6.99** ea.



854-04-923Explorer Single End #23 **\$6.45** ea



854-04-996Explorer Double End #DG16 **\$6.99** ea.



Explorer/Probe

854-16-520Expro #23/CP11 Color Coded
3, 6, 8, 11mm **\$11.45** ea.



854-16-521Expro #23/CP12 Color Coded
3, 6, 9, 12mm **\$11.45** ea.



854-16-522Explorer Williams Color Coded Probe
Double End Expro #23 **\$11.45** ea.



854-16-523Williams Periodontal Expro #23 **\$9.45** ea.



Mirrors/Mirror Handles

854-07-210Front Surface C.S. #4 Mirror
Chromium Coating **\$16.45** 12/pk



854-07-230Front Surface C.S. #5 Mirror
Chromium Coating **\$16.45** 12/pk



DIAGNOSTIC cont.

854-07-400Mirror Handle Oct. Simple Stem **\$2.99** ea.



854-07-405Mirror Handle Oct. Cone Socket **\$2.99** ea.



854-07-410Mirror Handle Round Simple Stem **\$2.99** ea.



854-07-420Mirror Handle Round Cone Socket **\$2.99** ea.



Pliers

854-08-200College Plier #317 6" Serrated Tips **\$4.45** ea.



854-08-210College Plier 6" Self Locking Serrated **\$5.99** ea.



854-08-220College Plier 6" Self Locking Grooved **\$5.99** ea.



854-08-255Meriam Plier 6.25" Serrated Tips **\$5.99** ea.



Probes

854-16-512Periodontal Probe #12-3, 6, 9, 12mm **\$8.99** ea.



854-16-559Periodontal Probe #CP12 Color Coded
3, 6, 9, 12mm **\$8.99** ea.



ORTHODONTIC

Explorers

854-09-114Johnson Contouring Pliers #114 **\$22.45** ea.



854-09-200Clasp Adjusting Plier #200 **\$22.45** ea.



PERIODONTAL

Curettes

854-02-285Lucas Curette #85 **\$12.99** ea.



854-02-286Lucas Curette #86 **\$12.99** ea.



854-02-291Molt Bone Curette #1 **\$16.45** ea.



854-95-612Gracey Curette #1/2 SILK Handle **\$20.45** ea.



854-95-694Gracey Curette #11/12 SILK Handle **\$20.45** ea.



854-95-695Gracey Curette #13/14 SILK Handle **\$20.45** ea.



Scalers

854-11-067Scaler #H6/H7 **\$11.99** ea.



854-95-067Scaler #S6/S7 SILK Handle **\$21.09** ea.



854-95-204Scaler #204S SILK Handle **\$21.09** ea.



RESTORATIVE

854-01-040Amalgam Carrier Double Ended
Regular/Regular **\$19.99** ea.



854-01-215Black Plugger #0/1 Smooth **\$8.99** ea.



854-02-540Cavity Liner 4.25" Regular Tip **\$4.45** ea.



RESTORATIVE cont.

854-07-300 Universal Tofflemire Adult
Adult Matrix Retainer **\$10.45** ea.



854-08-300 Articulating Paper Forceps Straight
Serrated **\$7.99** ea.



854-11-110 Friedman Rongeur 5.5" 30 degrees **\$38.45** ea.



Burnishers

854-02-504 Burnisher #25/26 13A Double Ended **\$8.99** ea.



854-02-510 Burnisher Double Ended #27/29 **\$8.99** ea.



854-02-602 Acorn Burnisher #21D Double Ended **\$8.99** ea.



Carvers

854-03-230 Hollenback Carver #3S (1/2) **\$8.99** ea.



854-03-335 Interproximal Carver **\$8.99** ea.



Cleoid Discoid

854-04-603 Cleoid Discoid #3/6 Large **\$8.99** ea.



854-04-689 Cleoid Discoid #89/92 **\$8.99** ea.



Excavators

854-04-717 Excavator #17 **\$8.99** ea.



854-04-718 Excavator #18 **\$8.99** ea.



RESTORATIVE cont.

Excavators

854-04-719Excavator #19 **\$8.99** ea.



854-04-731Excavator #31L **\$8.99** ea.



854-04-738Excavator #38/39 **\$8.99** ea.



Filling Instruments

854-08-000Gregg Plastic Filling Instrument #4/5 **\$8.99** ea.



854-08-003Ladmore Plastic Filling Instrument #3 **\$8.99** ea.



854-08-040Plastic Filling Instrument #4 **\$8.99** ea.



854-08-050Plastic Filling Instrument #5 **\$8.99** ea.



854-08-060Plastic Filling Instrument #13 **\$8.99** ea.



854-08-101Woodson Plastic Filling Instrument #1 **\$8.99** ea.



854-95-020Plastic Filling Instrument #2
Silk Handle **\$20.45** ea.




Rubber Dam Punches

854-11-500Ainsworth Rubber Dam Punch **\$58.99** ea.



854-11-550IV Type Rubber Dam Punch **\$58.99** ea.



RESTORATIVE cont.

Spatulas

854-13-122Cement Spatula #22 **\$6.45** ea.



854-13-124Cement Spatula #24 **\$6.45** ea.



SURGICAL

854-02-010Coupland Bone Chisel #1 **\$13.45** ea.



854-07-505University of Minnesota
Cheek Retractor **\$11.45** ea.



854-07-800Aspirating Syringe
Cook Weight 1.8cc **\$19.99** ea.



854-11-030Blumenthal Rongeur 6"
30 Degrees **\$37.99** ea.



854-19-3000Scalpel Handle #3
Fits Blades 12B, 15, 15C **\$2.45** ea.



854-25-6420Crile - Wood Needle Holder **\$10.99** ea.



Bone Files

854-02-120Howard Bone File #12 **\$12.99** ea.



854-02-210Miller Bone File #21 **\$12.99** ea.



854-02-364Miller Bone File #64 **\$12.99** ea.



SURGICAL cont.

Elevators

854-04-120.....Molt Periosteal Elevators #9 **\$14.99** ea.



854-04-130.....Goldman-Fox Periosteal Elevator #14 **\$14.99** ea.



854-04-331.....Cryer Elevator #31 - Large **\$16.99** ea.



854-04-332.....Cryer Elevator #32 - Large **\$16.99** ea.



854-04-334.....Elevator Gouge #34 **\$16.99** ea.



854-04-401.....Apical Elevator #301 **\$16.99** ea.



854-04-403.....Apical Elevator #303 **\$16.99** ea.



854-04-404.....Apical Elevator #304 **\$16.99** ea.



854-04-591.....Potts Cross-Bar Elevator #1 **\$18.99** ea.



854-04-592.....Potts Cross-Bar Elevator #2 **\$18.99** ea.



854-04-925.....Cryer Elevator #25 - Medium **\$16.99** ea.



854-04-926.....Cryer Elevator #26 - Medium **\$16.99** ea.



854-04-944.....Cryer Elevator #44 - Small **\$16.99** ea.



SURGICAL cont.

854-04-945.....Cryer Elevator #45 **\$16.99** ea.



854-04-961.....Schmeckeber Elevator #81 **\$16.99** ea.



854-04-964.....Schmeckeber Elevator #4 **\$16.99** ea.



854-04-965.....Schmeckeber Elevator #5 **\$16.99** ea.



Forceps

854-05-001.....Extracting Forceps #1 Standard
Upper Incisor Cuspid **\$32.99** ea.



854-05-105.....Extracting Forceps #10S Upper Molar
Universal Straight Handle **\$32.99** ea.




854-05-170.....Extracting Forceps #17 Lower 1st & 2nd Molar
Universal Straight Handle **\$32.99** ea.



854-05-230.....Extracting Forceps #23 Lower 1st & 2nd Molar
Horn Beak Straight Handle **\$32.99** ea.



854-05-650.....Extracting Forceps #65 Upper Bicuspid/
Incisor Root Bayonet **\$32.99** ea.



854-05-880.....Extracting Forceps #88L Upper Left
1st & 2nd Molar **\$32.99** ea.



854-05-885.....Extracting Forceps #88R Upper Right
1st & 2nd Molar **\$32.99** ea.



SURGICAL cont.

Forceps

854-06-150.....Extracting Forceps Upper Incisor Bicuspid Root Universal **\$32.99** ea.



854-06-151.....Extracting Forceps Universal Bicuspid Root Lower **\$32.99** ea.



854-06-152.....Extracting Forceps Universal Incisor Bicuspid Upper **\$32.99** ea.



854-06-153....Extracting Forceps #151A Universal Incisor Bicuspid Lower **\$32.99** ea.



854-06-154....Extracting Forceps #150S Universal Incisor Bicuspid Lower Pediatric **\$32.99** ea.



854-06-155....Extracting Forceps Universal Incisor Bicuspid Lower **\$32.99** ea.




854-06-156....Extracting Forceps #150 Universal **\$32.99** ea.



854-06-157....Extracting Forceps #151 Universal **\$32.99** ea.



854-06-158....Extracting Forceps #150XS Universal Bicuspid Root Upper Pediatric with Spring **\$36.45** ea.



SURGICAL cont.

854-06-159.....Extracting Forceps #151XS Universal Incisor Bicuspid Lower Pediatric with Spring **\$36.45** ea.



854-06-286.....Extracting Forceps #286 Upper Bicuspid Incisor Root Bayonet **\$32.99** ea.



854-06-323.....Extracting Forceps #MD3 Universal Lower Incisor Cuspid Bicuspid Root **\$32.99** ea.



854-24-1130...Kelly Forceps 5.5" Straight **\$8.99** ea.



854-1135.....Kelly Forceps 5.5" Curved **\$8.99** ea.



854-31-1000...Backhouse Towel Forceps 3.5" **\$8.45** ea.



Hemostats

854-24-1050...Mosquito Hemostat Straight FCP 5" **\$8.45** ea.



854-24-1055..Mosquito Hemostat Curved FCP 5" **\$9.99** ea.



PLEASE CHECK OUR WEBSITE
For many more Instrument types and sizes

SURGICAL cont.

Molt Mouth Gag

854-07-900Molt Mouth Gag w/ Silicon Tips
Pediatric **\$47.45** ea.



854-07-910Molt Mouth Gag w/ Silicon Tips
Adult **\$48.99** ea.




Scissors


854-12-020Crown Scissor Curved 4.5"
\$5.99 ea.



854-12-215Kelly Scissor 6.25" Straight **\$8.99** ea.



854-12-235Kelly Scissor 6.25" Curved **\$9.45** ea.



854-22-1500...Iris Scissor 4.5" Straight **\$6.99** ea.



854-22-1510...Iris Scissor 4.5" Curved **\$7.45** ea.



Miscellaneous

Front Surface Mirrors

Miltex

Simple Stem

540-49001010 ...#4

540-49001011 ...#5

\$14.95 12/pk

Cone Socket

540-49001014 ...#4

540-49001015 ...#5

\$14.79 12/pk



The ULTRA bioenzymatic cleaner is also the best priced.

PRO one SONICTAB™



- Biodegradable, non-chlorine formula
- Solution digests soils that may be resistant to other detergents
- Eliminates the handling and storage of heavy liquids
- Anti-Hard water mineralization
- **Non corrosive, neutral ph**
- **Anti-Rust**

80 Tablets **\$47.95**

Compare to VibraKleen at \$94.50

Will degrade **all** complex organics and soils that traditional enzymatic formulations cannot.




Matrix Bands & Strips

T-Bands

Pulpdent

- 590-BTSS/NStainless Steel - Straight/Narrow
 - 590-BTSS/WStainless Steel - Straight/Wide
 - 590-BTSS/A.....Stainless Steel - Straight/Assorted
 - 590-BTSC/NStainless Steel - Curved/Narrow
 - 590-BTSC/WStainless Steel - Curved/Wide
 - 590-BTSC/AStainless Steel - Curved/Assorted
 - 590-BTBS/NBrass - Straight/Narrow
 - 590-BTBS/WBrass - Straight/Wide
 - 590-BTBS/ABrass - Straight/Assorted
 - 590-BTBC/NBrass - Curved/Narrow
 - 590-BTBC/WBrass - Curved/Wide
 - 590-BTBC/ABrass - Curved/Assorted
- 0.002" / 0.05mm thick **\$15.75** 100/pk



Tofflemire Bands

WaterPik Technologies Inc

Regular

- 710-34577-0.....#1 .002
 - 710-34578-9.....#1 .0015 Universal
 - 710-34580-0#2 .0015
 - 710-34585-0#13 .0015
- \$6.49** 36/pk

Dead Soft .001

- 710-34603-6#1 Adult Universal
 - 710-34604-4#2 Adult MOD Wide
 - 710-34605-2#13 Pedo Universal
- \$8.75** 36/pk
- 710-061100-00 ..Retainer Universal
- \$33.95** 1/pk



H0 Bands

Young Dental

Regular

- 770-350110#1 Universal
- 770-350210#2 Mod
- 770-350310#13 Pedo

Thin Tofflemire Matrix Bands - 0.001 Gauge

\$32.99 100/pk

Dead Soft

- 770-351110#1 Universal
 - 770-351210#2 Mod
 - 770-351310#13 Pedo
- \$33.95** 100/pk




Wedges

Sycamore Wedges

Premier Dental

- 580-906110013XT - Neutral
- 580-906110113T - Green
- 580-906110213 - Yellow
- 580-906110315 - Blue
- 580-906110416 - Red
- 580-906110517 - Neutral



- 580-9061106Ultra Thin - Orange
 - 580-9061108Assorted
 - 580-9061109Mini - Light Blue
- \$39.49** 400/bx

Assorted Pack Contains: 100 each of Ultra Thin, 13T, 13 and 15

Wizard Wedges

WaterPik Technologies Inc

- 710-061211Small
- 710-061212Large
- 710-061213Slim Jim
- 710-061214Wide
- 710-061215Assorted

\$38.95 500/bx



Posts

MARK 3™

Screw Posts - Gold Plated

100-0613 Short S1	100-0612 Medium M6
100-0614 Short S2	100-0601 Long L1
100-0615 Short S3	100-0602 Long L2
100-0616 Short S4	100-0603 Long L3
100-0617 Short S5	100-0604 Long L4
100-0618 Short S6	100-0605 Long L5
100-0607 Medium M1	100-0606 Long L6
100-0608 Medium M2	100-0619 Extra Long XL4
100-0609 Medium M3	100-0620 Extra Long XL6
100-0610 Medium M4	\$10.95 12/pk
100-0611 Medium M5	

COMPARE TO
**Dentatus
Posts**

ParaPost Plus

Coltene/Whaledent

Titanium Refill Posts

830-P2843Brown .036
830-P2844Yellow .040
830-P2845Red .050
\$99.99 10/bx



Stainless Steel Refill Posts

830-P2443Brown .036
830-P2444Yellow .040
830-P2445Red .050
\$99.99 10/bx

ParaPost

Coltene/Whaledent

Stainless Steel Refill Posts

830-P443Brown .036
830-P444Yellow .040
830-P4445Blue .045
830-P445Red .050
\$99.99 10/bx



Flexi Post Refill

Essential Dental Systems

235-130-0 SS #0 Yellow
235-130-00 SS #00 White
235-130-01 SS #1 Red
235-130-02 SS #2 Blue
235-135-0 Titanium #0 Yellow
235-135-00 Titanium #00 White
235-135-01 Titanium #1 Red
235-135-02 Titanium #2 Blue
\$151.95 10/pk



Flexi-Flange

Essential Dental Systems

Titanium Refill Posts

235-435-0Size 0 - Yellow
235-435-01Size 1 - Red
235-435-02Size 2 - Blue
235-435-03Size 3 - Green



Contains: 10 serrated posts w/flange,
1 primary reamer and 1 countersink drill
\$153.95

Stainless Steel Refill Posts

235-430-0Size 0 - Yellow
235-430-01Size 1 - Red
235-430-02Size 2 - Blue
235-430-03Size 3 - Green

Contains: 10 serrated posts w/flange,
1 primary reamer and 1 countersink drill
\$151.95

IntegraPost

Premier Dental

580-3001105 Size 3
580-3001110 Size 4
580-3001115 Size 4.5
580-3001120 Size 4.5L
580-3001125 Size 5
580-3001130 Size 6



\$84.95 10/pk

PINS & POSTS

Parallel-Sided Posts (Whaledent-Type)

MARK 3™

Stainless Steel Posts

100-443-SS10Size 0.036
 100-444-SS10Size 0.040
 100-4445-SS10Size 0.045
 100-445-SS10Size 0.050
 100-4455-SS10Size 0.055
 100-446-SS10Size 0.060
 100-447-SS10Size 0.070

\$35.99 10/pk

100-443-SS25Size 0.036
 100-444-SS25Size 0.040
 100-4445-SS25Size 0.045
 100-445-SS25Size 0.050
 100-4455-SS25Size 0.055
 100-446-SS25Size 0.060
 100-447-SS25Size 0.070

\$79.99 25/pk

Compare to
Coltene Para Posts

P44-3 Brown
 P44-4 Yellow
 P44-4.5 Blue
 P44-5 Red
 P44-5.5 Purple
 P44-6 Black
 P44-7 Green

\$99.99 10/pk

P44-3B Brown
 P44-4B Yellow
 P44-4.5B Blue
 P44-5B Red
 P44-5.5B Purple
 P44-6B Black
 P44-7B Green

SAVE
OVER
\$69



Titanium Posts

100-843-T10Size 0.036
 100-844-T10Size 0.040
 100-8445-T10Size 0.045
 100-845-T10Size 0.050
 100-8455-T10Size 0.055
 100-846-T10Size 0.060
 100-847-T10Size 0.070

\$35.99 10/pk

100-843-T25Size 0.036
 100-844-T25Size 0.040
 100-8445-T25Size 0.045
 100-845-T25Size 0.050
 100-8455-T25Size 0.055
 100-846-T25Size 0.060
 100-847-T25Size 0.070

\$79.99 25/pk

Compare to
Coltene Para Posts

P84-3 Brown
 P84-4 Yellow
 P84-4.5 Blue
 P84-5 Red
 P84-5.5 Purple
 P84-6 Black
 P84-7 Green

\$99.99 10/pk

SAVE
OVER
\$69

P84-3B Brown
 P84-4B Yellow
 P84-4.5B Blue
 P84-5B Red
 P84-5.5B Purple
 P84-6B Black
 P84-7B Green

Post Drills

MARK 3™

System Drills

100-423-D2Size 0.036
 100-424-D2Size 0.040
 100-4245-D2Size 0.045
 100-425-D2Size 0.050
 100-4255-D2Size 0.055
 100-426-D2Size 0.060
 100-427-D2Size 0.070

\$35.99 2/pk

Compare to
Coltene Para Posts

P42-3 Brown
 P42-4 Yellow
 P42-4.5 Blue
 P42-5 Red
 P42-5.5 Purple

P42-6 Black

P42-7 Green



ParaPost System Drills

Coltene/Whaledent

830-P423Brown .036
 830-P424Yellow .040
 830-P4245Blue .045
 830-P425Red .050

\$64.95 3/bx



**FREE SHIPPING
NATIONWIDE**

With orders over \$300
Excludes Canada. Certain restrictions apply.

Compare these exclusive brands
to name brand products

Cut your supply
costs **IN HALF!**



SILMET ProFil Compules

Microhybrid Composite
With Advanced Nano Technology

\$53.95 20/pk



Compare to Filtek Supreme **\$109.45**

20/pk

**SAVE
\$55⁵⁰**

per box

Dental **PRODUCT SHOPPER**
RECOMMENDED
2015



MARK 3TM

VPS Impression Material

4, 50ml cartridges **\$31.95**

Compare to Genie **\$59.90**

2, 50ml cartridges x 2

**SAVE
\$61⁵⁰**

per box

"I switched from Kerr Extrude to MARK3 VPS. Not only did I save money, the quality of the MARK3 Impression Material was better."

Dr. Rafael Caniza

Piermont Dental Care
Nyack, New York



PREVENTATIVES

Floss

REACH Dental Floss Professional

J&J

370-002736\$4.95 ea.

Dental Floss Dispenser

370-002733\$3.95 200/yds

Waxed Mint

370-002749\$3.95 200/yds

Waxed



REACH Wild Flosser

J&J

370-002755\$43.95

Trial Size - 144/bx



Fluoride Products

Enamel Pro Fluoride Gel

Premier Dental

580-9007200Bubblegum

580-9007201Strawberry

580-9007202Cherry

580-9007203Orange

\$30.95 16 oz./btl



MARK3™

BEST VALUE

APF Thixotropic Gel 1.23% Fluoride

100-1611Cherry

100-1612Mint

100-1613Bubble Gum

100-1614Strawberry

100-1615Grape

100-1617Vanilla - Dye Free, Neutral ph

\$11.75 17 oz bottle



COMPARE TO
DentiGel

Fluoride Trays

MARK3™

Fluoride Trays - Dual Arch

100-8821Small - Yellow

100-8822Medium - White

100-8823Large - Blue

\$17.95 100/bx



- Hinged dual/double arch
- Anatomically designed for patient comfort

Oral Rinses

Chlorhexidine Gluconate 0.12%

RX Labs

900-200114\$3.25 4oz

900-200116\$6.69 16oz



Listerine

J&J

370-30635Total Care - Fresh Mint

1 Liter Bottles **\$40.95** 6/cs

370-42855Fresh Burst

1.5 Liter Bottles **\$55.59** 6/cs



Pit & Fissure Sealants

Embrace Pit & Fissure Sealant Kit

Pulpdent

590-EMSNatural

590-EMSWOff White

\$76.95

Kit Contains: 4 - 1.2ml Syringes and 20 Applicators



MARK3™

BEST VALUE

Pit & Fissure - Light Cure

100-4540.....Opaque
\$34.95 100/bx

Kit Contains: 3 - 1.5gm Pit & Fissure Sealant Syringes,
 1 - 1.2ml Etchant Gel Syringe and 20 Tips

COMPARE TO
Clinpro



Prophy Angles

Disposable Prophy Angles

Denticator

906-501214**\$76.95**
 Firm White - 144/bx

906-501314**\$76.95**
 Soft Green - 144/bx

906-500514**\$127.95**
 Taper - 144/bx



Disposable Prophy Angles - Webbed

Young Dental

770-131120Soft Grey

770-131320Firm White

\$136.95 200/pk



Care-Free Prophy Angle

Young Dental

770-018005 **\$199.95** 5/pk



Prophy Brushes - Standard

Young Dental

770-090101 **\$99.95** 144/bx
 Standard #0213 Black

770-090201 **\$99.95** 144/pk
 Pointed #0230



Prophy Cups - Webbed Screw Type

Young Dental

770-051101 **\$56.95** 144/pk
 Soft Grey



MARK3™

BEST VALUE

Disposable Prophy Angles

100-247-144**\$46.95**
 Soft Green - 144/bx

100-248-144**\$46.95**
 Firm White - 144/bx

100-247-500.....**\$139.95**
 Soft Green - 500/bx

100-245-500**\$139.95**
 Firm White - 500/bx



COMPARE TO
**Young
 Angles**

Disposable Prophy Cups

100-442144**\$20.95**
 Latch Type - 144/bx



MARK3™
THE SMART ALTERNATIVE



Special Formula makes it a true non-splatter paste

Non-Splatter Formula
Prophylaxis Paste

BEST VALUE

FEATURES & BENEFITS

- Contains 1.23% APF
- Maximizes stain removal
- Minimizes enamel loss and helps maintain an effective fluoride level
- Strengthens tooth enamel and prevents dental caries
- Disposable cups eliminate cross contamination
- Non-Splatter formula
- Great flavor with no bitter aftertaste



order#	grit/flavor	
100-1619	Extra Coarse - Mint	100-1623 Coarse - Mint
100-1620	Coarse - Assorted	100-1630 Medium - Assorted
100-1621	Coarse - Bubble Gum	100-1633 Medium - Mint
100-1622	Coarse - Cherry	100-1640 Fine - Bubble Gum

Each box includes a prophy ring for ease of application.

\$19.95 200 per box

Compare to Nupro, Enamel Pro, Topex, Kolorz and Sparkle

Disposable

BEST VALUE

Prophy Angles

Patented prophy angle is designed to ensure each prophy procedure will require only one angle to finish.



- Patented 2nd generation gear design
- Smooth running and no overheating
- Ergonomic grips reduce potential hand fatigue
- 100% latex-free cups
- Less splatter

FREE SAMPLES Available!

RISK FREE 100% Money Back Guarantee!

100-247-144Soft	} \$46.95	144 per box
100-248-144Firm		
100-247-500Soft	} \$139.95	500 per box
100-248-500Firm		

Compare to Sultan Upgrade Angles, Denticator Angles and Young Angles

Disposable

Prophy Cups



- Premium quality prophy cups that minimize chatter and heat build up
- 100% latex free

100-442144 **\$20.95** Latch Type - 144 per box



1.23% Fluoride

APF Thixotropic Gel

FEATURES & BENEFITS

- Superior viscosity making it a true thixotropic gel
- Exceptional flavoring
- Contains Xylitol - a proven anti-carries ingredient
- Also contains Vitamin E which helps with sensitivity
- For 1 to 4 minute treatments

- 100-1611Cherry
- 100-1612Mint
- 100-1613Bubble Gum
- 100-1614Strawberry
- 100-1615Grape
- 100-1617Vanilla - Dye Free, Neutral ph



\$11.75 17 oz. bottle



Compare to Topex, Enamel Pro and DentiGel

Preventatives

MARK 3™
THE SMART ALTERNATIVE



Light Cure

Pit & Fissure Sealant



FEATURES & BENEFITS

Provides excellent marginal seal for pits and fissures on teeth showing no decay lesions. Low water absorption assures resistance to staining and durability in the oral environment.

- Does not contain Bisphenol A, Bis-GMA and Bos-DMA
- Hard setting and wear resistant
- Ideal flowable viscosity for maximum coverage and protection
- Fluoride releasing

100-4540.....Opaque

\$34.95

BEST VALUE

3 x 1.5gm Pit & Fissure Sealant Syringes
1 x 1.2ml Etchant Gel Syringe and 20 Tips

Compare to Clinpro, Embrace and Delton

5% Sodium Fluoride with TCP

Varnish

BEST VALUE

FEATURES & BENEFITS

- Treats dentin hypersensitivity
- Contains Tri-Calcium Phosphate to help remineralize teeth
- Dries to natural enamel color
- Contains Xylitol for caries prevention
- Fast and easy application

100-7100Bubble Gum

100-7101Spearmint

100-7102Caramel



\$49.95

50 Unit-dose packages, each containing 0.4ml varnish and one applicator brush



Compare to Vanish, Duraflor, DuraShield, Embrace, Enamel Pro, DMG Kolorz, Nupro White and Prevident

Prophy Pastes

Glitter Prophy Paste

Premier Dental

With Fluoride

- 580-9007400Mint - Fine
- 580-9007402Mint - Medium
- 580-9007404Mint - Coarse
- 580-9007406Mint - Extra Coarse
- 580-9007410Cherry - Fine
- 580-9007412Cherry - Medium
- 580-9007414Cherry - Coarse
- 580-9007420Strawberry - Fine
- 580-9007422Strawberry - Medium
- 580-9007424Strawberry - Coarse
- 580-9007430Bubblegum - Fine
- 580-9007432Bubblegum - Medium



Unit-Dose Cups **\$29.95** 200/bx

Ziroxide Prophy Paste

Premier Dental

Jars With Fluoride

- 580-9007174Coarse
- 580-9007175Medium

Jars Without Fluoride

- 580-9007170Coarse
- 580-9007171Medium

Mint Flavored **\$32.49** 1 lb/jar

Classic Cups

- 580-9007172Coarse
- 580-9007173Medium/Fine

\$50.95 200/bx



MARK3™

Prophy Paste

- 100-1619Extra Coarse - Mint
- 100-1620Coarse - Assorted
- 100-1621Coarse - Bubble Gum
- 100-1622Coarse - Cherry
- 100-1623Coarse - Mint
- 100-1630Medium - Assorted
- 100-1633Medium - Mint
- 100-1640Fine - Bubble Gum

\$19.95 200/bx

Includes a prophy ring for ease of application

COMPARE TO
Topex

BEST VALUE



Enamel Pro Prophy Paste

Premier Dental

- 580-9007600Mint - Fine
- 580-9007601Mint - Medium
- 580-9007602Mint - Coarse
- 580-9007603Mint - Extra Coarse
- 580-9007613Mint - Coarse w/o Fluoride
- 580-9007604Cinnamon - Fine
- 580-9007605Cinnamon - Medium
- 580-9007606Cinnamon - Coarse
- 580-9007607Strawberry - Fine
- 580-9007608Strawberry - Medium
- 580-9007609Strawberry - Coarse
- 580-9007610Grape - Fine
- 580-9007611Grape - Medium
- 580-9007612Grape - Coarse
- 580-9007614Bubblegum - Fine
- 580-9007615Bubblegum - Medium
- 580-9007616Bubblegum - Coarse
- 580-9007620Raspberry Mint - Fine
- 580-9007621Raspberry Mint - Medium
- 580-9007622Raspberry Mint - Coarse



200 unit cups and 1 Finger Ring per box **\$52.95**

PREVENTATIVES

Toothbrushes

REACH Toothbrushes

J&J

Character

370-009534\$10.95 12/bag

Batman, Youth Toothbrush

Performance

370-007212Full, Soft

370-007213Compact, Extra Soft

370-007219.....Compact, Soft

\$9.95 12/pk



Toothpaste

Plaque HD - Plaque Identifying Toothpaste

Plaque HD

101-8671.....Bubble Gum

101-8670.....Mint

Minimum purchase 1 case/12 x 4.1 oz. tubes

\$131.40 case

Plaque HD is the first plaque-identifying toothpaste. TARGETOL TECHNOLOGY™ contains disclosing agents to color and identify plaque build-up on teeth.

This allows the patient to target any areas of plaque build-up that may have been missed during brushing to completely remove plaque from teeth.



100%

**MONEY BACK
GUARANTEE**
on every product we sell

See Terms and Conditions for details.

Varnish

MARK 3™

Varnish

5% Sodium Fluoride with TCP

BEST VALUE

100-7100.....Bubble Gum

100-7101.....Spearmint

100-7102.....Caramel

Contains: 50 Unit Dose Packages each containing 0.4ml Varnish and 1 applicator Brush

\$49.95

50/bx

COMPARE TO
Nupro
White



Enamel Pro Varnish

Premier Dental

Unit Dose - 0.4 ml Bulk 200 per Pack

580-9007543Strawberries 'N' Cream

580-9007544Bubblegum

580-9007545VanillaMint

Contains: 200 Single Unit Dose Application Sticks, 1 Comfy-Grip Varnish Holder, Instructions and 35 post-Treatment Patient Instructions

\$309.95

Unit Dose - 0.4 ml Standard 35 per Pack

580-9007540Strawberries 'N' Cream

580-9007541Bubblegum

580-9007547VanillaMint

Contains: 35 Single Unit Dose Application Sticks, 1 Comfy-Grip Varnish Holder, Instructions and 35 post-Treatment Patient Instructions

\$61.95

35/bx



Hemostatic Solutions

MARK 3™

Hemostat

100-3899 **\$22.95** 30cc/btl



COMPARE TO
**Premier
Hemodent**

Hemodent Liquid

Premier Dental

580-9007071 **\$16.25** 10cc/btl

580-9007072 **\$29.95** 20cc/btl

580-9007073 **\$49.95** 40cc/btl



FS Hemostatic

Premier Dental

580-9007057 **\$21.95**

Dropper - 1 oz./btl

580-9007058 **\$33.49**

Refill - 2 oz./btl



Retraction Materials

Traxodent Hemodent Paste

Premier Dental

580-9007093 **\$86.50** Starter Pack

Includes: 7 Syringes and 15 Tips

580-9007091 **\$239.95** Value Pack

Includes: 25 Syringes and 50 Tips



Traxodent Applicator Tips

Premier Dental

580-9007092 **\$34.49** 50/pk



Retraction Caps

Premier Dental

580-9048251 Size 1 - Small
\$26.95 120/pk

580-9048253 Size 3 - Medium
\$30.95 120/pk

580-9048255 Size 5 - Large
\$22.49 60/pk

580-9048260 Assorted Kit

Contains: 15 size 1, 15 size 3 and 30 size 5

\$20.95 60/pk



Knit Pak Retraction Cord

Premier Dental

580-9007547 Cord #000

580-9007554 Cord #1

260-9007552 Cord #00

260-9007555 Cord #2

260-9007553 Cord #0

260-9007556 Cord #3

Contains: 100 inches of cord. **\$11.49** each



Gingi-Pak Soft Twist

Gingi-Pak

260-10105M Cord #1 (Thin) with dl Epinephrine

260-10110M Cord #2 (Medium) with dl Epinephrine

260-10115M Cord #3 (Thick) with dl Epinephrine

260-11105M Cord #1 (Thin) Plain

260-11110M Cord #2 (Medium) Plain

Contains: 108 inches of cord. 0.5mg Epinephrine Solution per linear inch. Package is equipped with Kutter Kap.

\$12.95 each



RUBBER DAM

Rubber Dam

Hygienic Dental Dam

Coltene/Whaledent

Rubber Latex

830-H02142\$20.95 52/pk
5"x 5" Green Medium

830-H02147\$21.95 52/pk
6"x 6" Green Medium

830-H02148\$24.95 36/pk
6"x 6" Green Heavy

830-H04641\$22.95 52/pk
Fiesta 5"x 5" Medium

830-H04642\$22.49 36/pk
Fiesta 6"x 6" Medium



Non-Latex

830-H09928\$31.25 15/pk
5"x 5" Medium

830-H09105\$35.49 15/pk
6"x 6" Medium



House Brand

Rubber Dam Clamps

101-RDC2 #2

101-RDC2A #2A

101-RDC5 #5

101-RDC7 #7

101-RDC8 #8

101-RDC8A #8A

101-RDC9 #9

101-RDC12 #12

101-RDC12A #12A

101-RDC13A #13A

101-RDC14 #14

101-RDC212 #212

101-RDCW3 #W3

101-RDCW7 #W7

\$9.99 ea.



Small Equipment

Ledex Dental Curing Light WL070

Dentmate

725-LEDEX-BLLight Blue

725-LEDEX-BLK ...Black

725-LEDEX-GRApple Green

725-LEDEX-OROrange

725-LEDEX-PKPink

725-LEDEX-WHWhite

725-LEDEX-YELYellow

\$259.99



Optical Fiber Light Guide Rod

Dentmate

725-DC030\$34.50

8mm/Black Replacement - Universal



Curing Light Sleeves

Dentmate

725-DC051\$13.49 100/pk



Ledex WL-090+ Dental Curing Light & Caries Detector

Dentmate

725-WL-090+ \$499.99

OLED display, turbo light guide rod and 7 powerful extensive modes including: Low, Ramp, Standard, High, Fast Ortho, Plaque Inspect, Caries Inspect



Dentronix DDS 7000

Dentronix

165-DDS7000..... **\$4,600.00**

Rapid Dry Heat Sterilizer



Root ZX II Apex Locator

J. Morita

860-5336392 **\$979.00**

Contains: 1 Unit, 1 Probe, 3 File Holders, 5 Contrary Electrodes, 3 AA Batteries and 1 Tester

Slow speed handpiece sold separately



AED Plus - Automated External Defibrillator

Zoll's

101-839-ZLL **\$1,495.95**

With CPR-D Pads, Carrying Case and Batteries



 **FREE SHIPPING
NATIONWIDE**

With orders over \$300
Excludes Canada. Certain restrictions apply.

EZ Kwiklave

Tuttnauer

925-EZ10K **\$5,339.99**

Fully Automatic - 220 Volt



925-EZ10P **\$5,399.99**

With Internal Printer - Documents date, time, temperature and pressure



925-EZ9 **\$4,999.99**

Fully Automatic - 120 Volt



Dynamic Automatic Mixing Unit

Shanghai NSJ

101-3020 **\$1,399.00**

- Universal for 3M, DMG, MARK3, all brands of impression material
- Fully automatic mixer, no need to keep pressing. Easy and convenient
- Choose from 3 different mixing speeds
- Low noise
- Time reminding function



Handpieces

Prophy Pal Hygiene Handpieces

Denticator

906-750001Silver

906-751001Pink

906-752001Wild Green

\$499.99 1/pack

906-750003Silver

906-751003Pink

906-752003Wild Green

906-753003Variety

\$1,399.99 3/pack



Challenger High Speed Handpieces

Johnson Promident

Standard

898-CH-SHS22 Hole

898-CH-SHS44 Hole

\$139.99 ea.

Push Button

898-CH-PBHS22 Hole

898-CH-PBHS44 Hole

\$139.99 ea.

Miniature

898-CH-MHSStandard - 4 Hole

898-CH-MPBHSPush Button - 4 Hole

\$139.99 ea.

Fiber Optic (works with 4 hole cupler)

898-CH-FOHSStandard - 5 Hole

898-CH-PBFOHSPush Button - 5 Hole

\$199.99 ea.



Classic Highspeed Handpieces

Johnson Promident

898-CL-SHS2Standard - 2 Hole

898-CL-SHS4Standard - 4 Hole

\$139.99 ea.



Gemini High Speed Handpieces

Johnson Promident

High Speed - NSK Type

898-G-VSPB**\$259.99** Kinetic Viper Standard

898-G-VSFOPB**\$309.99** Kinetic Viper Stand. Fiber Optic

898-W-NQDS**\$83.99** Quick Disconnect Swivel

898-W-NFOQDS**\$99.99** Fiber Optic Quick Disconnect Swivel

898-W-N6PQDS**\$99.99** Optic 6 Pin Swivel w/Bulb

High Speed - Kavo Type

898-G-KSTPB**\$269.95**

Super Torque 625/630/640 (Large Head)

898-G-KSTFOPB**\$369.99**

Super Torque 625/630/640 (Large Head) Fiber Optic

898-G-KBTPB**\$269.95**

Bella Torque 632/642/643 (Small Head)

898-G-KBTFOPB**\$369.99**

Bella Torque 632/642/643 (Small Head) Fiber Optic

898-W-KQDS**\$129.99** Quick Disconnect Swivel

898-W-K6PQDS**\$183.99** Fiber Optic 6 Pin Swivel w/Bulb

898-W-KFOQDS**\$183.99** Fiber Optic QuickD Swivel 5 hole

898-W-K6PQDS-LED**\$226.99** LED 6 Pin Swivel w/Bulb

Midwest Tradition Type

898-G-MWTPB**\$229.99** Push Button

898-G-MWTFOPB**\$329.99** Push Button Fiber Optic



Clear Advantage Hygienist Handpiece

Johnson Promident


CA-RDH**\$395.00**

MW/RDH Type (5,000 RPM)





Midwest Type Slow Speed Handpieces


Johnson Promident


898-MWSH\$89.99
Contra Angle Sheath for Short/Rhino 

898-MWNC\$189.99
Nose Cone for Midwest Slow Speed Handpiece 

898-MW1.....\$53.99
30,000 RPM Latch Angle 


898-MWPB40.....\$149.99
40,000 RPM Ball Bearing PB Latch Angle 


898-MW2.....\$51.99
Sealed Screw Type Prophy Angle 


898-MWFGPB.....\$169.99
40,000 RPM Push Button Friction Grip Type 

NSK/"E" Type Slow Speed Handpieces

Johnson Promident


898-ECL-4\$209.99
Classic "E" Type Motor (25,000 RPM) 


898-ESH1:1\$63.99
1:1 Contra Angle E Type Sheath Only 

898-EL30\$39.99
30,000 RPM Latch Type 

898-ELPB40\$143.99
40,000 RPM Push Button Ball Bearing Latch Type 


898-EFGPB\$159.99
40,000 RPM Push Button Friction Grip 

898-ENC.....\$174.99
1:1 Straight Nose Cone 

898-ESGN4:1.....\$154.99
4:1 Short Gripper Nose Cone 


898-R20U\$25.49
Standard U Type 20,000 RPM Contra Angle 


898-PS-1.....\$11.99
Standard Screw-on Type Prophy Angle 

898-PK-2.....\$11.99
Standard Snap-on Type Prophy Angle 


Star Type Slow Speed Handpieces

Johnson Promident


898-TCL-4\$379.99
20K Replacement Motor for Star 


898-TCL5K-4.....\$449.99
5K Replacement Motor for Star 

898-TITMA\$105.99
Motor to Angle Adapter 

898-TITNC\$199.99
Nose Cone for Star Slow Speed 

898-TIT2.....\$34.99
Sealed Screw Type Prophy 

898-TIT3.....\$34.99
Sealed Snap On Prophy 

898-TIT4.....\$53.99
30,000 RPM Latch Type 


898-TITPB40.....\$149.99
40,000 RPM Ball Bearing Push Button Latch Type 

898-TITFGPB.....\$169.99
40,000 RPM Push Button Friction Grip Type 

Star Highspeed Handpieces

Star Dental

430 SWL Fiber Optic
902-264450\$869.99 F/O Lube Free
902-264451\$1045.99 F/O Lubricated

430 SW
902-262057\$599.99 4 Hole
902-261547\$155.99
360 Degree Quick Connect Swivel 4 Line 

430K Autochuck 4 Hole
902-262052\$589.99 Autochuck 4 Hole

Ultrasonic Cleaning Systems

Quantrex Cleaning Systems

L & R

Models	Input Watts @117 VAC 60 HZ	Amps	Overall Dimensions Inches Centimeters (LxWxH)	Tank Capacity Gallons Liters	Tank Internal Dimensions Inches Centimeters (LxWxH)	Price
Quantrex 140	95 150 w/heat	0.811 1.282	10.25 X 6.5 X 8.25 26.0 X 16.5 X 21.0	0.85 3.2	9.38 X 5.38 X 4.0 23.8 X 13.7 X 10.2	\$599.99
Quantrex 210	135 335 w/heat	1.154 2.863	12.75 X 7.0 X 11.0 32.4 X 17.8 X 27.9	1.51 5.7	11.75 X 6.0 X 6.0 29.8 X 15.2 X 15.2	\$999.99
Quantrex 310	320 520 w/heat	2.735 4.444	16.5 x 10.0 x 12.13 41.9 X 25.4 X 30.8	3.25 12.3	15.5 X 9.0 X 6.0 39.4 X 22.9 X 15.2	\$1,499.99
Quantrex 360	240 440 w/heat	2.051 3.760	12.75 X 10.5 X 14.25 32.4 X 26.7 X 36.2	3.59 13.6	11.5 X 9.25 X 8.0 29.2 X 23.5 X 20.3	\$1,599.99
Quantrex 650	415 815 w/heat	3.547 6.966	21.75 X 13.75 X 12.7 55.2 X 34.93 X 32.2	6.50 24.7	19.75 X 11.75 X 6.5 50.2 X 29.8 X 16.5	\$2,200.00

All Models Include: Timer, Drain and Cover (Quantrex 90 does not include drain).

Heater optional on all systems.

All systems operate @43KHz. Systems available with international voltages.



Top Quality Curing Lights by Dentmate

One of the best known names in the industry



Ledex WL070



Choose
Your
Favorite
Color

Simple to Operate

10 Second Fast Cure

Never Shutdown Within Continuous Use

Combo of Cordless & Cord Types

Radiant intensity: 1000mw/cm² above

THE NEW COMBINATION

Curing Light & Caries Detector Ledex WL-090+

Plaque Inspect & Caries Inspect

3 Second Fast Cure

7 Powerful Extensive Modes

Never Shutdown Within Continuous Use

Built-in Radiometer

Radiant intensity: up to 1800 mW/cm²



DemeTECH SUTURES

Chromic Gut Sutures

DemeTech's Chromic Catgut suture is an absorbable, sterile, surgical suture composed of highly purified connective tissue (mostly collagen) derived from either beef or sheep intestines. DemeTech's Chromic Catgut is treated with chromium salt solution to resist body enzymes, thereby prolonging the absorption time to over 90 days. This chromic covering not only prolongs the strand's integrity, but also minimizes tissue irritation, causing less internal reaction. The percentage of collagen in the catgut suture often determines the quality of the suture. Higher percentages of collagen allow for: superior tensile strength, longer absorption times, and lower reactions in vivo.

- DemeTech manufactures its Chromic Catgut sutures with the highest quality collagen, maintaining a composition of 97%-98% of pure collagen.
- DemeTech's Chromic Catgut sutures are covered in chromic salt to prolong integrity, minimize tissue reactivity, and minimize the chance of adverse reactions.
- Tensile strength is retained for up to 14 days, and sometimes as long as 21 days.
- Allows for smooth passage through tissue.
- Easy knot tying characteristics
- DemeTech's new technology of chromicizing strands minimizes the possibility of fraying and maximizes the potential and integrity of its Chromic Catgut.

Demetech Code	Look Code	Size	Needle	Length	Price
103-DT-525-1	525B (C6)	3/0	19mm	18"	\$27.95
103-DT-636-1	560B (C6)	3/0	19mm	30"	\$27.95
103-DT-558-1	558B (C6)	4/0	19mm	18"	\$27.95
103-DT-634-1	557B (C6)	5/0	19mm	30"	\$27.95
103-DT-635-1	559B (C6)	4/0	19mm	27"	\$27.95
103-DT-687-1	1248B (C3)	5/0	13mm	18"	\$79.95

12 per pack



PGA (Polyglycolic Acid) Coated Sutures

The Polyglycolic Acid suture is a synthetic, absorbable, sterile, surgical suture composed of 100% glycolide. The chemical formula of DemeTech's PGA is (C₂H₂O₂)_n. Clinical trials have shown that after two weeks, approximately 84% of DemeTech's PGA initial strength remained. At four weeks, approximately 23% of DemeTech's PGA linear tensile strength remained. Polyglycolic Acid Sutures have been found to be nonantigenic and elicit only mild tissue reactivity during the absorption process.

- Polyglycolic Acid, better known as PGA, is composed of polymers of glycolide acids.
- DemeTech's PGA is coated to allow for easier passage through tissues with only minimal drag.
- Coated PGA sutures facilitate:
 - Easy tissue passage
 - Precise knot placement
 - Smooth tie down
 - A decreased tendency to irritate tissue
 - Greater tensile strength

Demetech Code	Look Code	Size	Needle	Length	Price
103-DT-8683-1	1022B (C6)	4/0	19mm	18"	\$44.95
103-DT-421-1	421B (C6)	5/0	19mm	30"	\$44.95
103-DT-422-1	422B (C6)	4/0	19mm	30"	\$44.95
103-DT-493-1	493B (C3)	5/0	13mm	18"	\$44.95
103-DT-492-1	492B (C3)	6/0	13mm	18"	\$89.49

12 per pack



COMPARE TO LOOK

Plain Gut Sutures

DemeTech's Plain Catgut suture is an absorbable, sterile, surgical suture composed of highly purified connective tissue (mostly collagen) derived from either beef or sheep intestines. Due to an extremely short absorption period, DemeTech's Plain Catgut elicits only minimal reaction in tissues. The percentage of collagen in the catgut suture often determines the quality of the suture. Higher percentages of collagen allow for: superior tensile strength, longer absorption times, and lower reactions in vivo.

- DemeTech manufactures its Chromic Catgut sutures with the highest quality collagen, maintaining a composition of 97%-98% of pure collagen.
- DemeTech's Plain Catgut is rapidly absorbed.
- Tensile Strength is maintained for about 10 days.
- Absorption is complete within 70 days.
- DemeTech's Plain Catgut is ideal for tissues which heal rapidly and require minimal support

Demetech Code	Look Code	Size	Needle	Length	Price
103-DT-612-1	594B (C9)	3/0	22mm	18"	\$27.95
103-DT-591-1	591B (C6)	4/0	19mm	18"	\$27.95
103-DT-592-1	592B (C6)	3/0	19mm	18"	\$27.95
103-DT-822-1	554B (C6)	3/0	19mm	30"	\$27.95
103-DT-810-1	553B (C7)	3/0	24mm	30"	\$27.95
103-DT-821-1	551B (C6)	4/0	19mm	30"	\$27.95
103-DT-686-1	596B (C3)	5/0	13mm	18"	\$79.95

12 per pack



Silk Sutures

DemeTech's Silk suture is a non-absorbable, sterile, surgical suture composed of an organic protein called fibroin. These silk sutures are processed to remove the natural waxes and gums from its origin. As Silk is a natural fiber, it does illicit an acute inflammatory reaction, which is followed by gradual encapsulation of the suture by fibrous connective tissue (occurring usually in 14-21 days). While DemeTech's silk sutures are not considered absorbable sutures, progressive degradation of its fibers results in the gradual loss of tensile strength within the first year. The silk sutures are indicated for use in general soft tissue approximation and/or ligation.

- Raw silk is produced by a silkworm larva's cocoon. These natural waxes and gums are removed prior to the braiding process, allowing for a tighter, more compact braid.
- DemeTech's Silk suture is dyed black for easier visibility in tissue.
- DemeTech's Silk suture is non-absorbable and multifilament in nature.
- Although DemeTech classifies its silk suture as a non-absorbable suture, it loses most of its tensile strength in one year, and is completely absorbed in 2 years. Thus, DemeTech's Silk behaves like a very slow absorbing suture.
- DemeTech coats its Silk in silicon to limit absorption, increase ductility, and minimize tissue reactivity.

Demetech Code	Look Code	Size	Needle	Length	Price
103-DT-7762-1	778B (C26)	3/0	15mm	18"	\$26.49
103-DT-632-1	763B (C9)	3/0	22mm	18"	\$26.49
103-DT-7772-1	784B (C6)	3/0	19mm	18"	\$26.49
103-DT-684-1	766B (C7)	3/0	24mm	18"	\$26.49
103-DT-777-1	777B (C26)	4/0	15mm	18"	\$26.49
103-DT-736-1	761B (C17)	4/0	12mm	18"	\$26.49
103-DT-683-1	759B (C6)	4/0	19mm	18"	\$26.49
103-DT-719-1	756B (C17)	5/0	12mm	18"	\$26.49
103-DT-682-1	754B (C6)	5/0	19mm	18"	\$26.49

12 per pack



Automatic Processing Solutions

Developer & Fixer Set

Air Techniques

110-40270.....\$97.95 4/cs

2 Gallons each of Developer and Fixer



Peri-Pro Developer and Fixer

Air Techniques

110-90800\$52.95 6/cs

3 Quarts each of Developer and Fixer



Bite Blocks

MARK 3™

Bite Blocks

100-7000

Fits Midwest and Ritter Panorol and J. Morita Versaview

100-7010

Fits Siemens OP 10 and Instrumentarium

100-7030

Fits Phillips, Den-Tal-Ez Panolite and Soredex

100-7040

Fits GE Panolipse and Gendex

100-7060

Fits Planmeca

100-7090


Fits Belmont X-Caliber

100-7095

Fits Panoramic Corporation

\$37.99 100/pk

COMPARE TO
Dispos-a-Bite




Bite Wing Tabs

MARK 3™

Self Adhesive Bite Wing Tabs

100-0500.....\$7.99 500/bx



- Self-adhesive tabs attach to adult and children's X-ray films
- Holds securely for accurate, distortion-free bite wing radiographs

Film Flaps Original White Cushion Tabs

Microcopy

480-FLF.....\$39.95 500/bx



Digital X-Ray Products

Sensor Flaps Long Blue Digital Sensor Tabs

Microcopy

480-FLS\$39.95 375/bx



MARK 3™

Digital X-Ray Sensor Sleeves

Plastic Sleeves

100-7798.....#1 - 1 3/8 x 8

100-7799.....#2 - 1 5/8 x 8

\$12.95 500/bx



- Custom-fitted without baggy or wrinkled plastics for clear viewing and placement
- Made of soft plastic with no sharp edges for maximum patient comfort
- The paper-backing design ensures easy insertion of sensors
- 100% Latex-Free

Film

Comfee's Deluxe Sensor Sleeves

Flow X-Ray

915-80200.....Medium

915-80210.....Large
\$57.95 500/bx



Extraoral Film

Kopy Kat Duplicating Film

Flow X-Ray

915-18070.....**\$69.95** 50/bx
 5" x 12"

915-18071.....**\$119.95** 100/bx
 6" x 12"



Intraoral Film

Flow Silver D-Speed Film

Flow X-Ray

915-18060**\$71.25** 100/bx - Bitewing
 DX-42/ Size 3 - Single Packets



915-18200**\$44.99** 150/bx
 DV-58/ Size 2 - Single Packets



915-18201.....**\$57.95** 130/bx
 DV-57/ Size 2 - Double Packets

915-18235**\$37.95** 100/bx
 DV-54/ Size 0 - Single Packets



F-Speed Xpress Film

Flow X-Ray

915-18310.....**\$48.75** 150/bx
 FV-58/ Size 2 - Single Packets



915-18410.....**\$62.49** 130/bx
 FV-57/ Size 2 - Double Packets



915-18635**\$46.99** 100/bx
 FV-54/ Size 0 - Single Packs



Ultraspeed D-Speed Film

Kodak

Poly-Soft

460-122-8840.....**\$41.95** 100/bx
 DF-54/ Size 0 - Single Packets

Poly-Soft

E460-139-4543.....**\$62.95** 130/bx
 DF-57/ Size 2 - Double Packets

Poly-Soft

E460-881-8890.....**\$52.95** 150/bx
 DF-58/ Size 2 - Single Packets



PIX D Speed X-Ray Film

PIX X-Ray

109-88194.....**\$32.75** 150/bx
 Adult D-58 Single Film Packet #2

109-88210.....**\$45.99** 130/bx
 Adult D-57 Double Film Packet

109-88400**\$29.99** 100/pk
 Child X-Ray Film D-54 Size 0 Single

BEST VALUE



PIX F Speed X-Ray Film

PIX X-Ray

109-43401.....**\$49.95** 150/bx
Adult F-58 Single Film Packet #2

109-43402.....**\$57.99** 130/bx
Adult F-57 Double Film Packet

109-88400**\$33.99** 100/pk
Child X-Ray Film F-54 Size 0 Single



Film Mounts

EZ-Grip Plastic Mounts

Flow X-Ray

915-23203.....2H - Size 2**\$18.49** 100/pk
 915-23206.....2H-BT - Size 3**\$20.99** 100/pk
 915-23209.....4H - Size 2**\$26.49** 100/pk
 915-23215.....6H-2R - Size 2**\$32.49** 100/pk
 915-23216.....6-2V - Size 2**\$36.49** 100/pk
 915-23217.....7-3V - Size2**\$36.49** 100/pk
 915-23219.....8H - Size 2**\$39.99** 100/pk
 915-23225.....14-6V - Size 2**\$54.99** 100/pk
 915-23227.....16-6V-2BT.....**\$57.99** 100/pk
 915-23239.....4H-2R - Size 2**\$25.49** 100/pk
 915-23240.....10-6V - Size 2**\$42.99** 100/pk
 915-23260.....10HBW - Size 0**\$42.99** 100/pk
 915-23261.....10HBW - Size 2**\$43.99** 100/pk
 915-23262.....10HBW - Size 3**\$42.99** 100/pk
 915-23510.....4V - Size 2 (23209)**\$24.99** 100/pk
 915-23607.....1H"0"- Size 0.....**\$17.49** 100/pk
 915-23608.....2H"0"- Size 0.....**\$18.49** 100/pk
 915-23609.....4H - Size 0**\$26.49** 100/pk
 915-23625.....14-6V"0"- Size 2**\$56.49** 100/pk



Supa Dental X-Ray Film & PSP Positioner

Flow X-Ray

915-16087.....**\$23.95** 100/bx



Greylight X-Ray Mounts

Flow X-Ray

915-601-200G1H - Size 2**\$17.29** 100/pk
 915-602-203G2H - Size 2**\$21.99** 100/pk
 915-604-209G4H - Size 2**\$28.49** 100/pk
 915-6045-239G..4H-2R - Size 2**\$28.49** 100/pk
 915-604B-510G .4V - Size 2.....**\$27.49** 100/pk
 915-606-215G6H - Size 2**\$30.99** 100/pk
 915-606B-126G .6-2V - Size2**\$30.99** 100/pk
 915-607-217G7-3V - Size 2**\$33.99** 100/pk
 915-608-219G8H - Size 2**\$36.99** 100/pk
 915-610-261G10HBW - Size 2**\$52.99** 100/pk
 915-610-262G10HBW - Size 3**\$53.49** 100/pk
 915-614-225G14-6V - Size 2**\$51.27** 100/pk
 915-614B-625G .14-6V - Size 0**\$75.99** 100/pk
 915-618-230G20HBW - Size 2**\$68.99** 100/pk
 915-702-206G2H - Size 3**\$21.99** 100/pk
 915-716-227G16-6V-2BT.....**\$54.29** 100/pk



Klear-Vue Cardboard Mounts

Flow X-Ray

915-252001H - Size 2	\$21.29	100/pk
915-252032H - Size 2	\$21.29	100/pk
915-252062N-BT- Size 3 (bitewing)	\$21.29	100/pk
915-252094H - Size 2	\$28.49	100/pk
915-252156H - Size 2	\$32.49	100/pk
915-252166H-2V - Size 2.....	\$32.49	100/pk
915-252177-3V.....	\$37.49	100/pk
915-252198H - Size 2	\$47.99	100/pk
915-2522514-6V - Size 2	\$52.99	100/pk
915-2522616-6V - Size 2	\$68.99	100/pk
915-2522716-6V-2BT.....	\$68.99	100/pk
915-2522818-6V - Size 2	\$68.99	100/pk
915-252394H-2R - Size 2	\$28.49	100/pk
915-2524010-6V - Size 2.....	\$53.99	100/pk
915-2526010BT-“0” - Size 0	\$53.99	100/pk
915-2526110BT-“2” - Size 2	\$53.99	100/pk
915-2526210BT-“3” - Size 3	\$53.99	100/pk
915-256082H-“0”- Size 0.....	\$21.99	100/pk
915-256094H-“0”- Size 0.....	\$28.49	100/pk
915-2562514-6V “0” - Size 0	\$53.99	100/pk



Perfect Pocket Mounts

Flow X-Ray

915-81-234510BT#2.....	\$52.99	100/bx
915-81-616210H-6V#2	\$89.99	100/bx
915-81-618018V6V#2	\$79.99	100/bx
915-81-60111H#2.....	\$20.99	100/bx
915-81-456720H#2	\$93.99	100/bx
915-81-40212H#0.....	\$21.99	100/bx
915-81-60212H#2.....	\$14.99	100/bx
915-81-70212H#3.....	\$29.99	100/bx
915-81-40404H#0.....	\$25.99	100/bx
915-81-60404H#2.....	\$30.99	100/bx
915-81-60454H#2-2 Row	\$37.99	100/bx
915-81-60616H#2.....	\$47.99	100/bx
915-81-60818H#2.....	\$60.99	100/bx



MARK 3™

X-Ray Film Mounts Universal Roll

100-1249.....**\$27.95** 1000/pk

- Clear vinyl pocket film mounts are easy to handle and allows for a quick glance over x-ray
- Universal x-ray film mounts are packaged into rolls, perforated, and have a white strip for film records
- Dispensed through an opening in the box



Manual Processing Solutions

Insta Neg/Insta Fix

Microcopy

480-FGC.....\$21.99

Insta Fix Gallon

480-FQC.....\$11.99

Insta Fix Quart

480-NGC.....\$21.99

Insta Neg Gallon

480-NQC.....\$11.99

Insta Neg Quart



GBX Developer and Fixer

Kodak

460-1901859.....\$29.95 Twin Pack Set

Each Twin Pack Concentrate makes
1 Gallon Developer and 1 Gallon Fixer



Phosphor Imaging Plates & Envelopes

ScanX Extraoral Phosphor Plates (PSP)

Air Techniques

110-73578-5.....\$1099.99 1/pk

Panoramic Size 5" x 12"

110-73578-6.....\$1149.99 1/pk


Panoramic Size 6" x 12"

110-73578-8.....\$1229.99 1/pk

Cephalometric Size 8" x 10"

110-73578-57.....\$699.99 1/pk

TMJ Size 5" x 7"



ScanX Intraoral Phosphor Plates (PSP)

Air Techniques

110-73445-1.....\$86.99 2/pk

Size 1

110-73445-2.....\$169.99 4/pk

Size 2

110-73445-3.....\$99.99 2/pk

Size 3

110-73445-4.....\$129.99 2/pk

Size 4



ScanX Intraoral Phosphor Plates Guides

Air Techniques

110-73566-0.....\$38.99 ea.

Size 0

110-73566-1.....\$38.99 ea.

Size 1

110-73566-2.....\$38.99 ea.

Size 2

110-73566-3.....\$38.99 ea.

Size 3

110-73566-4.....\$38.99 ea.

Size 4



ScanX Phosphor Plate Cleaning Accessories

Air Techniques

ScanX Cleaning Sheets

110-B2010.....\$159.99 12/pk

110-B2020.....\$269.99 25/bx

Phosphor Plate (PSP) Cleaning Wipes

110-B8910.....\$74.49 50/bx



MARK 3™

Intraoral Phosphor Imaging Plates

100-7780\$69.99 2/pk
Size 0

100-7781\$69.99 2/bx
Size 1

100-7782\$129.99 4/pk
Size 2

100-7783\$124.95 2/pk
Size 3

100-7784\$99.99 1/pk
Size 4



These Direct Replacement Imaging Plates are compatible with the following units:

Air Tech-ScanX, Acteon PSPX, Apixia PSP, DenOptix, Durr Vista Scan, Gendex and Owandy

- Increased compliance and patient comfort
- Clear, high quality images in seconds
- Thin, flexible, easy to position
- Surface is 100% usable
- Reusable up to 3000 times
- Increased patient comfort

Phosphor Plate Barrier Envelopes

100-7790Size 0 - 22 x 35 mm

100-7791Size 1 - 24 x 40 mm

100-7792Size 2 - 31 x 41 mm

\$34.95 300/bx

- Designed for smooth and quick loading of phosphor plates
- The large paper tab provides for an easy, secure seal
- Notched vinyl tabs tear effortlessly for easy removal of phosphor plate
- Translucent on one side to ensure proper plate position
- Fits all Phosphor Plates including Scan-X



Safe & Sure PSP Barrier Envelopes

Flow X-Ray

915-80101.....Size 0

915-80111.....Size 1

915-80121.....Size 2

915-80131.....Size 3

\$44.95 300/bx



Processing Equipment

Safelight w/Swivel & Optical Sensor

Flow X-Ray

915-19298.....**\$144.99**



Safelight w/Swivel Base - Regular

Flow X-Ray

915-19198.....**\$135.99**



Processor Cleaners

Formula 2000 Plus Cleaner

Air Techniques

110-43945**\$59.95** 2 btls/pk



Spray 2000 Plus

Air Techniques

110-43965**\$48.95**

Contains:

- 1 liter Spray, 1 Sprayer Attachment, 4 Cleaning Webs



Miscellaneous

Denture Storage Boxes

Keystone

515-9576210**\$13.49**

1 3/4" Deep
Assorted Colors
12/case



**FREE SHIPPING
NATIONWIDE**

With orders over \$300
Excludes Canada. Certain restrictions apply.

NeoDry Absorbents

Microcopy

480-DYLLarge

480-DYSSmall

\$12.75 50/bx

480-DYRLReflective Large

480-DYRSReflective Small

\$13.95 50/bx



VISIT OUR NEW,
STATE-OF-THE-ART
WEBSITE

Faster, easier ordering
www.Nobledentalsupplies.com



12 O'Clock Rear Cabinet.....52	Applicators38	Caviwipes Towelettes.....69
38% Etch Microdose.....32	ARTICULATING PAPER/ACCESSORIES15	CEMENTS & LINERS.....20-24
Abrasives & Polishers56	Articulating Paper15	Chair Covers39
Acrylic Materials8	Aspirating Syringes14	Challenger High Speed Handpieces.....90
ACRYLICS & RELINE MATERIALS.....8-9	Aspirator Tips55	Chlorhexidine 0.12%.....80
Adper Prompt L-Pop31	Assistantant's North Carolina Dental Mobile Rolling Cart53	Chromic Gut Sutures94
Adper Single Bond Plus31	Automatic Processing Solutions96	Classic Highspeed Handpieces90
AED Plus - Automated External Defibrillator89	Barrier Film.....68	Clean-A-Diamond18
Air Water Syringe Sleeves.....68	Bausch.....15	Clear Advantage Hygenist Handpiece90
Alabama Dental Rolling	Bayonet Replacement Ring.....61	Clearfil Majesty Flow.....26
Assistantant's Cart53	Bib Holder - Autoclavable Silicone38	Clearfil SE Bond30
Alginate Alternative.....58	Birex SE Disinfectant69	Clearfil SE Protect.....30
Alginate Impression Material58	Bite Block Covers68	Clear Tips Plastic Air/Water Syringe Tips.....39
Alginates.....58	Bite Blocks96	Coe Comfort.....9
Alike.....8	Bite Registration.....58-59	Coe Pack.....9
ALLOYS10-11	Bite Registration Material58	Coe Pack Automix NDS9
Allure Complete Dental Package49	Bite Registration Trays.....59	Coe Soft9
Alpha Triple Trays59	Bite Wing Tabs96	Coe Tray Cleaner62
Amalgam.....10	Bonding Systems30	Comfee's Deluxe Sensor Sleeves97
Amalgam Separators11	Brush Tips25	CompCore.....29
Amalgam Separator Waste Containers ..11	Bur Accessories18	CompCore.....34
Amalsed Conversion Kit11	Bur Blocks18	CompCore AF Cartridge STACK Dual Cure/Automix29
Amalsed Direct Dry Vacuum Amalgam Separator11	BURS & DIAMONDS16-19	Compcore AF Mixing Tips60
Amalsed Direct Wet Vacuum Amalgam Separator11	C-Fold Towels.....38	CompCore AF SyringeMix STACK Dual Cure Automix29
Anesthetic Cartridges14	Canal Plus44	CompCore AF Twist STACK Auto Cure/Syringe Hand Mix.....29
Anesthetic Cartridge Warmer & Dispenser.....14	Carbide Burs16-17	Composite Accessories25
Anesthetic Gel.....15	Care-Free Prophy Angle81	Composite Finishing57
ANESTHETICS14	Cavicide69	Composite Gun/CompuLe Dispenser.....25
Apex Locators47	Cavit G Grey.....24	Composite Polishing System.....56
APF - Thixotropic Gel80, 83	Cavit Pink.....24	COMPOSITES & RESTORATIVES25-28
	Cavity Liner Liquid.....20	
	Cavity Liners & Bases.....20	

INDEX

Compo Strips.....	56	Disposable Dappen Dishes	8	Etch Royale.....	32
Consumables	38	Disposable Impression Trays	60	Evac Shock Evacuation Cleaner	55
Core Build-Up	29	Disposable Prophy Angles - Webbed	81	EVACUATION	55
CORE MATERIALS.....	29	Disposable Prophy Angles	81, 83	Evacuation System Cleaners	55
Core Materials.....	34	Disposable Prophy Cups	81, 83	Evacuation Traps	55
COSMETIC DENTISTRY.....	30-33	DISPOSABLES.....	38-39	Evacuator Tips.....	55
Cotton Rolls	38	Distinct Latex Gloves	65	Exabite II NDS	59
Cotton Tipped Applicators	38	Drying/Obturation.....	40	Exafast NDS	62
CROWN & BRIDGE PRODUCTS.....	34-35	Dura-White Stones	57	Examix NDS.....	62
CROWN FORMS	36-37	Durelon.....	23	Excite F Refill	31
Curing Lights.....	25, 26	Dynamic (Penta) Mixer Red.....	61	Express STD Putty.....	63
Curing Light Sleeve	26	Dynamic Automatic Mixing Unit	89	Extra-Safe Earloop Masks	65
Curing Light Sleeves	88	Dynamic Yellow Mixing Tips.....	61	Extra-Safe Sensitive Earloop Masks.....	65
Curved 412 Syringe	39	Econo Wrap - Pro Vision Safety Glasses	64	Extraoral Film	97
Demetech Sutures	94-95	EDTA Solution.....	44	Eyewear	64
Dental Bibs 2-Ply.....	38	Embrace Pit & Fissure Sealant Kit	80	EZ-Grip Plastic Mounts	98
Dental Bibs 2-ply.....	68	Embrace WetBond Repair Kit	31	EZ Kwiklave	89
DenTASTIC Uno-Duo.....	30	Embrace Wetbond Resin Cement.....	23	F-Speed Xpress Film	97
Dentin Desensitizer	20	Embrace WetBond Seal-N-Shine	31	Face Masks	64-65
Dentronix DDS 7000	89	Enamel Pro Fluoride Gel.....	80	Film	97
Denture Storage Boxes	102	Enamel Pro Prophy Paste	85	Film Flaps Original White Cushion Tabs.....	96
Desensitizers.....	20	Enamel Pro Varnish	86	Film Mounts.....	98-99
Developer & Fixer Set	96	Endo Aid Kit & Refills.....	47	Filtek P-60 Syringe.....	28
Diagnostic Instruments	70	ENDODONTIC PRODUCTS	40-47	Filtek Z-100 Restorative Capsules.....	28
Diamond Burs	19	Endo Ice Spray	44	Filtek Z-100 Syringes.....	28
Diamond Twist SCL Extra-Oral Polishing Kit	57	Endo Irrigation Combo Needles.....	44	Filtek Z-250 Caps.....	28
Digital X-Ray Products.....	96	Enzymatic Cleaners	64	FINISHING & POLISHING	56-57
Digital X-Ray Sensor Sleeves.....	68	EQUIPMENT	48-54	Finishing & Polishing System	56
Digital X-Ray Sensor Sleeves.....	96	Essence Complete Dental Package.....	48	Flexi-Flange	77
DisCide Ultra Disinfectant	69	Etchant Gel.....	31	Flexi-Flow Auto E	22
DisCide Ultra Towelettes.....	69	Etching Products	31-32	Flexi-Flow Cem	22
Dispensing Gun 10:1 - Type 50	35	Etch Rite Etching Gel.....	32	Flexi-Flow Natural	22

Flexi Post Refill.....77	Hedstrom Files.....43	Intra Oral Mixing Tips.....25
Floss.....80	Heliobond Liquid.....31	Intra Oral Mixing Tips.....61
Flowable Composites.....26	Heliomolar Cavifil.....28	Intraoral Phosphor Imaging Plates.....101
Flow Silver D-Speed Film.....97	Heliomolar Flow Refill.....26	Intraoral Tips XXF.....25
Fluoride Products.....80	Heliomolar Syringe.....28	Irrigation.....43
Fluoride Trays - Dual Arch.....80	Hemodent Liquid.....87	Iso-Form Bicuspid Crowns.....37
Fluoride Trays.....80	Hemostat.....87	Iso-Form Molar Crowns.....37
Formula 2000 Plus Cleaner.....101	Hemostatic Solutions.....87	Isofluid Earloop Masks.....64
Front Surface Mirrors.....75	High Performance Mixing Tips.....60	Isolation Gown - Elastic Cuffs.....67
FS Hemostatic.....87	High Volume Vented Evacuator Tips.....55	Isolation Gown - Knit Cuffs.....67
Fuji Cem 2 Automix Paste Pak.....20	HO Bands.....76	J&J Instruments.....70-75
FujiCem 2 Mixing Tips.....21	Hold Adhesive.....62	Jackets & Gowns.....66
Fuji I.....21	HP Mixing Tips.....61	Jacket X-Safe - Hip Length.....67
Fuji II LC Capsules.....21	Hydent Indicator Spray.....15	Jacket X-Safe - Knee Length.....66
Fuji IX GP Caps Fast Set.....21	Hygienic Dental Dam.....88	Jet Tooth Shade Acrylic.....8
Fuji Lining LC Paste Pak.....20	iBrite Automix Gel Whitening System.....33	K-Files.....42
Fuji Plus Capsules.....21	iBrite LED Whitening System.....32	Ketac Cem Aplicap.....21
G-Bond Kit.....30	II Island Unit.....50	Ketac Cem Radiopaque.....22
Galil Diamond Burs.....18	II Sterilization Center.....51	Ketac Silver Aplicap.....29
GBX Developer and Fixer.....100	Implant Cement.....23	Ketac Silver Aplicap.....29
Gemini High Speed Handpieces.....90	Impregum F.....61	Klear-Vue Cardboard Mounts.....99
Genie VPS.....62	Impregum Penta.....61	Knit Pak Retraction Cord.....87
Gingi Pak Soft Twist.....87	IMPRESSION MATERIAL.....58-63	Kooliner.....9
Glass Dappen Dish.....8	Impression Trays.....60	Kopy Kat Duplicating Film.....97
Glass Ionomer Cements.....20	INFECTION CONTROL.....64-69	Kromopan Alginate.....58
Glitter Prophy Paste.....85	Infusor Tips.....25	Latex Exam Gloves.....66
Gloves.....65-66	Insta Neg/Insta Fix.....100	Ledex WL-090+ Dental Curing Light & Caries Detector.....25, 88
Goldfogel Instrument.....33	INSTRUMENTS.....70-75	Ledex WL070 Curing Light.....26, 88
Greylight X-Ray Mounts.....98	Integra Bond.....30	Lidocaine HCL 2%.....14
Hand Files.....42, 43	IntegraCem.....22	Lighting Strips - Metal.....57
Handpieces.....90, 91	IntegraPost.....77	Lignospan Forte 2%.....14
HardCore.....34	Intermediate Restorative Material.....24	Lime-Light.....20
Headrest Covers.....39	Intraoral Film.....97, 98	

INDEX

Listerine.....	80	Mobile Dental Portable Cabinet.....	53	PacFil Thermal Obturators	42
LockTight Air/Water Syringe Tips.....	39	Monobond Plus	31	PacFil Verifiers	42
Luer-Lock Endo Irrigation Syringes	43	Monoject 412 Syringe	38	Parallel Sided Posts (Whaledent Type)...	78
Luminescence	57	Monoject Needles	14	ParaPost	77
LuxaCore Dual.....	34	Multi-Cal	20	ParaPost Plus	77
LuxaFlow Ultra.....	35	Multi-Cal	44	ParaPost System Drills	78
Luxatemp Automix Plus	34	Multi-Colored Air/Water Syringe Tips	39	Patient Bibs.....	68
Luxatemp Ultra.....	34	Multilink.....	23	Penta Red Mixing Tips	61
Manual Processing Solutions	100	Needles.....	14	Perfecta Bravo Tooth Whitening	33
Matrix Bands & Strips.....	76	Needles	44	Perfecta Rev Tooth Whitening.....	33
MATRIX MATERIALS.....	76	NeoDrys Absorbents.....	102	Perfect Pocket Mounts	99
Maxcem Elite Clear	22	Neo Trays.....	59	Perform Nitrile Gloves.....	65
Medentotainer Amalgam Waste Container	11	NeutraVAC	55	Peri-Pro Developer and Fixer.....	96
Medentotainer Cap Guard.....	11	New Truliner.....	9	Periodontal Instruments.....	71
Medentotainer Lead Foils	11	NexTemp.....	24	Permadyne Garant 2:1	61
Medentotainer Tooth Box	11	Nitrile Powder Free Gloves	66	Permanent Cements.....	22
Medicaments & Canal Prep Solutions	44	Nogama	10	Permanent Resin Cement	22
Mepivacaine HCL 3%.....	14	Non-Woven N/S Gauze Sponges	38	PGA Coated Sutures	94
Meta Gutta Percha Points	41	NSK/"E" Type Slow Speed Handpieces	91	Phosphor Imaging Plates & Envelopes	100
Meta Paper Points	40	Occlude Aerosol Powder.....	15	Phosphor Plate Barrier Envelopes.....	101
Micro Applicator Brushes	32	Occlusal & Pressure Indicators.....	15	Pic-n-Stic.....	25
Micro Applicators	32	OneGloss	57	PINS & POSTS.....	77-78
Microbrush Disposable Applicators.....	32	Optibond Solo Plus	30	Pit & Fissure - Light Cure.....	81
Midwest Type Slow Speed Handpieces	91	Optical Fiber Light Guide Rod.....	26, 88	Pit & Fissure Sealants.....	80-81, 84
Miracle Mix.....	29	Oral Rinses	80	PIX D Speed X-Ray Film	97
MISCELLANEOUS.....	102	Orthodontic Instruments	70	PIX F Speed X-Ray Film.....	98
Mixing Bowls.....	61	Pac-Dent Gutta Percha Points.....	41	Plain Gut Sutures.....	95
Mixing Pads.....	25	Pac-Dent Paper Points.....	40	Plaque HD -Plaque Identifying Toothpaste.....	86
Mixing Pads	61	Pac-File Flex Rotary Files.....	45	Plastic 5 oz. Cups	38
Mixing Tools.....	61	Pac-Taper Flex Rotary Files	46	Plastic Hub Needles	14
Mixing Wells	25	PacEndo Air/Water Irrigator	47	Poli Pro Disks	
Mixing Wells	61	PacEndo Pre-Filled Irrigation Kit.....	43	Polishing Pastes.....	57

Polycarbonate Crowns	36	Prophy Angles.....	81	Retraction Materials	87
Polyether.....	61	Prophy Brushes - Standard.....	81	RootPro Wireless Endo Motor	47
Pop-On Mandrels	56	Prophy Cups - Webbed Screw Type	81	Root ZX II Apex Locator	47
Porcelain Etch Ge	32	Prophylaxis Paste.....	85	Root ZX II Apex Locator.....	89
Porcelain Prep Kit	32	Prophy Pal Hygiene Handpieces	90	Rotary Files	45
Post Drills	78	Prophy Pastes	85	RUBBER DAM.....	88
Posts	77	Protege Nitrile Exam Gloves.....	65	Rubber Dam Clamps	88
Pre-Bent Needle Tips.....	25	Protomp 4 Garant	35	Safe & Sure PSP Barrier Envelopes.....	101
Premier Etch.....	32	ProVeneer Resin Based Blister Refills.....	33	Safelight w/ Swivel & Optical Sensor ...	101
Premium Bite Registration Trays.....	59	ProVeneer System Kit.....	33	Safelight w/Swivel Base - Regular.....	101
Premium Extra Soft Earloop Masks.....	65	Pulp Capping Paste.....	44	Saliva Ejectors.....	55
Prep-Rite RC.....	44	Pulp Vitality Endo Refrigerant Spray.....	44	ScanX Extraoral Phosphor Plates (PSP)	100
Pressure Indicating Paste	15	Purity Waterline Cleaning Tablets	64	ScanX Intraoral Phosphor Plate Guides.....	100
PREVENTATIVES.....	80-86	Quantrex 140	92	ScanX Intraoral Phosphor Plates (PSP)	100
Primary Molar Micro Crowns	37	Quantrex 210	92	ScanX Phosphor Plate Cleaning Accessories	100
ProBase Liner/Base	20	Quantrex 310	92	Screw Posts - Gold Plated.....	77
Procedural Earloop Masks	64	Quantrex 360	92	Sealapex.....	46
Processing Equipment	101	Quantrex 650.....	92	Sealers.....	46
Processor Cleaners.....	101	Quick Tips Metal Air/Water Syringe Tips.....	39	Self Adhesive Bite Wing Tabs	96
ProEtch	32	RC Prep.....	44	Self Seal Sterilization Pouches	68
ProFil Bulk Fill.....	27	Reach Character Toothbrush	86	Sensor Flaps Long Blue Digital Sensor Tabs.....	96
ProFil Flow	26	REACH Dental Floss Professional	80	Septocaine HCL 4%	14
ProFil Syringe Hybrid Composite	27	Reach Performance Toothbrush	86	Septoject Needles.....	14
ProFil Temp Grey	24	Reach Wild Flosser	80	Short HP Mixing Tips	60
ProFil Unidose Hybrid Composite.....	27	Reamers	43	Silane Bond Enhancer.....	30
ProGlass Nine Capsules	21	Reline Materials	9	Silk Sutures	95
ProGlass One Kit.....	21	RelyX Luting Cement.....	22	SMALL EQUIPMENT.....	88-92
ProGlass Plus Capsules.....	21	RelyX Luting Plus.....	22	Sof-Lex Discs.....	56
ProGlass Plus Kit.....	21	RelyX Unicem Aplicap	23	Softone	9
ProGlass Silver	29	Restorative Instruments.....	71-72	Spherodon +.....	10
ProLink Self Etch.....	30	Retraction Caps.....	87		
ProLink Universal Adhesive.....	30	RETRACTION MATERIALS.....	87		
Pro One Sonictab	64				

INDEX

Spherodon-M.....	10	Temporary Crown & Bridge Mixing Tips.....	61	Tubliseal	46
Spray 2000 Plus	101	Temporary Crowns	37	Ultra Sensitive Earloop Masks	65
Stainless Steel Crowns	36	Temporary Filling Material.....	24	Ultrasonic Bioenzymatic Cleaner.....	64
Stainless Steel Posts	78	Tetric EvoCeram Cavifil.....	28	Ultrasonic Cleaning Systems	92
Stain Out Tarter & Stain Remover	64	Tetric EvoCeram Syringe.....	28	Ultraspeed D-Speed Film.....	97
Star Highspeed Handpieces.....	91	Tetric EvoFlow Cavifill.....	27	Universal Composites	27
Star Type Slow Speed Handpieces	91	Tetric EvoFlow Syringe	26	Universal VPS Adhesive.....	62
Sterilization Pouches	68	Ti-Core Gray	34	Vac Attak Evacuation System Cleaner....	55
Strips.....	57	Ti-Core Natural	34	Vacuosol Ultra	55
Supa Dental X-Ray Film & PSP Positioner.....	98	Tissue Conditioner	9	Varnish	84, 86
Supergel	58	Titanium Posts	78	Veneer System	33
Supergel Fresh	58	Tofflemire Bands.....	76	Vibrant Micro-Textured Latex Exam Gloves.....	65
Surface Barriers.....	39	Toothbrushes	86	Vinyl Polysiloxane.....	62
Surface Barriers.....	68	Toothpaste.....	86	Vitrebond	20
Surface Disinfectants	69	Tooth Whitening	32-33	Vitribond Plus Clicker.....	20
Surgical Instruments.....	72-75	Topicale Gel Pump	15	Vivacaine HCL 0.5%.....	14
Surgical Tips.....	55	Topicale Xtra.....	15	VPS 380 Impression Material	63
SUTURES.....	94-95	Topical Gel Patch	15	VPS Adhesive	62
Sycamore Wedges.....	76	Topicals.....	15	VPS Impression Material	63
Syringe Mix Mixing Tips	60	Touch 200 Nitrile Exam Gloves.....	66	VPS Putty.....	63
Syringes & Tips - Disposable.....	39	Tra-Ten Trays.....	60	Wedges	76
Syringe Tip Covers	39	Traxodent Applicator Tips	87	White Stones	57
System Drills	78	Traxodent Hemodent Paste.....	87	Wizard Wedges.....	76
T-Bands.....	76	Tray-Aways	60	X-RAY	96-101
T-Style Sleeves	39	Tray Adhesives.....	62	X-Ray Film Mounts Universal Roll	99
TAC Adhesive Spray.....	62	Tray Cleaner	62	Ziroidex Prophy Paste	85
TempBond	23	Tray Covers Ritter "B"	38		
Temporary Cement NE	23	Tray Sleeves.....	38		
Temporary Cements	23	Trim II	8		
Temporary Crown & Bridge.....	35	Triple Tray	59		
Temporary Crown & Bridge.....	34-35	Truliner Original	9		

TERMS & CONDITIONS

Noble Dental Supplies is dedicated to providing quality brands at the best prices. While we strive to maintain the present prices in our catalog, we reserve the right to adjust them in accordance with manufacturer's increases without prior notice. We are not responsible for typographical errors.

Payment Methods:

Your orders can be paid via Mastercard, American Express, Visa, PayPal, or COD. COD checks should be made out to Noble Dental Supplies. We do not offer deferred payment terms.

Delivery Terms:

Continental U.S. - All orders are subject to a handling charge. This charge includes freight, except for additional carrier charges related to special delivery services and hazardous material shipments.

Outside the Continental U.S. - Standard UPS charges will apply. Any additional charges will be the responsibility of the recipient.

Controlled Substances:

Regulations allow us to sell controlled substances to registered, licensed healthcare professionals only. If you have recently changed your address or are a new customer, please supply us with a copy of your state registration. Mail, email or fax a copy of your DEA certificate.

Mail: Noble Dental Supplies
135 North Route 9W, Congers, New York 10920

Email: anthonyt@nobledentalsupplies.com

Fax: 845-510-3191

Returns:

PRIOR AUTHORIZATION IS NEEDED FOR RETURNS. Simply contact our Customer Service Department to arrange for a return. Please comply with the following conditions:

- Returns must be accompanied by a copy of your invoice and a reason for the return.
- All returns must come back in their original, unmarked packaging.
- Returns must have been purchased within the previous 30 days. Any returns over 30 days are subject to a restocking fee.
- Shortages or errors in shipments must be reported within 7 days of the invoice date to receive credit.
- Shipping charges apply on all returns.

SATISFACTION GUARANTEED

If the product you have ordered is defective or does not perform to satisfaction, you may choose to receive a credit, exchange or refund.



WE HELP YOU SAVE

Since 1978, we've brought dentists high quality dental supplies that are 30-50% below the expensive name-brands.

